







Hi, my name is Saffi, I'm the **SafMed** Product Specialist

My role is to assist anyone using **SafMed** Surgical products to identify & understand the value they bring into the Operating Theatre.

Symbols used in this document

	<p>Instructions for Use</p> <p>Always consult IFU's (equipment manuals or package inserts) for detailed manufacturer's instructions for use</p>
	<p>Quality Assurance</p> <p>Anywhere a standard is referred to or the product can form part of your quality process</p>
	<p>Did you know?</p> <p>Useful information</p>
	<p>Take Note!</p> <p>Take special note of interactions</p>

Surginox



The **Surginox** is designed to provide complete patient positioning and flexibility required

Why do I need this table?

- The SURGINOX Surgical table is a **mobile hydraulically and mechanically operated surgical table** that is designed to support all general surgical procedures including cardiac, endoscopic, gynaecology, urology, nephrectomy, neurology, ophthalmologic and orthopaedics

Features of the table

- The SURGINOX Surgical Table features hydraulically driven lateral tilt, Trendelenburg / reverse Trendelenburg, and adjustable height functions, and mechanically back, kidney, slide, head and legs sections
- The headrest, back section, kidney elevator, leg sections and slide are manually adjustable
- This table is designed to safely function with a 150kg patient

How do I use the Surginox?

- When adding the head rest, back rest and leg section, place the pins against the holes of the table until it clicks into place



The **headrest, back rest and leg section** is equipped with **gas spring cylinders** that **reduce the effort required to raise or lower** it while a patient is lying on the table

- Activate the table floor lock before positioning the patient and starting a procedure
- Do not exceed 150kg total weight of patient + accessories for the standard positions
- Always plug the table in at the end of the day
- Never leave anything on the foot piece/base or next to the column of the table
- Ensure proper cleaning between cases and at the end of the day
- Remember to clean the accessories before putting it on a trolley or into storage

- Do not use chlorine-based products or other products containing chlorine derivatives; these can cause excessive corrosion
- Ensure the theatre tables brakes are on when transferring a patient
- Service your theatre table at the recommended intervals to ensure continued functioning



Always refer to the instructions for use (manual)