







Hi, my name is Saffi, I'm the **SafMed** Product Specialist

My role is to assist anyone using **SafMed** Surgical products to identify & understand the value they bring into the Operating Theatre.

Symbols used in this document

	<p>Instructions for Use</p> <p>Always consult IFU's (equipment manuals or package inserts) for detailed manufacturer's instructions for use</p>
	<p>Quality Assurance</p> <p>Anywhere a standard is referred to or the product can form part of your quality process</p>
	<p>Did you know?</p> <p>Useful information</p>
	<p>Take Note!</p> <p>Take special note of interactions</p>

STERIS 4085



The **STERIS 4085** has a large sliding top for unrestricted radiological access

Why do I need this table?

- The STERIS 4085 is a mobile, electro-hydraulically operated table that is designed to support all general surgical procedures including cardiac and vascular, endoscopic, gynaecology, urology, nephrectomy, neurology, ophthalmology and orthopaedics

Features of the table

- This table features powered lateral tilt, Trendelenburg/ reverse Trendelenburg, a **movable table-top** and adjustable height functions
- This table is designed to safely function with a 272kg patient

How do I use the STERIS 4085?

- Accessories may have a specified weight limitation that is less than that of the table. Do not exceed the maximum patient weight
- Do not place the patient on the table unless the floor locks are engaged
- Do not disengage the floor locks when the patient is on the table
- When manually positioning the table sections, stay clear of potential pinching hazards
- **Centre the table-top over the column before applying chest compressions**
- Centre the table-top over the column **before transferring** the patient on or off the table
- Always plug the table in at the end of the day
- Ensure that the remote control is also properly connected to charge its battery
- Never leave anything on the foot piece/base or next to the column of the table



Ensure proper cleaning between cases and at the end of the day

- Remember to clean the accessories before putting it on a trolley or into storage
- Do not use chlorine-based products or other products containing chlorine derivatives; these can cause excessive corrosion
- Ensure the theatre tables brakes are on when transferring a patient
- Service your theatre table at the recommended intervals to ensure continued functioning



Always refer to the instructions for use (manual)