







Hi, my name is Saffi, I'm the **SafMed** Product Specialist

My role is to assist anyone using **SafMed** Surgical products to identify & understand the value they bring into the Operating Theatre.

Symbols used in this document

	<p>Instructions for Use</p> <p>Always consult IFU's (equipment manuals or package inserts) for detailed manufacturer's instructions for use</p>
	<p>Quality Assurance</p> <p>Anywhere a standard is referred to or the product can form part of your quality process</p>
	<p>Did you know?</p> <p>Useful information</p>
	<p>Take Note!</p> <p>Take special note of interactions</p>

Cmax X-Ray



The **Cmax X-Ray** offers a powered operating table for image assisted surgery

Why do I need this table?

- The Cmax X-Ray operating tables are designed to support and position adult patients during operations carried out in an operating theatre, and to transfer them from their bed
- The Cmax X-Ray is suitable for **imagery-assisted surgical procedures**

Features of the table

- The Cmax X-Ray operating table is equipped with a **carbon fibre** floating table-top and is operated in standalone mode and can be used for powered patient transfer
- This table is designed to safely function with a 230kg patient

How do I use the Cmax X-Ray?

- Activate the table floor lock before positioning the patient and starting a procedure
- It is strongly recommended that you use barriers or another system to hold the patient in place
- We advise the user to practise operating the table motorisation without a patient first
- Do not exceed 230 kg total weight of patient + accessories for the standard positions
- Always plug the table in at the end of the day
- It is recommended to re-connect the joystick every day after use, to ensure that it functions correctly and is **charged in Bluetooth mode**
- Never leave anything on the foot piece/base or next to the column of the table
- Ensure proper cleaning between cases and at the end of the day.
- Remember to clean the accessories before putting it on a trolley or into storage
- Do not use chlorine-based products or other products containing chlorine derivatives; these can cause excessive corrosion

- Ensure the theatre tables brakes are on when transferring a patient
- Service your theatre table at the recommended intervals to ensure continued functioning



Always refer to the instructions for use (manual)