







Hi, my name is Saffi, I'm the **SafMed** Product Specialist

My role is to assist anyone using **SafMed** Surgical products to identify & understand the value they bring into the Operating Theatre.

Symbols used in this document

	<p>Instructions for Use</p> <p>Always consult IFU's (equipment manuals or package inserts) for detailed manufacturer's instructions for use</p>
	<p>Quality Assurance</p> <p>Anywhere a standard is referred to or the product can form part of your quality process</p>
	<p>Did you know?</p> <p>Useful information</p>
	<p>Take Note!</p> <p>Take special note of interactions</p>

STERIS 3085



The **STERIS 3085** offers easy C-arm access and complete upper and lower body Image Amplification

Why do I need this table?

- The STERIS 3085 is a mobile, electro-hydraulically operated general surgical table providing flexible articulation of the patient
- This table accommodates all surgical procedures including Vascular, Cardiac, Thoracic, Endoscopy, ENT, Neuro, Ophthalmology, Arthroscopy, Spine Rectal, Gynaecology, Arm & Hand Orthopaedic Surgery

Features of the table

- A microprocessor controlled hydraulic system provides table-top articulation
- The base is cast iron with textured enamel finish
- The top is enclosed by a two-section, welded, stainless-steel cover that also forms a shroud for the lower portion of the column
- **Four swivel casters** facilitate easy table relocation and movement
- This table is designed to safely function with a 226kg patient

How do I use the STERIS 3085?

- Activate the table floor lock before positioning the patient and starting a procedure
- Do not exceed 226kg total weight of patient + accessories for the standard positions
- When performing surgery requiring a headrest accessory in a reversed patient orientation, do not exceed 181kg patient weight limit



You **should** always **plug the table in at the end of the day**

- Never leave anything on the foot piece/base or next to the column of the table
- Ensure proper cleaning between cases and at the end of the day

- Remember to clean the accessories before putting it on a trolley or into storage
- Do not use chlorine-based products or other products containing chlorine derivatives; these can cause excessive corrosion
- Ensure the theatre table brakes are on when transferring a patient
- Service your theatre table at the recommended intervals to ensure continued functioning



Always refer to the instructions for use (manual)