



Your solutions partner for **30 YEARS**

SUMA L4 ALKALINE DETERGENT

Instructions for use

Product description

Suma L4 is an Alkaline Detergent suitable for use in medium to hard water conditions and is designed to remove fatty and proteinaceous soiling. Suma L4 is particularly effective against staining.

Container size

Description	Size
Prolystica 2X Enzymatic Cleaner	3.8 litre

Features	Benefits
<ul style="list-style-type: none"> High quality sequestrants 	<ul style="list-style-type: none"> Provides water softening properties and enhances detergency
<ul style="list-style-type: none"> Concentrated cleaning action 	<ul style="list-style-type: none"> Prevents scale build-up
<ul style="list-style-type: none"> Automatically dispensed 	<ul style="list-style-type: none"> Removes organic material and leaves stainless steel surgical equipment clean.
<ul style="list-style-type: none"> Alpha numeric label 	<ul style="list-style-type: none"> Cost effective
	<ul style="list-style-type: none"> No caking or spilling.
	<ul style="list-style-type: none"> No waste
	<ul style="list-style-type: none"> Easy identification – assists staff training.

Wear personal protective equipment (gloves, goggles and aprons) when using Prolystica. Container size

Indications for use

Suma L4 is suitable for use on steel, stainless steel, copper, plastics, rubber, and ceramics. Suma L4 Alkaline Detergent is a highly alkaline liquid detergent, suitable for use in automated cleaning processes in the medical industry

Storage

Store in accordance with local and national regulations. Keep only in original container. Store in a closed container.

Direction for use

Step 1: Confirm Suma L4 can be used with the make and model of Washer Disinfector

Consult with the supplier of the Washer Disinfector and detergent to confirm that Suma L4 has been validated for use in the particular make and model of Washer Disinfector

Step 1: How much Suma must be used?

Request Washer Disinfector technician to dose the detergent at 1-3ml per litre of water.