



Surgical Table Accessories Catalog

Accessories by Category

Accessory Carts	5
Accessory Packages	7
Arm Supports	14
Body Restraints & Safety Rails	19
Body Supports	21
Clamps & Sockets	24
Disposables / Cleaning Supplies	26
Head Supports	27
Leg Supports	31
Patient Positioner Pads	35
Table Operating Controls	40
X-ray Tops	44

Accessories by Surgical Service

Anesthesia	47
Bariatric / Gastrointestinal	48
Neurosurgical	59
Ophthalmology/ENT/Plastic	63
Orthopedic	65
Radiology	74
Urology / Gynecology	77
Vascular	81

Indices

Index by Item Description	83
Index by Equipment Number	88

How to Order

On the following pages you will find specifications for accessories to be used on STERIS's complete line of surgical tables: STERIS® 5085 SRT, 5085, and 4085 tables; Cmax™ tables, Amsco® 3085 SP, 3080-R and 3080 tables; OrthoVision® tables, and the SurgiGraphic® 6000 and 1027 Image Guided Tables. For more information about STERIS's surgical accessories, contact your STERIS Account Manager, or call 888-478-3747. In Canada, call 800-661-3937.

Now you can place your order online at www.steris.com. Or you can call our toll-free number, 888-478-3747, Monday through Friday, between 7 a.m. and 5 p.m. (EST) and our Inside Sales Representatives will be glad to assist you.

Thank you for your business.

Total Patient Support

Each STERIS surgical table has its own total patient weight limitation. **Accessories and accessory systems also have total patient weight limitations based on the proportional weight of the intended body part(s) supported by the accessory.**

When more than one accessory is used to support a patient body part, the overall rated load will be that of the weakest accessory in the system or the surgical table to which it is applied. Known total patient weight limitations for surgical tables and accessories are identified throughout this catalog. All personnel should be informed of these limitations.

Refer to the Operator's Manual supplied with each table for more information and proper use instructions.

STERIS Warranty

STERIS Corporation is committed to customer satisfaction. STERIS will correct any defects(s) in material or workmanship occurring during the warranty period in accordance with the following terms:

- For 12 months after the warranty start date, STERIS Corporation will provide free replacement parts for all equipment (except surgical table pads and hand controls) that does not perform according to its specifications.
- For 90 days after the warranty start date, STERIS Corporation will provide free replacement parts for surgical table pads and surgical table hand controls that do not conform to STERIS's specifications.

All work on STERIS equipment during the warranty period must be performed by a qualified technician. For more information, see STERIS's Warranty and Customer Assistance Information booklet, P-056399-116.

STERIS® HLT Pads

Prolonged exposure to high pressure during surgery, particularly over bony prominences, can lead to ischemia and necrosis of the skin. Damage occurs at the cellular level and is often irreversible. Because ulcers present themselves hours or days after patients leave surgery, the OR staff is often unaware of the damage.

The prevalence of OR-acquired pressure ulcers is an increasing concern for hospitals. There are many risk factors contributing to patient susceptibility to pressure related injury including the trend toward increases in weight and age. As patient risk factors increase, the potential for OR-acquired pressure ulcers becomes a greater clinical concern. Generally, table pads feature the same configuration: uniform foam layers of differing densities. That means every part of the body receives the same protection.

Through extensive research, clinical partnerships and anthropometric patient scenarios, STERIS has developed an innovative Hybrid Layer Technology, engineered to redistribute pressure in specific areas commonly associated with OR-acquired pressure ulcers. The table pad's unique geometry optimizes load distribution in areas such as the head, sacrum and heels, redistributing weight to protect these bony prominences. This clinically-inspired table pad construction provides hospitals with a different solution.

STERIS® TLT Pads

The TLT has been developed as an overall table surface and a STERIS standard since 1990. The TLT has a history of proven performance using the Triple Layer Technology in providing overall patient protection.

TLT pads protect patients from both shear and direct pressure forces. The Fusion Two-Way Stretch vinyl fabric Herculite cover moves relative to the patient skin surface to reduce shear pressure. The cover material is fluid-proof, bacteria-resistant, flame-retardant, self-deodorizing, and electrically conductive. It is also non-toxic and non-allergenic. The cover material is on Underwriters Laboratories approval list for conductive sheeting and meets National Fire Prevention Association (NFPA) requirements. Seams are sealed in normally wet areas to prevent absorption. The bottom of the pad is durable anti-static Herculite Lectrolite cover stock.

The foam core is comprised of a proprietary construction of three layers of specialized foam:

- The top layer of soft slow-recovery foam conforms to patient weight and skin temperature to assure patient comfort.
- The middle layer of cushioning foam spreads forces for even distribution of weight.
- The bottom layer of high density foam suspends the patient safely above the hard surface of the surgical table top.

STERIS's entire pad assembly does not contain Latex material.

General Surgical Tables		2" TLT	4" TLT	HLT
5085 SRT	Complete Set, 3 Pcs	BF095	BF295	BF383
	Head Section	BF096	BF296	BF384
	Back/Seat Sections	BF582	BF293	BF387
	Back/Seat Sections w/o Perineal Cutout	BF83400206	N/A	N/A
	Leg Section	BF583	BF294	BF388
	Leg Section w/o Perineal Cutout	BF83300339	N/A	N/A
5085 4085 / Cmax	Complete Set, 3 Pcs.	BF58	BF291	BF385
	Complete Set, 3 Pcs with Rectangular Leg Section <i>(for use with Featherweight Leg Section)</i>	BF65	N/A	N/A
	Head Section	BF581	BF292	BF386
	Back/Seat Sections	BF582	BF293	BF387
	Back/Seat Sections w/o Perineal Cutout	BF83400206	N/A	N/A
	Leg Section	BF583	BF294	BF388
3085SP 3080	Complete Set, 3 Pcs	BF321	BF112	BF389
	Head Section	BF321001	BF113	BF391
	Back/Seat Sections	BF321002	BF114	BF392
	Back/Seat Section w/o Perineal Cutout	BF83300206	N/A	N/A
	Leg Section	BF321003	BF115	BF393

Specialty Tables	2" TLT	4" TLT	HLT
Ortho 2, OrthoVision, and OrthoView	BF323	N/A	N/A
SurgiGraphic 1027	TA001	N/A	N/A
SurgiGraphic 6000	TA00007	N/A	N/A

Table Accessory Cart

#BF7

This stainless steel cart features four shelves to accommodate storage of a full variety of accessories and is essential for an organized, productive OR. The four casters allow it to be moved from room to room easily to adapt to your surgical schedule. (Accessories not included.)



Orthopedic Extension Cart

#BF24


This cart was specifically designed with "risers" on the top of the cart to hold the abductor bar assembly of the Orthopedic Extension. These small stainless steel platforms allow OR staff to wheel the cart right up to the table for quick, easy and safe attachment of the Orthopedic Accessory. The cart also features four shelves to organize all related accessories. (Accessories not included.)



A-Frame Accessory Cart

#BF27

This mobile peg board, constructed of 16-gauge stainless steel, is extremely durable and easily maintained and cleaned. It comes equipped with 24 single, 12 double, and 4 "J" hooks, enabling accessories of various sizes and shapes to be attached for easy access or storage. The cart also contains a storage bin and two storage shelves located between the cart walls. (Accessories not included.)



Bariatric Accessory Storage Cart

#BF551

This cart is a compact storage and transport device for the full range of STERIS bariatric accessories. The eight side rails were specifically designed to accommodate each of the unique Bariatric Width Extenders. The two back rails easily hold STERIS's Split Leg Positioner or other leg holder accessories. Other accessories may be attached to interior rails or stored in the deep bin / weighted base. (Accessories not included.)



Power-lift / Bariatric Power-lift Stirrup Storage Cart

#BF46

This cart provides stable, lightweight, portable storage for your surgical stirrups, table clamps and other leg holder accessories; the pre-configured handles allow safe and easy transport. The cart was designed with an integral set of support prongs that accommodate either the Power-lift Stirrups or the Bariatric Power-lift Stirrups. (Accessories not included.)



Shoulder Table Accessory Cart

#BF461

This dual purpose storage cart was specifically engineered to hold the Shoulder Table or the Shoulder Chair and related accessories, and has a weighted base that makes this cart tip-resistant. The ergonomic design helps ensure that devices are always at a convenient height for the user to reduce potential back strain. The optional accessory hooks allow quick and easy organization and storage of components and posturing accessories, and the deep bin in the base provides additional valuable storage space. (Accessories not included.)



Standard Accessory Package with Foot Extension, 3-inch Armboard (to be used with X-ray tops), and Standard Length Restraint Strap 66-1/2"

#BF564

Includes:

- (2) Anesthesia Armboards with 3-inch TLT Pads (BF804)
- Perineal Cutout Filler Piece (BF205)
- Cable Wrap (BF316)
- Foot Extension with Siderail Locks (BF2)
- Standard Length Restraint Strap 66-1/2" (BF45)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Standard Accessory Package with Foot Extension, 3-inch Armboard (to be used with X-ray tops), and Extra-long Restraint Strap 126"

#BF565

Includes:

- (2) Anesthesia Armboards with 3-inch TLT Pads (BF804)
- Perineal Cutout Filler Piece (BF205)
- Cable Wrap (BF316)
- Foot Extension with Siderail Locks (BF2)
- Extra-long Restraint Strap 126" (BF451)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Standard Accessory Package with Foot Extension, 2-inch Armboard (to be used without X-ray tops), and Standard Length Restraint Strap 66-1/2"

#BF566

Includes:

- (2) Anesthesia Armboards with 2-inch TLT Pads (BF803)
- Perineal Cutout Filler Piece (BF205)
- Cable Wrap (BF316)
- Foot Extension with Siderail Locks (BF2)
- Standard Length Restraint Strap 66-1/2" (BF45)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Standard Accessory Package with Foot Extension, 2-inch Armboard (to be used without X-ray tops), and Extra-long Restraint Strap 126"

#BF567

Includes:

- (2) Anesthesia Armboards with 2-inch TLT Pads (BF803)
- Perineal Cutout Filler Piece (BF205)
- Cable Wrap (BF316)
- Foot Extensions with Siderail Locks (BF2)
- Extra-long Restraint Strap 126" (BF451)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Standard Accessory Package, 3-inch Armboard (to be used with X-ray tops), and Standard Length Restraint Strap 66-1/2"

#BF413

Includes:

- (2) Anesthesia Armboards with 3-inch TLT Pads (BF804)
- Perineal Cutout Filler Piece (BF205)
- Cable Wrap (BF316)
- Standard Length Restraint Strap 66-1/2" (BF45)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Standard Accessory Package, 3-inch Armboard (to be used with X-ray tops), and Extra-long Restraint Strap 126"

#BF414

Includes:

- (2) Anesthesia Armboards with 3-inch TLT Pads (BF804)
- Perineal Cutout Filler Piece (BF205)
- Cable Wrap (BF316)
- Extra-long Restraint Strap 126" (BF451)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Standard Accessory Package, 2-inch Armboard (to be used without X-ray tops), and Standard Length Restraint Strap 66-1/2"

#BF424

Includes:

- (2) Anesthesia Armboards with 2-inch TLT Pads (BF803)
- Perineal Cutout Filler Piece (BF205)
- Cable Wrap (BF316)
- Standard Length Restraint Strap 66-1/2" (BF45)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Standard Accessory Package, 2-inch Armboard (to be used without X-ray tops), and Extra-long Restraint Strap 126"

#BF425

Includes:

- (2) Anesthesia Armboards with 2-inch TLT Pads (BF803)
- Perineal Cutout Filler Piece (BF205)
- Cable Wrap (BF316)
- Extra-long Restraint Strap 126" (BF451)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Standard Accessory Package with Foot Extension, and 3-inch Armboard (to be used with X-ray tops)

#BF093

Includes:

- (2) Anesthesia Armboards with 3-inch TLT Pads (BF804)
- Perineal Cutout Filler Piece (BF205)
- Cable Wrap (BF316)
- Foot Extensions with Siderail Locks (BF2)

Table Compatibility: 5085 SRT

Standard Accessory Package with Foot Extension, and 2-inch Armboard Pad (to be used without X-ray tops)

#BF099

Includes:

- (2) Anesthesia Armboards with 2-inch TLT Pads (BF803)
- Perineal Cutout Filler Piece (BF205)
- Cable Wrap (BF316)
- Foot Extension with Siderail Locks (BF2)

Table Compatibility: 5085 SRT

Bariatric Procedure Package with Split Leg

#BF852

Includes:

- Split Leg Positioner (BF351)
- (2) Multi-Axis Armboards RC (BF71239)
- Left Foot Rest (MTPA09)
- Right Foot Rest (MTPA10)
- Pair of Knee Restraint Straps (BF552)
- Extra-long Restraint Strap 126" (BF451)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Special Note(s): Side Rail Locks (BF133) and Clark Sockets (BF083) are required (sold separately). Bariatric Table Extenders also recommended.

Gastric Bypass Package

#BF851

Includes:

- (2) Multi-Axis Armboards RC (BF71239)
- Bariatric Foot Extension (BF542)
- Extra-long Restraint Strap 126" (BF451)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Special Note(s): Side Rail Locks (BF133) and Clark Sockets (BF083) are required (sold separately). Bariatric Table Extenders are also required.

Cardiovascular Accessory Package

#BF412

Includes:

Fem/Pop Board (BF287)

Perineal Cutout Filler Piece (BF205)

Table Compatibility: 3080, 3080-R, 3085 SP

Special Note(s): See 36" Carbon Fiber Fluoro Extension (BF642) for alternate product option.

Split Leg Package

#BF537

Includes:

Split Leg Positioner (BF351)

(2) Clark Sockets (BF083)

Left Foot Rest (MTPA09)

Right Foot Rest (MTPA10)

Pair of Knee Restraint Straps (BF552)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Trendelenburg Accessory Package

#BF607

Includes:

TrenStop (BF728)

(2) Rail Clamps (BF462)

(2) Multi-Axis Armboards RC with Pads (BF71239)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Uro/Cysto Accessory Package

#BF574

Includes:

Lightweight Transfer Board (BF588)

Power-Lift Stirrup Set (BF453) with Rail Clamps (BF462)

Power-Lift Stirrup Storage Cart (BF46)

Fluid Collector System (BF457)

Box of 10 Non-Sterile Fluid Collector Bags (BF459)

Foot Control (BF586)

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Uro/Endo Accessory Package with Fluid Collection System

#BF576

Includes:

9 Inch Urology Extension (BF578)

Patient Transfer Board with 2" Pad (BF204)

Power-lift Stirrup Set (BF453) with Rail Clamps (BF462)

Power-lift Stirrup Storage Cart (BF46)

Fluid Collector System (BF457)

Box of 10 Non-Sterile Fluid Collector Bags (BF459)

Foot Control (BF586)

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Uro/Endo Accessory Package with Drain Pan

#BF609

Includes:

- 9 Inch Urology Extension (BF578)
- Patient Transfer Board with 2" Pad (BF204)
- Power-lift Stirrup Set (BF453) with Siderail Clamps
- Power-lift Stirrup Storage Cart (BF46)
- Drain Pan (BF186)
- Foot Control (BF586)

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

3085 Uro/Endo Accessory Package

#BF402

Includes:

- Foot Control (BF187)
- Drain Pan (BF186)
- Patient Transfer Board with 2-inch Pad (BF204)
- IA Extender with 2-inch Pad and X-ray Top (BF188)
- Cable Wrap (BF316)

Table Compatibility: 3080, 3080-R, 3085 SP

Bariatric Urology Package

#BF853

Includes:

- Extra-long Restraint Strap 126" (BF451)
- (2) Multi-Axis Armboard RC (BFR71239)
- Bariatric Power-lift Stirrup (BF454)
- Power-lift Stirrup Storage Cart (BF46)
- (2) Bariatric Reverse Trendelenburg Strap (BF563)

Special Note(s): Fits Bariatric Table Extenders (sold separately) applied to 5085 SRT; 5085; 4085; Cmax; 3085 SP; 3080-R and 3080 tables.

Shoulder Chair Package

#BF571

Includes:

- Shoulder Chair (BF594)
- Set of 10 Head Restraint Straps (BF039)
- Multi-Posture Armboard with 2-inch Pad (BF432)
- Pair of Clark Sockets (BF083)
- Standard Length Restraint Strap 66-1/2" (BF45)

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

3085 Shoulder Table Accessory Package

#BF382

Includes:

- Shoulder Table (BF437)
- (2) Lateral Brace Assemblies (BF602)
- Multi-Posture Armboard with 2-inch Pad (BF432)
- Standard Length Restraint Strap 66-1/2" (BF45)

Table Compatibility: 3080, 3080-R, 3085 SP

Orthopedic Extension Accessory Package

#BF806

Includes:

- Accessory Cart (BF24)
- Adult Wedge-shaped Sacral Rest with Pad (BF00066)
- Adult Perineal Post with Pad (BF00065)
- (2) L-shaped Traction Extensions (BF0003)
- Anesthesia Armboard with 2-inch Pad (BF803)
- (2) TenZor™ Carbon Fiber Traction Unit Assemblies with Padded Boots (BF72)
- Universal Legholder (BF266)
- (2) Traction Accessory Clamps (BF0002)
- Multi-Posture Armboard with 2-inch Pad (BF432)
- Patient Transfer Board with 4-inch Pad (BF214)
- Clark Socket (BF083)
- Ortho Extension Setup Guide (P150832316)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Enables hip pinning and supine femoral nailing procedures.

OrthoVision & Orthopedic Extension Supplemental Accessory Package

#BF808

Includes

- Tibia/Lateral Counteraction Support (BF00062)
- Hip Rest with Pad (BF00072)
- Surgical Drape Support with Clark Socket and Traction Accessory Clamp (BF00018)
- Straight Traction Extension (BF268)
- Skeletal Traction Bow/Block (BF724)
- Skeletal Traction Block (BF725)
- Traction Accessory Clamp (BF0002)

Table Compatibility: Orthopedic Extension, OrthoVision

Usage/Procedures: This package extends procedure capability of the Orthopedic Extension and OrthoVision table to include lateral patient positioning and skeletal traction.

OrthoVision & Orthopedic Extension Pediatric Accessory Package

#BF809

Includes:

- Child Wedge-shaped Sacral Rest with Pad (BF00048)
- Child Perineal Post without Pad (BF41)
- Child Foot Traction Boots - Pair (BF539)
- Child Hip Spica Assembly (BF00038)

Table Compatibility: Orthopedic Extension, OrthoVision

SurgiGraphic 1027 & 6000 Standard Accessory Package

#TA005

Includes:

- Patient Safety Straps - Pair (TA00001)
- Standard Siderail Adaptors - Set of 3 (TA00005)
- Universal Ether Screen with Clark Socket (BF074)
- Cable Wrap (BF316)
- Hand Control Bracket (BF318)
- Anesthesia Armboard with 1-inch Pad (BF448)

Table Compatibility: SurgiGraphic 1027, SurgiGraphic 6000CE

SurgiGraphic 6000 Vascular Accessory Package

#TA00501

Includes:

- Catheter Extension (TA00002)
- Radial AV Fistula Board including Pad and Pair of Siderail Locks (TA00003)

Table Compatibility: SurgiGraphic 6000CE

SurgiGraphic 1027 Vascular Accessory Package

#TA00502

Includes:

- Radial AV Fistula Board including pad and pair of Siderail Locks (TA00003)

Table Compatibility: SurgiGraphic 1027

Standard Accessory Package

#BF807

Includes:

- TLT Tabletop Pad (BF323)
- Wedge-shaped Sacral Rest with Pad - Adult (BF00066)
- Perineal Post with Pad - Adult (BF00065)
- (2) L-shaped Traction Extensions (BF0003)
- Anesthesia Armboard with 2-inch Pad (BF803)
- (2) TenZor Carbon Fiber Traction Unit Assembly with Boots (BF72)
- Universal Legholder (BF266)
- (2) Traction Accessory Clamp (BF0002)
- Multi-Posture Armboard with 2-inch Pad (BF432)
- Patient Transfer Board with 2-inch Pad (BF213)
- Clark Socket (BF083)

Table Compatibility: OrthoVision

Usage/Procedures: For Orthopedic procedures using OrthoVision table

Anesthesia Armboard without Pad

#BF431

The radiolucent carbon fiber armboard is easily attached and removed from surgical table siderail and is easy to clean. Components are enclosed within molded body, protecting against denting, breakage, and damage when dropped. Includes a 30-inch velcro arm restraint strap.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Replacement Pad(s): BF345 (1" TLT); BF326 (2" TLT); BF327 (3" TLT); BF394 (4" TLT); BF395 (5" HLT)

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).

Anesthesia Armboard with 1" Pad

#BF448

The radiolucent carbon fiber armboard with 1-inch Pad is easily attached and removed from surgical table siderail and is easy to clean. Components are enclosed within molded body, protecting against denting, breakage, and damage when dropped. Includes a 30-inch velcro arm restraint strap.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Replacement Pad(s): BF345

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



Anesthesia Armboard with 2" Pad

#BF803

The radiolucent carbon fiber armboard with 2-inch Pad is easily attached and removed from surgical table siderail and is easy to clean. Components are enclosed within molded body, protecting against denting, breakage, and damage when dropped. Includes a 30-inch velcro arm restraint strap.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Replacement Pad(s): BF326

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



Anesthesia Armboard with 3" Pad

#BF804

The radiolucent carbon fiber armboard with 3-inch Pad is easily attached and removed from surgical table siderail and is easy to clean. Components are enclosed within molded body, protecting against denting, breakage, and damage when dropped. Includes a 30-inch velcro arm restraint strap.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Replacement Pad(s): BF327

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



IV Armboard without Pad

#BF141

Supports arm at a comfortable angle. Locks by gravity to table siderail. Adjusts 180 degrees horizontally. Supplied with 15-inch velcro arm restraint strap.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Replacement Pad(s): BF345 (1"); BF326 (2"); BF327 (3")

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).

IV Armboard with 2" Pad

#BF363

Supports arm at a comfortable angle. Locks by gravity to table siderail. Adjusts 180 degrees horizontally. Supplied with 15-inch velcro arm restraint strap.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Replacement Pad(s): BF326

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



IV Armboard with 3" Pad

#BF364

Supports arm at a comfortable angle. Locks by gravity to table siderail. Adjusts 180 degrees horizontally. Supplied with 15-inch velcro arm restraint strap.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Replacement Pad(s): BF327

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



Multi-axis Armboard, RC


#BF71239

The Multi-axis Armboard allows the arm to be postured to meet the anatomical requirements of each patient. In Trendelenburg / Reverse Trendelenburg postures, the armboard adjusts so the arm can be positioned level to the floor. With the squeeze of a single handle, the armboard permits flexion, extension, abduction / adduction and rotation of the arm. This armboard is also height adjustable. It also features a rotating clamp that allows attachment to the back section of the table to accommodate the sitting position. In the seated posture, the patient's flexed arm can be adjusted to an ergonomically correct position. A contoured pad is included.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Replacement Pad(s): BF00029

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



Multi-posture Armboard with 2-inch Pad

#BF432

The Multi-posture Armboard with 2-inch Pad attaches to table siderail with a Clark Socket and is mounted on a ball joint for positioning flexibility. Provides upper arm support during lateral positioning, as well as general arm support during supine, sitting, and "park-bench" positioning; especially useful during prone positioning for preventing hyperextension.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Replacement Pad(s): BF81

Special Note(s): Clark Sockets (BF083) or Multi-Accessory Clamps (BF43) are required (sold separately). Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).

Krause Arm Support

#BF1

An arm length stockinette supports the arm without noticeable pressure. The support is appropriate for unilateral, prone, supine, and seated positioning of patient. The frame is adjustable for arm length with a universal ball joint and side-rail mounting allowing complete positioning flexibility. A Clark Socket is furnished with each support.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



Pediatric Arm and Hand Table

#BF434

Small, lightweight carbon fiber arm and hand table perfect for pediatric procedures and arterial line access. Installation and transportation easy due to light weight and size.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Required Accessories: BF083

Replacement Pad(s): BF814

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



Adult Arm and Hand Table

#BF435

Lightweight carbon fiber arm and hand table provides unobstructed C-arm and surgeon access for major and minor procedures on arms and hands. Size accommodates both major and minor vascular or orthopedic procedures.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Required Accessories: BF083

Replacement Pad(s): BF812

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Sidrail Adaptors (TA00012).



Micro-Surgical Arm and Hand Table

#BF436

Hourglass-shaped table provides stability and surgeon access with ample working space for when performing vascular and orthopedic procedures on hands or arms. Table top is constructed of carbon fiber material allowing for imaging of the extremity and lightweight placement on the surgical table. Elbow rests reduce surgeon fatigue during lengthy procedures. Support leg for micro-surgical stability. Distal end rail enables arm and hand traction.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Required Accessories: BF083

Replacement Pad(s): BF811

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Sidrail Adaptors (TA00012).



AQUAGEL Arm Protector, Short

#BF475

This AQUAGEL Pressure Management Positioner is made with dual density construction that most closely resembles human tissue. AQUAGEL arm protectors are ergonomically designed to conform to each patient's arm. The 3-dimensional construction manages and re-distributes pressure as well as provides protection for skin tissue and nerves.

Specifications: 13 x 5-1/2 x 1"



AQUAGEL Arm Protector, Long

#BF476

This AQUAGEL Pressure Management Positioner is made with dual density construction that most closely resembles human tissue. AQUAGEL arm protectors are ergonomically designed to conform to each patient's arm. The 3-dimensional construction manages and re-distributes pressure as well as provides protection for skin tissue and nerves.

Specifications: 18 x 7 x 1-1/2"



OASIS Armboard, Short

#BF426

Flat Armboard Pads are available in two sizes to accommodate the requirements of any armboard. Flat AQUAGEL armboard pads contain a silicone gel that most closely replicates human tissue to protect against skin shear and nerve injury.

Specifications: 15 x 4 x 1/2"



OASIS Armboard, Long

#BF427

Flat Armboard Pads are available in two sizes to accommodate the requirements of any armboard. Flat AQUAGEL armboard pads contain a silicone gel that most closely replicates human tissue to protect against skin shear and nerve injury.

Specifications: 23-1/2 x 5 x 1/2"



Patient Safety Rail with Cover


#BF092

This padded safety rail is used to provide safe patient transport throughout the perioperative area while in transfer on the 5085 SRT table. Sold individually with integral clamp.

Table Compatibility: 5085 SRT

Usage/Procedures: Patient Transport

Replacement Pad(s): BF094



Restraint Strap; Standard Length 66-1/2"

#BF45

Designed with quick-release buckles on both ends for easy, positive adjustment. Latex free extra-wide webbing is soil-resistant and washable.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax



Restraint Strap; Extra-long 126"

#BF451

Designed with quick-release buckles on both ends for easy, positive adjustment. Latex free extra-wide webbing is soil-resistant and washable.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax



Standard Reverse Trendelenburg Strap (Pair)

#BF455

The broad pads on these straps safely encircle the inner thigh without binding or pinching, making them ideal for Nissen Fundoplication or other extreme Reverse Trendelenburg procedures. On this standard model, the safety belt buckle ensures a snug fit for patients up to 500 lbs. These straps attach easily and securely to the table with integral rail clamps.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Reverse Trendelenburg

Optional Accessories: BF453



Bariatric Reverse Trendelenburg Strap (Pair)

#BF563

The broad pads on these straps safely encircle the inner thigh without binding or pinching, making them ideal for Nissen Fundoplication or other extreme Reverse Trendelenburg procedures. On this bariatric model, the safety belt buckle ensures a snug fit for patients up to 1000 lbs. These straps attach easily and securely to the table with integral rail clamps.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Reverse Trendelenburg

Optional Accessories: BF454



Knee Restraint Strap (Pair)

#BF552

The broad pads on these straps safely prevent the knee from flexing during surgery and feature a full-length Velcro closure for a snug fit. Used in conjunction with the Split Leg Positioner.

Usage/Procedures: Laparoscopic; Colo-rectal; Vascular; Extreme Reverse Trendelenburg

Special Note(s): Compatible with tables using the Split Leg Positioner (BF351) including 5085 SRT, 5085, 4085, Cmax, 3085 SP, 3080-R, and 3080 tables.



Lateral Safety Strap

#BF562

This padded safety strap allows you to secure the operative side of the patient safely and comfortably. Attaches to the side rail with a post that can be inserted into a Clark Socket or Multi-accessory Clamp.

Usage/Procedures: Open Shoulder; Closed Shoulder

Special Note(s): Clark Sockets (BF083) or Multi-accessory Clamps (BF43) are required (sold separately). Compatible with 5085 SRT, 5085, 4085, Cmax, 3085 SP, 3080-R, and 3080 tables using the Shoulder Chair (BF594) or Shoulder Table (BF437).




Wristlets with Straps and Brackets (Pair)

#BF111

These support wrists and forearms while patient is in supine operative position. Padded leather straps secure wrists quickly and comfortably. Supporting brackets, with locks, mount to table siderails for sliding adjustment. Straps, brackets, and siderail locks are included.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax



Patient Safety Straps (Pair)

#TA0001

Radiolucent strap wraps around patient and tabletop. Fits all size patients with Velcro securement.

Table Compatibility: SurgiGraphic 1027, SurgiGraphic 6000

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.



TrenStop

#BF728

The TrenStop combines self-adjusting shoulder support wings with a unique cervical notch pad. Shoulder support wings rotate freely (and locked securely) to meet each patient's unique anatomical requirements. Extra wide pads distribute weight across a large area for patient safety. The cervical support pad is placed under patient's cervical spine curvature and attached to the horizontal support bar. This pad prevents funneling effect, reducing the risk of anatomical distortion.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity 450 lbs. (204 kg).

Usage/Procedures: Trendelenburg up to 45 degrees and Robotics procedures.

Replacement Pad(s): BF72876 (Shoulder); BF72877 (Cervical)

Special Note(s): Two Rail Clamps (BF462) are required (sold separately).



Perineal Cutout Filler Piece

#BF205

This filler piece covers perineal opening in seat of table. Provides additional patient support.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Special Note(s): Also covers perineal opening for 3080/3085 accompanying seat section X-ray top. Compatible with most Cmax versions. Consult sales rep for additional information.



Lateral Brace Assemblies (Pair)

#BF602

Each brace assembly attaches to table using a Siderail Lock and consists of a flat-bar stock attached to an angular control socket and ball joint. Each brace can be positioned horizontally or vertically and can be angled 7-inches internally. Supplied with a 3x7" pad.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Lateral; Orthopedic

Required Accessories: BF133



Solid Lateral Braces (Pair)

#BF331

Anatomically curved for comfortable lateral body support. Mount easily on kidney elevator (kidney bridge).

Table Compatibility: 3080, 3080-R, 3085 SP

Usage/Procedures: Chest; Lateral



GEL Supine Package, Large

#BF486

Includes:

Supine Headrest, Large (BF466)

(2) Arm Protectors, Long (BF476)

Sacral Protector (BF474)

(2) Leg Protectors (BF478)

GEL Supine Package, Medium

#BF487

Includes:

- Supine Headrest, Medium (BF467)
- (2) Arm Protectors, Short (BF475)
- Sacral Protector (BF474)
- (2) Leg Protectors (BF478)

GEL Lateral Package, Large

#BF488

Includes:

- Head Protector (BF463)
- Lateral Positioner (BF484)
- (3) Arm Protector, Long (BF476)

GEL Lateral Package, Medium

#BF489

Includes:

- Head Protector (BF463)
- Lateral Positioner (BF484)
- (3) Arm Protector, Short (BFR475)

AQUAGEL Sacral Protector, Small

#BF474

This AQUAGEL Pressure Management Positioner is made with dual density construction that most closely resembles human tissue. The Sacral Protector is an essential positioner for all supine patients, especially those in the Lithotomy position. This 3-dimensional positioner is contoured specifically to re-distribute pressure away from the sacrum. Recommended companion products are either the 20" Profile Table Overlay (BF471) or the 45" Profile Table Overlay (BF472).

Specifications: 15-3/4 x 15 x 1"




AQUAGEL Patient Positioner

#BF481

This AQUAGEL Pressure Management Positioner is made with dual density construction that most closely resembles human tissue. The Patient Positioner is a versatile volumetric product that may be used to posture and protect patients while in supine, lateral, and prone positions.

Specifications: 16 x 5-3/4 x 2-1/2"




ELITE Flat Bottom Chest Roll, Medium

#BF482

Used for prone and lateral postures and may be used in conjunction with the Prone Positioner (BF485) and Prone Headrest (BF464 large/BF465 medium) as shown.

Specifications: 16 x 6 x 4"



ELITE Flat Bottom Chest Roll, Large

#BF483

Used for prone and lateral postures and may be used in conjunction with the Prone Positioner (BF485) and Prone Headrest (BF464 large/BF465 medium) as shown.

Specifications: 20 x 6 x 4"



ELITE Lateral Positioner

#BF484

Use this posturing device for all lateral positions. The Lateral Positioner permits a safe, quick and easy transition from supine to lateral. The patient's lower shoulder fits into the integral channel to provide exceptional patient comfort and nerve protection. The built-in headrest and axillary support assists in proper spinal alignment and surgical site access.

Specifications: 28 x 20 x 6"



ELITE Prone Positioner

#BF485

The unique shape of the Prone Positioner is ideal for spinal procedures. Use this positioner in conjunction with the Prone Headrest (BF464 large or BF465 medium), Flat-bottom Chest Rolls (BF483 large or BF482 medium) and the Patient Positioner (BF481) for proper patient alignment (as shown in photo).

Specifications: 24 x 4-1/2 x 20-1/2"




Clark Sockets (Pair)

#BF083

Ideal for mounting accessories that have round posts anywhere along table siderail. These sockets provide angular control and positive locking for any accessory in a plane parallel to siderail, and are easily removable.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Special Note(s): Also fits SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



Multi-Accessory Clamps (Pair)

#BF43

Allows independent release and positioning of horizontal, vertical, and rotational functions, making it easy for one person to adjust the position of accessories. Multi-Accessory Clamps can be applied anywhere along table siderail, and are designed for easy cleaning with no exposed teeth. These clamps are used for securing accessories with either round, square or rectangular posts.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Special Note(s): Also fits SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012). Maximum patient weight capacity is equivalent to those STERIS table accessories where clamp or lock is specified to be used to support the accessory.



Siderail Locks (Pair)

#BF133

Notched siderail locks attach to table siderails for securing flat-bar stock accessories. Useful with Foot Extensions (BF2 or BF542) either horizontally or vertically.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



Rail Clamp (Each)

#BF462

These Rail Clamps come standard with both models (standard and bariatric) of our Power-lift Stirrups, but if needed, individual extra (or replacement) Rail Clamps can be ordered separately. They can be used to support any flat-bar stock accessory perpendicular to the table surface.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision


Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).

Standard Siderail Adaptor (Each)

#TA00012

The Standard Siderail Adaptors provide attachment of accessories requiring standard siderails. Moveable and removable anywhere along side of table.

Table Compatibility: SurgiGraphic 1027, SurgiGraphic 6000, 36" Carbon Fiber Fluoro Extension



Standard Siderail Adaptor (Set of 3)

#TA00005

The Standard Siderail Adaptors provide attachment of accessories requiring standard siderails. Moveable and removable anywhere along side of table.

Table Compatibility: SurgiGraphic 1027, SurgiGraphic 6000, 36" Carbon Fiber Fluoro Extension



Perineal Post Pad - Disposable

#NO001

Relieves pressure from the perineal area. Fits standard perineal posts and can be trimmed to fit child perineal posts. Made of latex-free disposable foam. Case of 10.

Table Compatibility: Orthopedic Extension, OrthoVision

Usage/Procedures: Orthopedic procedures in supine position



Drain Pan Tubing

#R003500835

This replacement tubing can be ordered by the foot. Indicate the length desired, otherwise, the tubing will come in foot long segments.

Usage/Procedures: Urologic

Special Note(s): Compatible with Drain Pan (BF186) (sold separately).



Fluid Collector Bags; Sterile Box of 10


#BF458

These sterile disposable bags slide over the flexible drape support. Each bag has a tissue retention filter as well as a flexible, expandable drain hose that can be connected to suction or used with a kick bucket. These bags fit both the Standard and Bariatric drape supports.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Urologic; Cystoscopic; Hysteroscopic

Special Note(s): Fluid Collection System: Standard (BF457) or Bariatric (BF56) is required (sold separately).



Fluid Collector Bags; Non-Sterile Box of 10

#BF459

These non-sterile disposable bags slide over the flexible drape support. Each bag has a tissue retention filter as well as a flexible, expandable drain hose that can be connected to suction or used with a kick bucket. These bags fit both the Standard and Bariatric drape supports.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Urologic; Cystoscopic; Hysteroscopic

Special Note(s): Fluid Collection System: Standard (BF457) or Bariatric (BF56) is required (sold separately).



Narrow Headrest

#BF592

This headrest is ideal for neurological procedures. The Wrist Rest (BF593) attaches to the Narrow Headrest to support the surgeon's hands during long procedures.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Neurological; Ophthalmic; ENT; Plastic Surgery

Optional Accessories: BF593

Replacement Pad(s): BF491

Special Note(s): The Eye-ENT-Neuro Headrest Adaptor is required (sold separately). Order BF185 for 3085 tables, BF173 for 5085, 4085, and Cmax tables, or BF309 for 5085SRT.



Cone-Shaped Headrest

#BF215

This headrest fits directly into a table's standard headrest area and provides easier access to head area for eye, ENT, plastic surgery, and neurosurgical applications. Includes a 2-inch pad and a 4-inch pad for use with X-ray tops.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Eye; ENT; Neurosurgical; Plastic Surgery

Replacement Pad(s): BF83 (2"); BF831 (4")



Eye-ENT-Neuro Headrest Adaptor - 3080, 3080-R, 3085

#BF185

This specialized adaptor accommodates all standard Eye-ENT-Neuro attachments. Conveniently inserts into either end of the 3085 SP surgical table.

Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is 500 pounds (227 kg) for 3085 and 300 (136 kg) for 3080 and 3080-R.

Usage/Procedures: Neurosurgical; Ophthalmic; ENT; Plastic Surgery



Eye-ENT-Neuro Headrest Adaptor - 5085, 4085, Cmax

#BF173

This specialized adaptor accommodates all standard Eye-ENT-Neuro attachments. Conveniently inserts into either end of the 5085, 4085, or Cmax surgical tables.

Table Compatibility: 4085, 5085, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Neurosurgical; Ophthalmic; ENT; Plastic Surgery



Eye-ENT-Neuro Headrest Adaptor - 5085 SRT

#BF309

This specialized adaptor accommodates all standard Eye-ENT-Neuro attachments. Conveniently inserts into either end of the 5085 SRT surgical table.

Table Compatibility: 5085 SRT

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Neurosurgical; Ophthalmic; ENT; Plastic Surgery



Frontal Sitting Attachment

#BF01


This accessory mounts neurosurgical headrests to head or foot end of 3085 SP, 3080-R or 3080 tables or to head end of 5085 SRT, 5085, 4085 or Cmax tables. Ball joints allow full angular positioning of headrest.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 300 lbs. (136 kg).

Usage/Procedures: For support of patient's head from behind when patient is in a seated position for temporal, parietal, and frontal surgical procedures.

Special Note(s): Eye-ENT-Neuro Headrest Adaptor and either Multi-Poise Headrest (BF038) or Horseshoe Headrest (BF036) are required (sold separately).



Horizontal Attachment #1800

#BF008

Mounts neurosurgical headrests to head or foot end of 3085 SP, 3080-R and 3080 tables or to head end of 5085 SRT, 5085, 4085 or Cmax tables. Ball joints allow full angular positioning of headrest.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 300 lbs. (136 kg).

Usage/Procedures: For supporting patient's head in supine, prone, or lateral positions.

Special Note(s): Eye-ENT-Neuro Headrest Adaptor and either Horizontal Attachment #1800 (BF008) or Frontal Sitting Attachment (BF01) are required (sold separately).



Horseshoe Headrest

#BF036

For surgical procedures involving nose, paranasal, sinuses, trachea, and larynx approaches, as well as some neurological procedures.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 300 lbs. (136 kg).

Usage/Procedures: Cradles patient's head for eye, ENT, neuro, and plastic surgery.

Replacement Pad(s): BF492

Special Note(s): Eye-ENT-Neuro Headrest Adaptor and either Horizontal Attachment (BF008) or Frontal Sitting Attachment (BF01) are required (sold separately).



Posterior Sitting Attachment #1801

#BF011

Mounts neurosurgical headrests to table siderails. Ball joints allow lateral and vertical adjustment.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 300 lbs. (136 kg).

Usage/Procedures: For support of patient's head from the front when patient is in a seated position for suboccipital and parietal surgical procedures.

Special Note(s): Multi-Poise Headrest (BF038) and either Clark Sockets (BF083) or Multi-Accessory Clamps (BF43) are required (sold separately).



Multi-Poise Headrest

#BF038

This headrest permits posturing flexibility for neurosurgical techniques. Three sterilizable spring-laden pins assure maximum rigidity of patient's head without discomfort.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 300 lbs. (136 kg).

Usage/Procedures: Used with patient in supine, prone, lateral (park bench), or sitting position. Neurosurgical Craniotomies

Special Note(s): Eye-ENT-Neuro Headrest Adaptor and either Horizontal Attachment #1800 (BF008), Frontal Sitting Attachment #1802 (BF01), or Posterior Sitting Attachment #1801 (BF011) are required (sold separately).



AQUAGEL Head Protector

#BF463

This AQUAGEL Pressure Management Positioner is made with dual density construction that most closely resembles human tissue. The comfortable 3-dimensional Contoured Head Protector provides stability and pressure management for supine and lateral positioning.

Specifications: 10 x 10 x 2"



ELITE Prone Headrest, Large

#BF464

The unique molded design of the lightweight AQUAGEL Prone Headrest permits easy endotracheal tube access as well as excellent head support and stability. This positioner also protects the forehead, eyes and chin. Companion products include Flat-bottom Chest Rolls (BF483 large or BF482 medium) and the Prone Positioner (BF485).

Specifications: 11 x 9 x 6"



ELITE Prone Headrest, Medium

#BF465

The unique molded design of the lightweight AQUAGEL Prone Headrest permits easy endotracheal tube access as well as excellent head support and stability. This positioner also protects the forehead, eyes and chin. Companion products include Flat-bottom Chest Rolls (BF483 large or BF482 medium) and the Prone Positioner (BF485).

Specifications: 9 x 7-1/2 x 5"



ELITE Supine Headrest, Large

#BF466

The ergonomically molded AQUALGEL Supine Headrest is ideal for positioning the patient during eye surgery as well as for other procedures requiring the supine position.

Specifications: 10-3/4 x 9-1/2 x 3-1/2"



ELITE Supine Headrest, Medium

#BF467

The ergonomically molded AQUALGEL Supine Headrest is ideal for positioning the patient during eye surgery as well as for other procedures requiring the supine position.

Specifications: 10½ x 9 x 2½"



Wrist Rest for Narrow Headrest

#BF593

This Wrist Rest attaches to the Narrow Headrest (BF592) to support the surgeon's hands during long procedures.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Neurological; Ophthalmic; ENT; Plastic Surgery

Required Accessories: BF592



Featherweight Leg Section - 5085, 4085, Cmax


#BF009

This carbon fiber foot section offers a lightweight design engineered to allow for easier removal and installation from 5085, 4085, or Cmax table. Great replacement option for clinical procedures requiring frequent removal and installation of leg section since lightweight design supports easier handling with less strain.

Table Compatibility: 4085, 5085, Cmax

Replacement Pad(s): BF83300339

Special Note(s): If Cmax table uses BF584 X-ray tops, BF744 will need to be ordered for the Featherweight Leg Section X-ray top.




Featherweight Leg Section - 5085 SRT

#BF012

This carbon fiber foot section offers a lightweight design engineered to allow for easier removal and installation from the 5085 SRT table. Great replacement option for clinical procedures requiring frequent removal and installation of leg section since lightweight design supports easier handling with less strain.

Table Compatibility: 5085 SRT



Power-lift Stirrups (Pair)

#BF453

The integral power assist mechanism in these stirrups provides excellent surgeon access and virtually effortless patient positioning of patients up to 500 lbs. The unique grip-and-twist control handles allow staff to adjust the angle of Lithotomy and degree of ab/adduction during set up and during the procedure without lifting the sterile drape. Includes Rail Clamps and pads.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Endourology; Gynecology; Perineal; Other mid-line access procedures

Optional Accessories: BF455, BF46

Replacement Pad(s): BF85000381 (Left); BF85000382 (Right)



Bariatric Power-lift Stirrups (Pair)

#BF454

These stirrups feature the same grip-and-twist control handles and use the same power assist technology as the standard Power-lift Stirrups (BF453), but can accommodate patients up to 720 lbs., depending on the table used. The extra large boots and pads provide comfort and security. Includes rail clamps and pads.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 720 lbs. (377 kg) for 5085 SRT, 5085, 4085, Cmax and 3085 SP; 600 lbs for 3080-R and 400 lbs. (181 kg) for 3080 tables.

Usage/Procedures: Endourology; Gynecology; Perineal; Bariatric; Colo-rectal; Other mid-line access procedures

Optional Accessories: BF46, BF563

Replacement Pad(s): BF191 (Left); BF192 (Right)



Split Leg Positioner

#BF351

The Split Leg Positioner includes both left and right leg supports and provides a safe system for positioning patients for laparoscopic procedures as well as lower extremity orthopedic and vascular procedures. It is radiolucent allowing both the surgeon and imaging equipment excellent access to the operative site. Patient's legs can be adducted, abducted, raised or lowered for low Lithotomy procedures.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg) for 5085 SRT, 5085, 4085, Cmax, 3085 SP and 3080-R, and 400 lbs. (181 kg) for 3080 tables.

Usage/Procedures: Lower extremity procedures; Vascular; Orthopedic; Bariatric; Mid-line access to abdomen

Optional Accessories: BF552, MTPA09, MTPA10

Replacement Pad(s): BF351000337 (Left); BF351000338 (Right)



Foot Rest; Left

#MTPA09

The foot rests assist in ergonomic positioning of the patient's feet. They are especially useful when using the Split Leg Positioner. The foot rests can be attached using Clark Sockets or Multi-Accessory Clamps (sold separately).

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Lower extremity procedures; Vascular Orthopedic; Mid-line access

Required Accessories: BF351

Special Note(s): Clark Sockets (BF083) or Multi-Accessory Clamps (BF43) are required (sold separately).



Foot Rest; Right

#MTPA10

The foot rests assist in ergonomic positioning of the patient's feet. They are especially useful when using the Split Leg Positioner. The foot rests can be attached using Clark Sockets or Multi-Accessory Clamps (sold separately).

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Lower extremity procedures; Vascular Orthopedic; Mid-line access

Required Accessories: BF351

Special Note(s): Clark Sockets (BF083) or Multi-Accessory Clamps (BF43) are required (sold separately).



Lithotomy Legholders - Adult (Pair)

#BF054


The Lithotomy Legholders provide cantilevered suspension-strap support for legs. Height and angulation are adjustable for full positioning flexibility.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 300 lbs. (136 kg).

Usage/Procedures: OB/GYN; Urologic

Special Note(s): Clark Sockets (BF083) or Multi-Accessory Clamps (BF43) are required (sold separately).



Multi-Position Legholders (Pair)

#BF433

The Multi-position Legholders are designed to accommodate varying positioning requirements. Boot design prevents uncontrolled abduction, even when unlocked. Telescoping knee joints provide added flexibility. Mount to table side rails. Includes one set of disposable foam insert pads.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax
Usage/Procedures: Gynecologic; Urologic; Laparoscopic Nissen Fundoplication
Special Note(s): Clark Sockets (BF083) or Multi-Accessory Clamps (BF43) are required (sold separately).



Universal Legholder (Each)


#BF266

Each Universal Legholder is contoured and padded for comfortable leg support, and is mounted on a ball joint for positioning flexibility. Includes strap.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Usage/Procedures: Orthopedic fracture procedures to support patient's unaffected leg

Special Note(s): Clark Sockets (BF083) or Multi-Accessory Clamps (BF43) for Cmax, 3080-R, 3085SP are required (sold separately). Traction-Accessory Clamp (BF0002) for the OrthoVision table and Orthopedic Extension is required (sold separately).



Bariatric Foot Extension

#BF542


The Bariatric Foot Extension is 28-inches wide to accommodate wider tables with Bariatric Extenders installed. It adds 8-inches of table length for taller patients and swings 90 degrees upward for foot support when in reverse Trendelenburg.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax
Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Required Accessories: BF133

Replacement Pad(s): BF493

Special Note(s): Bariatric Table Extensions are required (sold separately).



Foot Extension with Siderail Locks

#BF2

This Foot Extension extends table 8-inches to accommodate taller patients and swings 90 degrees upward for foot support when patient is in reverse Trendelenburg. This full-width (20-inches) extension locks to the side rails of the table foot section with Siderail Locks (pair included).

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax
Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Replacement Pad(s): BF816



AQUAGEL Heel Protector

#BF477

These AQUAGEL Pressure Management Positioners are made with dual density construction that most closely resembles human tissue. These 3-dimensional Contoured Heel Protectors provide heel protection and positioning stability.

Specifications: 7-1/2 x 7 x 1"



AQUAGEL Leg Protector

#BF478

This AQUAGEL Pressure Management Positioner is made with dual density construction that most closely resembles human tissue. The 3-dimensional Contoured Leg Protector provides ultimate pressure relief for the heel as it stabilizes the lower leg and protects the knee from hyperextension.

Specifications: 20 x 7 x 1-1/2"



OASIS Crutch Stirrup

#BF48

The Crutch Stirrup pads were designed for use with various types of legholders during Lithotomy. Lightweight and flexible, these flat gel pads contain a silicone gel that most closely replicates human tissue to protect against skin shear and nerve injury. They come in pairs and may be used with patients of all sizes.

Specifications: 15-1/2 x 9-1/2 x 1/2"



GEL Supine Package, Large

#BF486

Includes:

- Supine Headrest, Large (BF466)
- (2) Arm Protectors, Long (BF476)
- Sacral Protector (BF474)
- (2) Leg Protectors (BF478)

GEL Supine Package, Medium

#BF487

Includes:

- Supine Headrest, Medium (BF467)
- (2) Arm Protectors, Short (BF475)
- Sacral Protector (BF474)
- (2) Leg Protectors (BF478)

GEL Lateral Package, Large

#BF488

Includes:

- Head Protector (BF463)
- Lateral Positioner (BF484)
- (3) Arm Protector, Long (BF476)

GEL Lateral Package, Medium

#BF489

Includes:

- Head Protector (BF463)
- Lateral Positioner (BF484)
- (3) Arm Protector, Short (BFR475)

AQUAGEL Head Protector

#BF463

This AQUAGEL Pressure Management Positioner is made with dual density construction that most closely resembles human tissue. The comfortable 3-dimensional Contoured Head Protector provides stability and pressure management for supine and lateral positioning.

Specifications: 10 x 10 x 2"

ELITE Prone Headrest, Large

#BF464

The unique molded design of the lightweight AQUAGEL Prone Headrest permits easy endotracheal tube access as well as excellent head support and stability. This positioner also protects the forehead, eyes and chin. Companion products include Flat-bottom Chest Rolls (BF483 large or BF482 medium) and the Prone Positioner (BF485).

Specifications: 11 x 9 x 6"



ELITE Prone Headrest, Medium

#BF465

The unique molded design of the lightweight AQUAGEL Prone Headrest permits easy endotracheal tube access as well as excellent head support and stability. This positioner also protects the forehead, eyes and chin. Companion products include Flat-bottom Chest Rolls (BF483 large or BF482 medium) and the Prone Positioner (BF485).

Specifications: 9 x 7-1/2 x 5"



ELITE Supine Headrest, Large

#BF466

The ergonomically molded AQUAGEL Supine Headrest is ideal for positioning the patient during eye surgery as well as for other procedures requiring the supine position.

Specifications: 10-3/4 x 9-1/2 x 3-1/2"



ELITE Supine Headrest, Medium

#BF467

The ergonomically molded AQUAGEL Supine Headrest is ideal for positioning the patient during eye surgery as well as for other procedures requiring the supine position.

Specifications: 10½ x 9 x 2½"



AQUAGEL Sacral Protector, Small

#BF474

This AQUAGEL Pressure Management Positioner is made with dual density construction that most closely resembles human tissue. The Sacral Protector is an essential positioner for all supine patients, especially those in the Lithotomy position. This 3-dimensional positioner is contoured specifically to re-distribute pressure away from the sacrum. Recommended companion products are either the 20" Profile Table Overlay (BF471) or the 45" Profile Table Overlay (BF472).

Specifications: 15-3/4 x 15 x 1"



AQUAGEL Patient Positioner

#BF481

This AQUAGEL Pressure Management Positioner is made with dual density construction that most closely resembles human tissue. The Patient Positioner is a versatile volumetric product that may be used to posture and protect patients while in supine, lateral, and prone positions.

Specifications: 16 x 5-3/4 x 2-1/2"




ELITE Flat Bottom Chest Roll, Medium

#BF482

Used for prone and lateral postures and may be used in conjunction with the Prone Positioner (BF485) and Prone Headrest (BF464 large/BF465 medium) as shown.

Specifications: 16 x 6 x 4"




ELITE Flat Bottom Chest Roll, Large

#BF483

Used for prone and lateral postures and may be used in conjunction with the Prone Positioner (BF485) and Prone Headrest (BF464 large/BF465 medium) as shown.

Specifications: 20 x 6 x 4"



ELITE Lateral Positioner

#BF484

Use this posturing device for all lateral positions. The Lateral Positioner permits a safe, quick and easy transition from supine to lateral. The patient's lower shoulder fits into the integral channel to provide exceptional patient comfort and nerve protection. The built-in headrest and axillary support assists in proper spinal alignment and surgical site access.

Specifications: 28 x 20 x 6"



ELITE Prone Positioner

#BF485

The unique shape of the Prone Positioner is ideal for spinal procedures. Use this positioner in conjunction with the Prone Headrest (BF464 large or BF465 medium), Flat-bottom Chest Rolls (BF483 large or BF482 medium) and the Patient Positioner (BF481) for proper patient alignment (as shown in photo).

Specifications: 24 x 4-1/2 x 20-1/2"



AQUAGEL Arm Protector, Short

#BF475

This AQUAGEL Pressure Management Positioner is made with dual density construction that most closely resembles human tissue. AQUAGEL arm protectors are ergonomically designed to conform to each patient's arm. The 3-dimensional construction manages and re-distributes pressure as well as provides protection for skin tissue and nerves.

Specifications: 13 x 5-1/2 x 1"



AQUAGEL Arm Protector, Long

#BF476

This AQUAGEL Pressure Management Positioner is made with dual density construction that most closely resembles human tissue. AQUAGEL arm protectors are ergonomically designed to conform to each patient's arm. The 3-dimensional construction manages and re-distributes pressure as well as provides protection for skin tissue and nerves.

Specifications: 18 x 7 x 1-1/2"



OASIS Armboard, Short

#BF426

Flat Armboard Pads are available in two sizes to accommodate the requirements of any armboard. Flat AQUAGEL armboard pads contain a silicone gel that most closely replicates human tissue to protect against skin shear and nerve injury.

Specifications: 15 x 4 x 1/2"



OASIS Armboard, Long

#BF427

Flat Armboard Pads are available in two sizes to accommodate the requirements of any armboard. Flat AQUAGEL armboard pads contain a silicone gel that most closely replicates human tissue to protect against skin shear and nerve injury.

Specifications: 23-1/2 x 5 x 1/2"



AQUAGEL Heel Protector

#BF477

These AQUAGEL Pressure Management Positioners are made with dual density construction that most closely resembles human tissue. These 3-dimensional Contoured Heel Protectors provide heel protection and positioning stability.

Specifications: 7-1/2 x 7 x 1"



AQUAGEL Leg Protector

#BF478

This AQUAGEL Pressure Management Positioner is made with dual density construction that most closely resembles human tissue. The 3-dimensional Contoured Leg Protector provides ultimate pressure relief for the heel as it stabilizes the lower leg and protects the knee from hyperextension.

Specifications: 20 x 7 x 1-1/2"



OASIS Crutch Stirrup

#BF48

The Crutch Stirrup pads were designed for use with various types of legholders during Lithotomy. Lightweight and flexible, these flat gel pads contain a silicone gel that most closely replicates human tissue to protect against skin shear and nerve injury. They come in pairs and may be used with patients of all sizes.

Specifications: 15-1/2 x 9-1/2 x 1/2"



Hand Control (Advanced) - 5085, 5085 SRT

#BF329

An easy-to-read LCD window illustrates table positioning data, making it perfect for clinicians behind the drape. Pre-programmed, backlit buttons allow quick, one-touch patient posturing, even in darkened ORs. Programmable and manual control features provide additional flexibility to articulate virtually every patient posture.

Table Compatibility: 5085, 5085 SRT



Hand Control (Universal) - 5085, 4085, Cmax

#BF319

The 4085's intuitive hand control instantly indicates table positioning data to caregivers with an easy-to-read digital display without lifting the sterile drape. A return-to-level button quickly and safely returns the patient to level mode. Its Auto Limit Sensor(TM) (ALS) allows you to operate the table without having to worry about unintended sectional collisions. This hand control is backwards compatible with Cmax tables.

Table Compatibility: 4085, 5085, Cmax



Hand Control - Cmax

#BF61

No longer available, see BF319



Hand Control - Amsco 3085 SP ACT-enabled

#BF374

This replacement hand control was designed to work exclusively with the Automated Control Technologies (ACT) digital integrated OR software.

Table Compatibility: 3085 SP



Hand Control - 3085 SP, 3080 RL/SP

#BF336

Lightweight, removable, fluid-tight, and shock-resistant, this remote control device includes automatic "return-to-level" positioning. Also provides easy patient orientation and a complete array of tabletop articulation.

Table Compatibility: 3080-R, 3080RL/SP, 3085 SP



Hand Control - 3080 RC

#BF201

Lightweight, removable, fluid-tight, and shock resistant. Also provides for easy patient orientation and a complete array of tabletop articulations.

Table Compatibility: 3080 RC



Hand Control - SurgiGraphic 6000 CE

#BF35

Lightweight, removable, fluid-tight, and shock-resistant, this remote control device includes a complete array of tabletop articulations for the anesthesiologist, including four programmable memory function buttons. The "Return-to-Level" features return table to level position with the push of a button.

Table Compatibility: SurgiGraphic 6000CE



Hand Control - SurgiGraphic 6000 (old style)

#BF33

Lightweight, removable, fluid-tight, and shock-resistant, this remote control device includes a complete array of tabletop articulations for the anesthesiologist, including four programmable memory function buttons. The "Return-to-Level" features return table to level position with the push of a button.

Table Compatibility: SurgiGraphic 6000 (Old Style)



Palm Control - Surgigraphic 6000 CE

#BF348

This unique controller gives the surgeon virtually limitless positioning of the patient relative to the C-arm with lateral and longitudinal table float, as well as all discreet motions in between. The manual Palm Control offers the surgeon effortless, precise, and remarkably smooth glide of the table top. This fluid-tight, shock-resistant control can be positioned anywhere along the table sides and has a complete array of table top articulations for the surgeon.

Table Compatibility: SurgiGraphic 6000CE, SurgiGraphic 6000FB



Palm Control - SurgiGraphic 6000 (old style)

#BF347

This unique controller gives the surgeon virtually limitless positioning of the patient relative to the C-arm with lateral and longitudinal table float, as well as all discreet motions in between. The manual Palm Control offers the surgeon effortless, precise, and remarkably smooth glide of the table top. This fluid-tight, shock-resistant control can be positioned anywhere along the table sides and has a complete array of table top articulations for the surgeon.

Table Compatibility: SurgiGraphic 6000CE, SurgiGraphic 6000FB



Foot Control - SurgiGraphic 6000

#TA00004

Provides surgeon with Trendelenburg/reverse, raise/lower and tilt control during surgical procedures. Unit plugs into receptacle next to hand control on column of table.

Table Compatibility: SurgiGraphic 6000

Special Note(s): Also compatible with 6000CE if bypassing palm control (plugging in foot control instead of palm); Substitute alternative to similar control functions on SurgiGraphic 6000 Palm Control.



Hand Control - SurgiGraphic 1027

#BF346

Lightweight, removable, fluid-tight, and shock-resistant, this remote control device includes a complete array of tabletop articulations for the anesthesiologist, including Trendelenburg/reverse, raise/lower and tilt functions.

Table Compatibility: SurgiGraphic 1027



Foot Control - 5085 SRT, 5085, 4085 and Cmax

#BF586

The foot control provides the surgeon with Trendelenburg/reverse, raise/lower and tilt control during surgical procedures.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax



Foot Control - 3080 SP, 3080-R and 3080

#BF187

Provides surgeon with Trendelenburg/reverse, raise/lower and tilt control during surgical procedures. Unit plugs into grey color-coded receptacle next to hand control on table column.

Table Compatibility: 3080, 3080-R, 3085 SP



Hand Control Bracket

#BF318

This bracket provides an optimal viewing angle of the hand control for the anesthesiologist as well as a sturdy place for using the hand control. It is attached to either the table's siderail or Side Rail Adapter by means of a Siderail Lock.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Special Note(s): Also compatible with SurgiGraphic 6000 with Siderail Adaptors (TA00012) sold separately.




Cable Wrap

#BF316

Attaches the power cord and/or hand control cord to siderail of table for storage

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax



X-ray Tops Complete Set - 5085SRT, 5085, 4085, Cmax


#BF219

Removable, radiolucent X-ray tops mount onto surgical tables. The 2-inch standoffs add elevation to the tabletop, providing space to insert flat X-ray cassettes from head, foot, or side of table. For use with existing table pads. These X-ray tops are available in Complete Sets – Head, Back, Seat, & Leg Section (BF219); Head Section only (BF174); Back Section only (BF175), Seat Section only (BF176), or Leg Section only (BF744).

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Special Note(s): Compatible with tables with Standard or Featherweight Leg Sections. For Cmax tables, the tabletop upgrade may need to be purchased in order to apply the X-ray tops listed above. For more information, call 888-478-3747.



X-ray Tops Head Section - 5085 SRT, 5085, 4085, Cmax

#BF174

Removable, radiolucent X-ray tops mount onto surgical tables. The 2-inch standoffs add elevation to the tabletop, providing space to insert flat X-ray cassettes from head, foot, or side of table. For use with existing table pads. These X-ray tops are available in Complete Sets – Head, Back, Seat, & Leg Section (BF219); Head Section only (BF174); Back Section only (BF175), Seat Section only (BF176), or Leg Section only (BF744).

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Special Note(s): For Cmax tables, the tabletop upgrade may need to be purchased in order to apply the X-ray tops listed above. For more information, call 888-478-3747.

X-ray Tops Back Section - 5085 SRT, 5085, 4085, Cmax

#BF175

Removable, radiolucent X-ray tops mount onto surgical tables. The 2-inch standoffs add elevation to the tabletop, providing space to insert flat X-ray cassettes from head, foot, or side of table. For use with existing table pads. These X-ray tops are available in Complete Sets – Head, Back, Seat, & Leg Section (BF219); Head Section only (BF174); Back Section only (BF175), Seat Section only (BF176), or Leg Section only (BF744).

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Special Note(s): For Cmax tables, the tabletop upgrade may need to be purchased in order to apply the X-ray tops listed above. For more information, call 888-478-3747.

X-ray Tops Seat Section - 5085 SRT, 5085, 4085, Cmax

#BF176

Removable, radiolucent X-ray tops mount onto surgical tables. The 2-inch standoffs add elevation to the tabletop, providing space to insert flat X-ray cassettes from head, foot, or side of table. For use with existing table pads. These X-ray tops are available in Complete Sets – Head, Back, Seat, & Leg Section (BF219); Head Section only (BF174); Back Section only (BF175), Seat Section only (BF176), or Leg Section only (BF744).

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Special Note(s): For Cmax tables, the tabletop upgrade may need to be purchased in order to apply the X-ray tops listed above. For more information, call 888-478-3747.

X-ray Tops Leg Section - 5085 SRT, 5085, 4085, Cmax

#BF744

Removable, radiolucent X-ray tops mount onto surgical tables. The 2-inch standoffs add elevation to the tabletop, providing space to insert flat X-ray cassettes from head, foot, or side of table. For use with existing table pads. These X-ray tops are available in Complete Sets – Head, Back, Seat, & Leg Section (BF219); Head Section only (BF174); Back Section only (BF175), Seat Section only (BF176), or Leg Section only (BF744).

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Compatible with tables with Standard or Featherweight Leg Sections. Maximum patient weight capacity is equivalent to that of the table.

Special Note(s): For Cmax tables, the tabletop upgrade may need to be purchased in order to apply the X-ray tops listed above. For more information, call 888-478-3747.

X-ray Tops Complete Set - 3080, 3080-R, 3085 SP

#BF202

Removable, radiolucent X-ray tops mount onto surgical tables. The 2-inch standoffs add elevation to the tabletop, providing space to insert flat X-ray cassettes from head, foot, or side of table. For use with existing table pads.

Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Special Note(s):



Cassette Tray with Handle - 14" x 17" Film

#BF137

Positions hard X-ray cassettes at any section of the surgical table. Detachable handle is graduated in 1-inch increments from 11 inches to 36 inches to facilitate cassette placement. Handle is approximately 30 inches long.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax



X-ray Top for 9 Inch Urology Extension

#BF579

Removable, radiolucent X-ray top mounts onto the 9 Inch Urology Extension. The 2-inch standoffs add elevation to the extension providing space to insert flat X-ray cassettes. Provides even surface when X-ray tops are mounted to the surgical table.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum weight capacity is the equivalent to that of the 9 Inch Urology Extension.

Usage/Procedures: Endourologic; KUB; Lithotomy

Required Accessories: BF578


Universal Ether Screen with Clark Socket

#BF074

Supports surgical drape and is mounted on a ball joint for positioning flexibility (parallel or perpendicular to patient body as the case dictates). Screen is height adjustable and locks or releases with a half-turn of control knob. Includes Clark Socket for mounting to table siderail.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



Winged Ether Screen #1401

#BF075

Features two extension arms to support surgical drapes. Control knob for each extension arm allows height and lateral adjustments. Each extension is independently height adjustable and locks or releases with a half-turn of the control knob. Includes mounting sockets.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



Transducer Holder

#BF263

Provides a rigid, stable post for securing a transducer at the appropriate height for neuro and cardiac surgery. Curved design makes it easy for anesthesiologist to position and view the transducer.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Required Accessories: BF083, BF43

Special Note(s): Also compatible with SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



Bariatric Accessory Storage Cart

#BF551

This cart is a compact storage and transport device for the full range of STERIS bariatric accessories. The eight side rails were specifically designed to accommodate each of the unique Bariatric Width Extenders. The two back rails easily hold STERIS's Split Leg Positioner or other leg holder accessories. Other accessories may be attached to interior rails or stored in the deep bin / weighted base. (Accessories not included.)



Bariatric Reverse Trendelenburg Strap (Pair)

#BF563

The broad pads on these straps safely encircle the inner thigh without binding or pinching, making them ideal for Nissen Fundoplication or other extreme Reverse Trendelenburg procedures. On this bariatric model, the safety belt buckle ensures a snug fit for patients up to 1000 lbs. These straps attach easily and securely to the table with integral rail clamps.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Reverse Trendelenburg

Optional Accessories: BF454



Knee Restraint Strap (Pair)

#BF552

The broad pads on these straps safely prevent the knee from flexing during surgery and feature a full-length Velcro closure for a snug fit. Used in conjunction with the Split Leg Positioner.

Usage/Procedures: Laparoscopic; Colo-rectal; Vascular; Extreme Reverse Trendelenburg

Special Note(s): Compatible with tables using the Split Leg Positioner (BF351) including 5085 SRT, 5085, 4085, Cmax, 3085 SP, 3080-R, and 3080 tables.



Bariatric Power-lift Stirrups (Pair)

#BF454

These stirrups feature the same grip-and-twist control handles and use the same power assist technology as the standard Power-lift Stirrups (BF453), but can accommodate patients up to 720 lbs., depending on the table used. The extra large boots and pads provide comfort and security. Includes rail clamps and pads.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 720 lbs. (377 kg) for 5085 SRT, 5085, 4085, Cmax and 3085 SP; 600 lbs for 3080-R and 400 lbs. (181 kg) for 3080 tables.

Usage/Procedures: Endourology; Gynecology; Perineal; Bariatric; Colo-rectal; Other mid-line access procedures

Optional Accessories: BF46, BF563

Replacement Pad(s): BF191 (Left); BF192 (Right)



Split Leg Positioner

#BF351

The Split Leg Positioner includes both left and right leg supports and provides a safe system for positioning patients for laparoscopic procedures as well as lower extremity orthopedic and vascular procedures. It is radiolucent allowing both the surgeon and imaging equipment excellent access to the operative site. Patient's legs can be adducted, abducted, raised or lowered for low Lithotomy procedures.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg) for 5085 SRT, 5085, 4085, Cmax, 3085 SP and 3080-R, and 400 lbs. (181 kg) for 3080 tables.

Usage/Procedures: Lower extremity procedures; Vascular; Orthopedic; Bariatric; Mid-line access to abdomen

Optional Accessories: BF552, MTPA09, MTPA10

Replacement Pad(s): BF351000337 (Left); BF351000338 (Right)



Bariatric Foot Extension

#BF542

The Bariatric Foot Extension is 28-inches wide to accommodate wider tables with Bariatric Extenders installed. It adds 8-inches of table length for taller patients and swings 90 degrees upward for foot support when in reverse Trendelenburg.


Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Required Accessories: BF133

Replacement Pad(s): BF493

Special Note(s): Bariatric Table Extensions are required (sold separately).



Bariatric Table Extenders without Pads; Complete Set - 5085SRT, 5085, 4085, Cmax

#BF661

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This complete set includes 2 extenders each for the Head, Back, Seat, and Leg sections. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Bariatric Table Extenders with 2" Pads; Complete Set - 5085SRT, 5085, 4085, Cmax

#BF643

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This complete set includes 2 extenders each for the Head, Back, Seat, and Leg sections. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.


Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Replacement Pad(s): BF652



Bariatric Table Extenders with 4" Pads; Complete Set - 5085SRT, 5085, 4085, Cmax

#BF644

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This complete set includes 2 extenders each for the Head, Back, Seat, and Leg sections. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Replacement Pad(s): BF653



Bariatric Table Extenders without Pads; Back and Seat Sections - 5085SRT, 5085, 4085, Cmax

#BF662

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This set includes 2 extenders each for the Back and Seat sections. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Bariatric Table Extenders with 2" Pads; Back and Seat Sections - 5085SRT, 5085, 4085, Cmax

#BF645

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This set includes 2 extenders each for the Back and Seat sections. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Replacement Pad(s): BF654

Bariatric Table Extenders with 4" Pads; Back and Seat Sections - 5085SRT, 5085, 4085, Cmax

#BF646

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This set includes 2 extenders each for the Back and Seat sections. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Replacement Pad(s): BF655

Bariatric Table Extenders without Pads; Head Section - 5085SRT, 5085, 4085, Cmax

#BF663

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This set includes 2 extenders each for the Head section. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Bariatric Table Extenders with 2" Pads; Head Section - 5085SRT, 5085, 4085, Cmax**#BF647**

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This set includes 2 extenders each for the Head section. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Replacement Pad(s): BF656

Bariatric Table Extenders with 4" Pads; Head Section - 5085SRT, 5085, 4085, Cmax**#BF648**

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This set includes 2 extenders each for the Head section. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Replacement Pad(s): BF657

Bariatric Table Extenders without Pads; Leg Section - 5085SRT, 5085, 4085, Cmax**#BF664**

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This set includes 2 extenders each for the Leg section. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Special Note(s): Compatible with Standard or Featherweight Leg Sections.

Bariatric Table Extenders with 2" Pads; Leg Section - 5085SRT, 5085, 4085, Cmax

#BF649

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This set includes 2 extenders each for the Leg section. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Replacement Pad(s): BF658

Special Note(s): Compatible with Standard or Featherweight Leg Sections.

Bariatric Table Extenders with 4" Pads; Leg Section - 5085SRT, 5085, 4085, Cmax

#BF651

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This set includes 2 extenders each for the Leg section. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Replacement Pad(s): BF659

Special Note(s): Compatible with Standard or Featherweight Leg Sections.

Bariatric Table Extender Pads Only 2"; Complete Set - 5085 SRT, 5085, 4085, Cmax

#BF652

This complete set of 2" pads for the Bariatric Table Extenders includes 2 pads each for the Head, Back, Seat, and Leg sections.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF661

Bariatric Table Extender Pads Only - 4"; Complete Set - 5085 SRT, 5085, 4085, Cmax

#BF653

This complete set of 4" pads for the Bariatric Table Extenders includes 2 pads each for the Head, Back, Seat, and Leg sections.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF661



Bariatric Table Extender Pads Only - 2"; Back and Seat Sections - 5085 SRT, 5085, 4085, Cmax

#BF654

This set of 2" pads for the Bariatric Table Extenders includes 2 pads each for the Back and Seat sections.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF662

Bariatric Table Extender Pads Only - 4"; Back and Seat Sections - 5085 SRT, 5085, 4085, Cmax

#BF655

This set of 4" pads for the Bariatric Table Extenders includes 2 pads each for the Back and Seat sections.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF662

Bariatric Table Extender Pads Only - 2"; Head Section - 5085 SRT, 5085, 4085, Cmax

#BF656

This set of 2" pads for the Bariatric Table Extenders includes 2 pads for the Head section.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF663

Bariatric Table Extender Pads Only - 4"; Head Section - 5085 SRT, 5085, 4085, Cmax

#BF657

This set of 4" pads for the Bariatric Table Extenders includes 2 pads for the Head section.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF663

Bariatric Table Extender Pads Only - 2"; Leg Section - 5085 SRT, 5085, 4085, Cmax

#BF658

This set of 2" pads for the Bariatric Table Extenders includes 2 pads for the Leg Section.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF664

Bariatric Table Extender Pads Only - 4"; Leg Section - 5085 SRT, 5085, 4085, Cmax

#BF659

This set of 4" pads for the Bariatric Table Extenders includes 2 pads for the Leg section.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF664

Bariatric Table Extenders with 2" Pads; Complete Set - 3080, 3085

#BF541

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This complete set includes 2 extenders each for the Head, Back, Seat, and Foot sections. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.


Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Replacement Pad(s): BF543 (2"); BF556 (4")



Bariatric Table Extenders with 2" Pads; Back and Seat - 3080, 3085

#BF553

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This set includes 2 extenders each for the Back and Seat sections. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Replacement Pad(s): BF544 (2"); BF557 (4")

Bariatric Table Extenders with 2" Pads; Head Section - 3080, 3085**#BF554**

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This set includes 2 extenders each for the Back section. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Replacement Pad(s): BF545 (2"); BF558 (4")

Bariatric Table Extenders with 2" Pads; Foot Section - 3080, 3085**#BF555**

These radiolucent modular extenders increase the width of your table to a full 28-inches to accommodate bariatric patients. This set includes 2 extenders each for the Foot section. The Bariatric Table Extenders can be conveniently stored on the Bariatric Accessory Cart (BF551); sold separately.

Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Optional Accessories: BF542, BF551

Replacement Pad(s): BF546 (2"); BF559 (4")


Bariatric Table Extender Pads Only; 4" Complete Set - 3080, 3085SP**#BF556**

This set of 4" pads for the Bariatric Table Extenders includes 2 pads each for the Head, Back, Seat, and Leg sections.

Table Compatibility: 3080, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

**Bariatric Table Extender Pads Only; 4" Back and Seat sections - 3080, 3085SP****#BF557**

This set of 4" pads for the Bariatric Table Extenders includes 2 pads each for the Back and Seat sections.

Table Compatibility: 3080, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Bariatric Table Extender Pads Only; 4" Head section - 3080, 3085SP

#BF558

This set of 4" pads for the Bariatric Table Extenders includes 2 pads for the Head section.

Table Compatibility: 3080, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Bariatric Table Extender Pads Only; 2" Complete Set - 3080, 3085SP

#BF543

This set of 2" pads for the Bariatric Table Extenders includes 2 pads each for the Head, Back, Seat, and Leg sections.

Table Compatibility: 3080, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Bariatric Table Extender Pads Only; 2" Back and Seat sections - 3080, 3085SP

#BF544

This set of 2" pads for the Bariatric Table Extenders includes 2 pads each for the Back and Seat sections.

Table Compatibility: 3080, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Bariatric Table Extender Pads Only; 2" Head section - 3080, 3085SP

#BF545

This set of 2" pads for the Bariatric Table Extenders includes 2 pads for the Head section.

Table Compatibility: 3080, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Bariatric Table Extender Pads Only; 2" Foot section - 3080, 3085SP

#BF546

This set of 2" pads for the Bariatric Table Extenders includes 2 pads for the Foot section.

Table Compatibility: 3080, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Bariatric Fluid Collection System

#BF56

The Bariatric Fluid Collection System includes adapters that extend to accommodate the additional width created by the Bariatric Table Extenders applied to the table. It is ideal for use in cystoscopic or hysteroscopic procedures, and comes with the flexible drape support and a fluid collection bag. The flexible drape support device quickly attaches to the Seat Section of the Bariatric Table Extender during set up, and is easily manipulated intra-procedurally to obtain maximum unobstructed surgical site access. Both the Sterile and Non-sterile Fluid Collector Bags fit either the standard or Bariatric model drape support. The built-in flexible drain hose works with a kick bucket or a suction canister.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Urologic; Cystoscopic; Hysteroscopic

Special Note(s): Also compatible with 3080 with Bariatric Table Extenders applied. Set of 10 Fluid Collection Bags: Sterile (BF458) or Non-sterile (BF459) are required (sold separately).



Narrow Headrest

#BF592

This headrest is ideal for neurological procedures. The Wrist Rest (BF593) attaches to the Narrow Headrest to support the surgeon's hands during long procedures.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Neurological; Ophthalmic; ENT; Plastic Surgery

Optional Accessories: BF593

Replacement Pad(s): BF491

Special Note(s): The Eye-ENT-Neuro Headrest Adaptor is required (sold separately). Order BF185 for 3085 tables, BF173 for 5085, 4085, and Cmax tables, or BF309 for 5085SRT.



Cone-Shaped Headrest

#BF215

This headrest fits directly into a table's standard headrest area and provides easier access to head area for eye, ENT, plastic surgery, and neurosurgical applications. Includes a 2-inch pad and a 4-inch pad for use with X-ray tops.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Eye; ENT; Neurosurgical; Plastic Surgery

Replacement Pad(s): BF83 (2"); BF831 (4")



Eye-ENT-Neuro Headrest Adaptor - 3080, 3080-R, 3085

#BF185

This specialized adaptor accommodates all standard Eye-ENT-Neuro attachments. Conveniently inserts into either end of the 3085 SP surgical table.

Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is 500 pounds (227 kg) for 3085 and 300 (136 kg) for 3080 and 3080-R.

Usage/Procedures: Neurosurgical; Ophthalmic; ENT; Plastic Surgery



Eye-ENT-Neuro Headrest Adaptor - 5085, 4085, Cmax

#BF173

This specialized adaptor accommodates all standard Eye-ENT-Neuro attachments. Conveniently inserts into either end of the 5085, 4085, or Cmax surgical tables.

Table Compatibility: 4085, 5085, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Neurosurgical; Ophthalmic; ENT; Plastic Surgery



Eye-ENT-Neuro Headrest Adaptor - 5085 SRT

#BF309

This specialized adaptor accommodates all standard Eye-ENT-Neuro attachments. Conveniently inserts into either end of the 5085 SRT surgical table.

Table Compatibility: 5085 SRT

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Neurosurgical; Ophthalmic; ENT; Plastic Surgery



Frontal Sitting Attachment

#BF01


This accessory mounts neurosurgical headrests to head or foot end of 3085 SP, 3080-R or 3080 tables or to head end of 5085 SRT, 5085, 4085 or Cmax tables. Ball joints allow full angular positioning of headrest.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 300 lbs. (136 kg).

Usage/Procedures: For support of patient's head from behind when patient is in a seated position for temporal, parietal, and frontal surgical procedures.

Special Note(s): Eye-ENT-Neuro Headrest Adaptor and either Multi-Poise Headrest (BF038) or Horseshoe Headrest (BF036) are required (sold separately).



Horizontal Attachment #1800

#BF008

Mounts neurosurgical headrests to head or foot end of 3085 SP, 3080-R and 3080 tables or to head end of 5085 SRT, 5085, 4085 or Cmax tables. Ball joints allow full angular positioning of headrest.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 300 lbs. (136 kg).

Usage/Procedures: For supporting patient's head in supine, prone, or lateral positions.

Special Note(s): Eye-ENT-Neuro Headrest Adaptor and either Horizontal Attachment #1800 (BF008) or Frontal Sitting Attachment (BF01) are required (sold separately).



Horseshoe Headrest

#BF036

For surgical procedures involving nose, paranasal, sinuses, trachea, and larynx approaches, as well as some neurological procedures.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 300 lbs. (136 kg).

Usage/Procedures: Cradles patient's head for eye, ENT, neuro, and plastic surgery.

Replacement Pad(s): BF492

Special Note(s): Eye-ENT-Neuro Headrest Adaptor and either Horizontal Attachment (BF008) or Frontal Sitting Attachment (BF01) are required (sold separately).



Posterior Sitting Attachment #1801

#BF011


Mounts neurosurgical headrests to table siderails. Ball joints allow lateral and vertical adjustment.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 300 lbs. (136 kg).

Usage/Procedures: For support of patient's head from the front when patient is in a seated position for suboccipital and parietal surgical procedures.

Special Note(s): Multi-Poise Headrest (BF038) and either Clark Sockets (BF083) or Multi-Accessory Clamps (BF43) are required (sold separately).



Multi-Poise Headrest

#BF038

This headrest permits posturing flexibility for neurosurgical techniques. Three sterilizable spring-laden pins assure maximum rigidity of patient's head without discomfort.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 300 lbs. (136 kg).

Usage/Procedures: Used with patient in supine, prone, lateral (park bench), or sitting position. Neurosurgical Craniotomies

Special Note(s): Eye-ENT-Neuro Headrest Adaptor and either Horizontal Attachment #1800 (BF008), Frontal Sitting Attachment #1802 (BF01), or Posterior Sitting Attachment #1801 (BF011) are required (sold separately).



Morgan Pelvic Rest

#BF135

Supports patient knee-chest position with knees on extended foot section. Adjustable lateral patient support mounts to table siderails. Includes a 4-inch thick pad with cutout, two padded supports, and an extended foot section with pad.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 300 lbs. (136 kg).

Usage/Procedures: Spinal; Laminectomy

Required Accessories: BF133, BF45

Replacement Pad(s): BF129 (Table pad w/cutout); BF132 (Foot section pad); BF131 (Lateral pad)

Special Note(s): Clark Sockets (BF083) or Multi-Accessory Clamps (BF43) are required (sold separately).



Andrews Frame Adaptor Kit

#BF305

Adapts 3085 SP, 3080-R and 3080 table to accommodate older Andrews Frame models supplied by Mizuho OSI. An adaptor is not required for newer Andrews Frame models. Check Mizuho OSI model for compatibility. Includes siderail extenders, six flat head screws, and a new leg section with spaldite top.

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg). This weight capacity does not apply to the Andrews Frame, which is not supplied by STERIS.

Usage/Procedures: Neurologic procedures involving spinal surgery

Special Note(s): Andrews Frame Compatibility: Model 914 is compatible with 1080 and 2080 tables and also compatible with 3085 SP, 3080-R, and 3080 with adapter kit. Model 916 is compatible with 3085 Sp, 3080-R, and 3080 tables. And, Model 919 is compatible with 2080 IA tables.



Wrist Rest for Narrow Headrest

#BF593

This Wrist Rest attaches to the Narrow Headrest (BF592) to support the surgeon's hands during long procedures.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Neurological; Ophthalmic; ENT; Plastic Surgery

Required Accessories: BF592



Narrow Headrest

#BF592

This headrest is ideal for neurological procedures. The Wrist Rest (BF593) attaches to the Narrow Headrest to support the surgeon's hands during long procedures.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Neurological; Ophthalmic; ENT; Plastic Surgery

Optional Accessories: BF593

Replacement Pad(s): BF491

Special Note(s): The Eye-ENT-Neuro Headrest Adaptor is required (sold separately). Order BF185 for 3085 tables, BF173 for 5085, 4085, and Cmax tables, or BF309 for 5085SRT.



Cone-Shaped Headrest

#BF215

This headrest fits directly into a table's standard headrest area and provides easier access to head area for eye, ENT, plastic surgery, and neurosurgical applications. Includes a 2-inch pad and a 4-inch pad for use with X-ray tops.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Eye; ENT; Neurosurgical; Plastic Surgery

Replacement Pad(s): BF83 (2"); BF831 (4")



Eye-ENT-Neuro Headrest Adaptor - 3080, 3080-R, 3085

#BF185

This specialized adaptor accommodates all standard Eye-ENT-Neuro attachments. Conveniently inserts into either end of the 3085 SP surgical table.

Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is 500 pounds (227 kg) for 3085 and 300 (136 kg) for 3080 and 3080-R.

Usage/Procedures: Neurosurgical; Ophthalmic; ENT; Plastic Surgery



Eye-ENT-Neuro Headrest Adaptor - 5085, 4085, Cmax

#BF173

This specialized adaptor accommodates all standard Eye-ENT-Neuro attachments. Conveniently inserts into either end of the 5085, 4085, or Cmax surgical tables.

Table Compatibility: 4085, 5085, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Neurosurgical; Ophthalmic; ENT; Plastic Surgery



Eye-ENT-Neuro Headrest Adaptor - 5085 SRT

#BF309

This specialized adaptor accommodates all standard Eye-ENT-Neuro attachments. Conveniently inserts into either end of the 5085 SRT surgical table.

Table Compatibility: 5085 SRT

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Neurosurgical; Ophthalmic; ENT; Plastic Surgery



Wrist Rest for Narrow Headrest

#BF593

This Wrist Rest attaches to the Narrow Headrest (BF592) to support the surgeon's hands during long procedures.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Neurological; Ophthalmic; ENT; Plastic Surgery

Required Accessories: BF592



Orthopedic Extension Accessory Package

#BF806

Includes:

- Accessory Cart (BF24)
- Adult Wedge-shaped Sacral Rest with Pad (BF00066)
- Adult Perineal Post with Pad (BF00065)
- (2) L-shaped Traction Extensions (BF0003)
- Anesthesia Armboard with 2-inch Pad (BF803)
- (2) TenZor™ Carbon Fiber Traction Unit Assemblies with Padded Boots (BF72)
- Universal Legholder (BF266)
- (2) Traction Accessory Clamps (BF0002)
- Multi-Posture Armboard with 2-inch Pad (BF432)
- Patient Transfer Board with 4-inch Pad (BF214)
- Clark Socket (BF083)
- Ortho Extension Setup Guide (P150832316)

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Enables hip pinning and supine femoral nailing procedures.

Universal Legholder (Each)

#BF266

Each Universal Legholder is contoured and padded for comfortable leg support, and is mounted on a ball joint for positioning flexibility. Includes strap.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Usage/Procedures: Orthopedic fracture procedures to support patient's unaffected leg

Special Note(s): Clark Sockets (BF083) or Multi-Accessory Clamps (BF43) for Cmax, 3080-R, 3085SP are required (sold separately). Traction-Accessory Clamp (BF0002) for the OrthoVision table and Orthopedic Extension is required (sold separately).



Shoulder Chair

#BF594

This Shoulder Chair, used for open or arthroscopic procedures, features a unique, multi-axis head positioning capability, which allows for easy intubation and head posturing. Powered positioning of the back section inclines from 0 to 90 degrees. Removable shoulder panels provide excellent surgical access to each shoulder. Includes head restraint device and an integral lateral brace for left or right side patient positioning.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Shoulder Arthroscopy

Optional Accessories: BF039, BF461

Replacement Pad(s): BF442 (Body Section), BF443 (Right Shoulder Panel), BF445 (Lateral Brace), BF446 (Headrest)



Head Restraint Straps

#BF039

These single-use straps provide chin and head restraint when used with the Shoulder Chair (BF594). Each box includes 12 head straps and 12 chin straps.

Table Compatibility: Shoulder Chair (BF594)

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Shoulder Arthroscopy



Shoulder Table

#BF437

Used for arthroscopic or open shoulder surgery. Powered positioning of back section inclines from 0 degrees to 90 degrees. Removable shoulder panels provide an 8-inch access area behind each shoulder.

Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Arthroscopic; Open Shoulder

Optional Accessories: BF039, BF452, BF461

Special Note(s): For lateral support, use Lateral Brace Assemblies (BF602) or Lateral Safety Strap (BF562) sold separately.



Shoulder Table Chin Cups

#BF452

These chin cups are used in conjunction with the Shoulder Table and provide stability for the head and chin during shoulder procedures. Each box contains 10 chin cups.

Table Compatibility: Shoulder Table (BF437)

Usage/Procedures: Arthroscopic; Open Shoulder



Patient Transfer Board with 4-inch Pad for Tables with Orthopedic Extension

#BF214

Supports patient's legs during transfer to table. Recommended for leg support only. Includes a 4-inch removable foam pad. Pad attaches to transfer board with Velcro strips.

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Replacement Pad(s): BF219

Special Note(s): Compatible with Orthopedic Extension applied to 5085 SRT; 5085; 4085; Cmax; 3085 SP; 3080-R; and 3080 tables.



Patient Transfer Board with 2-inch Pad - OrthoVision

#BF213

Lightweight, duotone-covered, aluminum frame supports patient's legs during transfer to and from table. Recommended for leg support only. Fits into access holes in sacral area of OrthoVision tables.

Table Compatibility: OrthoVision

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Replacement Pad(s): BF818



Hip Spica Assembly - Adult

#BF00031

Elevated sacral rest provides added height for casting in sacrolumbar area. Includes elevated back board, elevated sacral rest, body slat, and perineal post.

Table Compatibility: OrthoVision

Weight Capacity: Maximum weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Body Casting

Special Note(s): Must be custom ordered to fit the Orthopedic Extension when applied to general surgical tables. Contact your sales rep for details.



Hip Spica Assembly - Child

#BF00038

Elevated sacral rest provides added height for casting in sacrolumbar area. Includes elevated back board: 22"L x 20"W; elevated sacral rest: 6"L x 4"W; body slat: 20"L x 2"W, 2" pad.

Table Compatibility: OrthoVision

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Body Casting



Hip Spica Assembly for Ortho Extension - Child

#BF177

For use on the Orthopedic Extension. Elevated sacral rest provides added height for casting in sacrolumbar area. Includes elevated back board: 22"L x 20"W; elevated sacral rest: 6"L x 4"W; body slat: 20"L x 2"W, 2" pad.

Usage/Procedures: Body Casting

Required Accessories: BF133

Special Note(s): Orthopedic Extension applied to 5085 SRT, 5085, 4085, Cmax, 3085 SP, 3080-R, 3080

Straight Traction Extension

#BF268

When patient is positioned on the OrthoVision, this extension is used to adjust the level of the traction unit assembly and the patient's leg relative to the abductor bars. The extension is in a position that lowers the patient's leg relative to the abductor bar for tibia nailing.

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Lower Extremity Fractures

Required Accessories: BF72



L-Shaped Traction Extension

#BF0003


Holds the Traction Unit Assembly to provide greater positioning flexibility. For example, it can be used to provide additional length or greater abduction during orthopedic procedures.

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Orthopedic Fractures

Required Accessories: BF72



TenZor™ Carbon Fiber Traction Unit Assembly

#BF72

This lightweight carbon fiber traction assembly attaches to the abductor bars of orthopedic tables by means of a Traction Accessory Clamp to hold foot traction boot and pad (included) in place to pull traction. The tension-adjusting handle at the distal end allows precise adjustment of traction being pulled, and easy-to-read numbers inform amount of torque being used. It also holds the optional traction bow/block and the optional skeletal traction block. Includes: Carbon Fiber Support Plate (BF729), Adult Foot Pad (BF538), and Traction Unit. For use with coupler connection. Sold as single unit.

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Orthopedic Traction

Required Accessories: BF0002

Special Note(s): Foot Traction Boot Pads Child (BF539) or Obese Adult (BF548) can also be ordered as needed.



Foot Traction Boot Pads with Straps - Adult

#BF538

This replaceable traction boot pads fit snugly against the Carbon Fiber Support Plate (BF729) used in the TenZor™ Carbon Fiber Traction Unit Assembly (BF72). Adult pads come standard with the TenZor™ Carbon Fiber Traction Unit Assembly (BF72). Sold as pair. Straps included.



Foot Traction Boot Pads with Straps - Obese Adult

#BF548

This replaceable traction boot pads fit snugly against the Carbon Fiber Support Plate (BF729) used in the TenZor™ Carbon Fiber Traction Unit Assembly (BF72). Adult pads come standard with the TenZor™ Carbon Fiber Traction Unit Assembly (BF72). Sold as pair. Straps included.



Foot Traction Boot Pads with Straps - Child

#BF539

This replaceable traction boot pads fit snugly against the Carbon Fiber Support Plate (BF729) used in the TenZor™ Carbon Fiber Traction Unit Assembly (BF72). Child-sized boots are adjustable to accommodate sizes from pre-adolescents to smaller adults. Sold as pair. Straps included.



Traction Bow/Block, Skeletal

#BF725

Inserts into the TenZor™ Carbon Fiber Traction Unit Assembly (BF72) to allow use of Steinman or Kirschner holders/bows for skeletal traction. For use with coupler connection. Sold individually.

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Orthopedic Fractures with Skeletal Traction

Required Accessories: BF72



Traction Bow, Skeletal

#BF724


This lightweight carbon fiber device, sold individually, inserts into the TenZor™ Carbon Fiber Traction Unit Assembly (BF72) for skeletal traction when performing procedures in supine or lateral positions. The flexibility of bow permits use with patients of various sizes. The Carbon Fiber Traction Bow/Block allows for more precise and stable positioning of a fracture and is radiolucent. For use with coupler connection.

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Orthopedic Fractures with Skeletal Traction

Required Accessories: BF72



Traction Unit Assembly

#BF0006


Attaches to abductor bars of orthopedic tables by means of a Traction Accessory Clamp to hold traction boots in place to pull traction. Also holds optional traction bow/block and optional traction bow. This version of the Traction Unit Assembly is the former style with the flat plate connection and brown Foot Traction Boots.

Table Compatibility: Orthopedic Extension, OrthoVision

Usage/Procedures: Orthopedic Traction

Required Accessories: BF0002

Special Note(s): Foot Traction Boots Adult (BF00023), Child (BF00032), or Obese Adult (BF0008) are required (sold separately).



Traction Block, Skeletal

#BF00028

Inserts into traction unit assembly to allow use of Steinman or Kirschner bows for skeletal traction. For use with older style plate connection.

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Patient weight limitation is 400 lbs. (181 kg).

Usage/Procedures: Orthopedic Fractures

Required Accessories: BF0006

Special Note(s): Two each Traction Assembly Clamps (BF0002) are required (sold separately).



Traction Bow/Block

#BF371

Inserts into traction unit assembly for skeletal traction when performing procedures in supine or lateral positions. Flexibility of bow permits use with patients of various sizes. Fits into traction unit assembly at foot plate slot, allowing for more precise and stable positioning of fracture. For use with older style plate connection.

Table Compatibility: Orthopedic Extension, OrthoVision

Usage/Procedures: Orthopedic Fractures

Required Accessories: BF0006

Special Note(s): For use with the BF0006 former style Traction Unit Assembly with the flat foot plate.




Traction Accessory Clamp

#BF0002

Can be attached to either side of abductor bars for mounting orthopedic accessories.

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).



Carbon Fiber Tibia/Lateral Countertraction Support

#BF0062

This radiolucent support provides comfortable support of the leg in a lateral position, as well as countertraction support for intramedullary and tibia nailing procedures. It features a padded "J-bar," a perineal post adaptor, and a pivot arm (perineal post not included).

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Femoral Intramedullary; Tibia Nailing Procedures; Lateral Positioning

Required Accessories: BF0002

Optional Accessories: BF00065



Arthroscopy Knee Holder

#BF717

The Arthroscopy Knee Holder attaches to the seat section of the table with a Clark Socket (not included) and is typically used for the "off the end" of the table arthroscopy technique. The diameter of the knee holder can be adjusted as surgical needs dictate. A three piece cushion set is included.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 218 lbs. (99 kg).

Usage/Procedures: Knee Arthroscopy

Replacement Pad(s): BF718

Special Note(s): Clark Sockets (BF083) or Multi-Accessory Clamps (BF43) are required (sold separately).



Carbon Fiber T-Shaped Sacral Rest with Pad

#BF00073

Attaches to orthopedic table to support the patient's sacral area. Maximum width measures 13-inches.

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Obese Patients and Pelvic Fractures

Required Accessories: BF00065

Replacement Pad(s): BF824



Carbon Fiber Wedge-Shaped Sacral Rest with Pad - Adult

#BF00066

Attaches to orthopedic table to support the patient's sacral area.


Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Ortho procedures using the Fracture table in supine position

Required Accessories: BF00065

Replacement Pad(s): BF825



Wedge-Shaped Sacral Rest with Pad - Child

#BF00048

Attaches to orthopedic table to support the patient's sacral area. Maximum width measures approximately 4 inches.

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Maximum patient weight capacity is 150 lbs. (68 kg).

Usage/Procedures: Pediatric ortho procedures using the Fracture table in supine position

Required Accessories: BF41

Optional Accessories: NO001

Replacement Pad(s): BF826



Carbon Fiber Hip Rest with Pad

#BF00072

Attaches to sacral area of table for support of patient's hip in lateral position. Total carbon composite. Includes 2-inch pad for pressure management.

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Lower extremity orthopedic procedures in lateral position

Replacement Pad(s): BF827



Perineal Post with Pad - Adult

#BF00065

Post provides countertraction and is made of radiolucent carbon fiber to enhance C-arm imaging. The post pad relieves pressure on the perineal area and is made of disposable foam.

Table Compatibility: Orthopedic Extension, OrthoVision

Usage/Procedures: Orthopedic procedures in supine position

Replacement Pad(s): NO001 (Case of 10)

Special Note(s): Wedge-shaped Sacral Rest - Adult (BF00066) or T-shaped Sacral Rest (BF00073) are required (sold separately).



Perineal Post without Pad - Child

#BF41

Provides countertraction and is used in supine positions with a sacral rest. Measures approximately 7-inches including the insert.

Table Compatibility: Orthopedic Extension, OrthoVision

Usage/Procedures: Orthopedic procedures in supine position

Required Accessories: BF00048, NO001



Perineal Post Pad - Disposable

#NO001

Relieves pressure from the perineal area. Fits standard perineal posts and can be trimmed to fit child perineal posts. Made of latex-free disposable foam. Case of 10.

Table Compatibility: Orthopedic Extension, OrthoVision

Usage/Procedures: Orthopedic procedures in supine position



Surgical Drape Support

#BF0018

Consists of two flat drape-supports on round, angled, stainless steel bars that attach to the head and foot end of the table. A Clark Socket is used to attach the drape support to the head end and the traction assembly clamp is used to attach the drape support to the abductor bar. Sterile field can be maintained if table is raised and/or lowered during the surgical procedure. Flat drape-support bars can be positioned in a "V" shape to enhance C-arm imaging. Clark socket and traction accessory clamp included.

Table Compatibility: Orthopedic Extension, OrthoVision



Foot Traction Boots (Old Style) - Adult (Pair)

#BF0023

Made of flexible material with foam inserts that secure to patient's feet with Velcro straps, these boots slide onto the foot plate of the Traction Unit Assembly. For use with the ratchet style stainless steel traction component.

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Patient weight limitation is 400 lbs. (181 kg).

Usage/Procedures: Orthopedic Fractures; Hip Spica

Required Accessories: BF0006

Special Note(s): These boots are replacements for use with the BF0006 former style Traction Unit Assembly with the flat foot plate.



Foot Traction Boots (Old Style) - Child (Pair)

#BF0032

Made of flexible material with foam inserts that secure to patient's feet with Velcro straps, these boots slide onto the foot plate of the Traction Unit Assembly. For use with the ratchet style stainless steel traction component. Boots are adjustable to accommodate sizes from pre-adolescents to smaller adults.

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Patient weight limitation is 400 lbs. (181 kg).

Usage/Procedures: Orthopedic Fractures; Hip Spica

Required Accessories: BF0006

Special Note(s): These boots are replacements for use with the BF0006 former style Traction Unit Assembly with the flat foot plate.



Foot Traction Boots - Obese Adult (Pair)

#BF0008

Made of flexible material with foam inserts that secure to patient's feet with Velcro straps, these boots slide onto the foot plate of the Traction Unit Assembly. For use with the ratchet style stainless steel traction component.

Table Compatibility: Orthopedic Extension, OrthoVision

Weight Capacity: Patient weight limitation is 400 lbs. (181 kg).

Usage/Procedures: Orthopedic Fractures; Hip Spica

Required Accessories: BF0006

Special Note(s): These boots are replacements for use with the BF0006 former style Traction Unit Assembly with the flat foot plate.



Standard Accessory Package

#BF807

Includes:

TLT Tabletop Pad (BF323)

Wedge-shaped Sacral Rest with Pad - Adult (BF00066)

Perineal Post with Pad - Adult (BF00065)

(2) L-shaped Traction Extensions (BF0003)

Anesthesia Armboard with 2-inch Pad (BF803)

(2) TenZor Carbon Fiber Traction Unit Assembly with Boots (BF72)

Universal Legholder (BF266)

(2) Traction Accessory Clamp (BF0002)

Multi-Posture Armboard with 2-inch Pad (BF432)

Patient Transfer Board with 2-inch Pad (BF213)

Clark Socket (BF083)

Table Compatibility: OrthoVision

Usage/Procedures: For Orthopedic procedures using OrthoVision table

X-ray Tops Complete Set - 5085SRT, 5085, 4085, Cmax

#BF219

Removable, radiolucent X-ray tops mount onto surgical tables. The 2-inch standoffs add elevation to the tabletop, providing space to insert flat X-ray cassettes from head, foot, or side of table. For use with existing table pads. These X-ray tops are available in Complete Sets – Head, Back, Seat, & Leg Section (BF219); Head Section only (BF174); Back Section only (BF175), Seat Section only (BF176), or Leg Section only (BF744).

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Special Note(s): Compatible with tables with Standard or Featherweight Leg Sections. For Cmax tables, the tabletop upgrade may need to be purchased in order to apply the X-ray tops listed above. For more information, call 888-478-3747.



X-ray Tops Head Section - 5085 SRT, 5085, 4085, Cmax

#BF174

Removable, radiolucent X-ray tops mount onto surgical tables. The 2-inch standoffs add elevation to the tabletop, providing space to insert flat X-ray cassettes from head, foot, or side of table. For use with existing table pads. These X-ray tops are available in Complete Sets – Head, Back, Seat, & Leg Section (BF219); Head Section only (BF174); Back Section only (BF175), Seat Section only (BF176), or Leg Section only (BF744).

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Special Note(s): For Cmax tables, the tabletop upgrade may need to be purchased in order to apply the X-ray tops listed above. For more information, call 888-478-3747.

X-ray Tops Back Section - 5085 SRT, 5085, 4085, Cmax

#BF175

Removable, radiolucent X-ray tops mount onto surgical tables. The 2-inch standoffs add elevation to the tabletop, providing space to insert flat X-ray cassettes from head, foot, or side of table. For use with existing table pads. These X-ray tops are available in Complete Sets – Head, Back, Seat, & Leg Section (BF219); Head Section only (BF174); Back Section only (BF175), Seat Section only (BF176), or Leg Section only (BF744).

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Special Note(s): For Cmax tables, the tabletop upgrade may need to be purchased in order to apply the X-ray tops listed above. For more information, call 888-478-3747.

X-ray Tops Seat Section - 5085 SRT, 5085, 4085, Cmax

#BF176

Removable, radiolucent X-ray tops mount onto surgical tables. The 2-inch standoffs add elevation to the tabletop, providing space to insert flat X-ray cassettes from head, foot, or side of table. For use with existing table pads. These X-ray tops are available in Complete Sets – Head, Back, Seat, & Leg Section (BF219); Head Section only (BF174); Back Section only (BF175), Seat Section only (BF176), or Leg Section only (BF744).

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Special Note(s): For Cmax tables, the tabletop upgrade may need to be purchased in order to apply the X-ray tops listed above. For more information, call 888-478-3747.

X-ray Tops Leg Section - 5085 SRT, 5085, 4085, Cmax

#BF744

Removable, radiolucent X-ray tops mount onto surgical tables. The 2-inch standoffs add elevation to the tabletop, providing space to insert flat X-ray cassettes from head, foot, or side of table. For use with existing table pads. These X-ray tops are available in Complete Sets – Head, Back, Seat, & Leg Section (BF219); Head Section only (BF174); Back Section only (BF175), Seat Section only (BF176), or Leg Section only (BF744).

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Compatible with tables with Standard or Featherweight Leg Sections. Maximum patient weight capacity is equivalent to that of the table.

Special Note(s): For Cmax tables, the tabletop upgrade may need to be purchased in order to apply the X-ray tops listed above. For more information, call 888-478-3747.

X-ray Tops Complete Set - 3080, 3080-R, 3085 SP

#BF202

Removable, radiolucent X-ray tops mount onto surgical tables. The 2-inch standoffs add elevation to the tabletop, providing space to insert flat X-ray cassettes from head, foot, or side of table. For use with existing table pads.

Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Special Note(s):



36" Carbon Fiber Fluoro Extension

#BF642

The Carbon Fiber Fluoro Extension features metal-free radiolucent beveled edges for 360 degrees imaging and increases upper or lower body imaging of up to 55 inches. This 36-inch long accessory broadens the surgical capabilities of STERIS general surgical tables for cardiovascular, AV fistulas, pain management, and other procedures requiring the enhanced imaging ability typically found on image guided surgical tables. Includes 2-inch pad.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Replacement Pad(s): BF171



Fem / Pop Board

#BF287

Increases lower body imaging capability to 55 inches of continuous imaging coverage. Radiolucent top extends approximately 22-inches and inserts into leg section of table. Includes a 2-inch pad and X-ray top. No support leg needed. For support of patient's legs only.

Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Cardiovascular procedures requiring C-arm

Replacement Pad(s): BF817



IA Extender with 2" Pad and X-ray Top

#BF188

This accessory serves two purposes: (1) for patients in Normal Orientation, the IA Extension is attached between the head and back section of the table to add 8" of imaging length to accommodate taller patients and extends upper body fluoroscopic coverage, and (2) for patients in Reverse Orientation, the IA Extension is attached to the back section of the table to lengthen the imaging window to 24" - ideal for endourology applications.

Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Endourologic; KUB; Lithotomy

Replacement Pad(s): BF829

Special Note(s): Patient Transfer Board (BF204) is recommended for lithotomy positioning or taller patients.



Cassette Tray with Handle - 14" x 17" Film

#BF137

Positions hard X-ray cassettes at any section of the surgical table. Detachable handle is graduated in 1-inch increments from 11 inches to 36 inches to facilitate cassette placement. Handle is approximately 30 inches long.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax



9 Inch Urology Extension

#BF578

This lightweight 9-inch long radiolucent extension provides additional fluoro imaging length for imaging of kidney, ureter and bladder (KUB) anatomy of tall patients. It easily attaches to the seat section of the table once the leg section is removed, and includes standard side rails for applying leg holders, or Fluid Collector System. The perineal cutout will accommodate a metal drain pan.


Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 600 lbs. (272 kg) with full slide.

Required Accessories: BF204

Replacement Pad(s): BF83400053

Special Note(s): X-ray top (BF579) sold separately.



IA Extender with 2" Pad and X-ray Top

#BF188

This accessory serves two purposes: (1) for patients in Normal Orientation, the IA Extension is attached between the head and back section of the table to add 8" of imaging length to accommodate taller patients and extends upper body fluoroscopic coverage, and (2) for patients in Reverse Orientation, the IA Extension is attached to the back section of the table to lengthen the imaging window to 24" - ideal for endourology applications.

Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Usage/Procedures: Endourologic; KUB; Lithotomy

Replacement Pad(s): BF829

Special Note(s): Patient Transfer Board (BF204) is recommended for lithotomy positioning or taller patients.



Power-lift Stirrups (Pair)

#BF453

The integral power assist mechanism in these stirrups provides excellent surgeon access and virtually effortless patient positioning of patients up to 500 lbs. The unique grip-and-twist control handles allow staff to adjust the angle of Lithotomy and degree of ab/adduction during set up and during the procedure without lifting the sterile drape. Includes Rail Clamps and pads.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 500 lbs. (227 kg).

Usage/Procedures: Endourology; Gynecology; Perineal; Other mid-line access procedures

Optional Accessories: BF455, BF46

Replacement Pad(s): BF85000381 (Left); BF85000382 (Right)



Bariatric Power-lift Stirrups (Pair)

#BF454

These stirrups feature the same grip-and-twist control handles and use the same power assist technology as the standard Power-lift Stirrups (BF453), but can accommodate patients up to 720 lbs., depending on the table used. The extra large boots and pads provide comfort and security. Includes rail clamps and pads.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 720 lbs. (377 kg) for 5085 SRT, 5085, 4085, Cmax and 3085 SP; 600 lbs for 3080-R and 400 lbs. (181 kg) for 3080 tables.

Usage/Procedures: Endourology; Gynecology; Perineal; Bariatric; Colo-rectal; Other mid-line access procedures

Optional Accessories: BF46, BF563

Replacement Pad(s): BF191 (Left); BF192 (Right)



Lithotomy Legholders - Adult (Pair)

#BF054


The Lithotomy Legholders provide cantilevered suspension-strap support for legs. Height and angulation are adjustable for full positioning flexibility.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum patient weight capacity is 300 lbs. (136 kg).

Usage/Procedures: OB/GYN; Urologic

Special Note(s): Clark Sockets (BF083) or Multi-Accessory Clamps (BF43) are required (sold separately).



Multi-Position Legholders (Pair)

#BF433

The Multi-position Legholders are designed to accommodate varying positioning requirements. Boot design prevents uncontrolled abduction, even when unlocked. Telescoping knee joints provide added flexibility. Mount to table side rails. Includes one set of disposable foam insert pads.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Gynecologic; Urologic; Laparoscopic Nissen Fundoplication

Special Note(s): Clark Sockets (BF083) or Multi-Accessory Clamps (BF43) are required (sold separately).



Drain Pan

#BF186

The drain pan attaches to the perineal edge of the seat section of the table, the IA Extender (BF188), or the 9-inch Urology Extension (BF578). Provides fast drainage. Includes a stainless steel pan measuring 12 1/2"L x 10 1/2"W, reusable screen insert, and drain pan tubing (8.3 feet long).

Table Compatibility: See Special Note(s)

Special Note(s): Compatible with the 9-inch Urology Extension attached to 5085 SRT, 5085, 4085, and Cmax. Will install directly onto Amsco 3080/3085 tables or IA Extender.



Drain Pan Tubing

#R003500835

This replacement tubing can be ordered by the foot. Indicate the length desired, otherwise, the tubing will come in foot long segments.

Usage/Procedures: Urologic

Special Note(s): Compatible with Drain Pan (BF186) (sold separately).



Transfer Board

#BF588

This lightweight board supports the patient's legs during transfer to table. It is recommended for leg support only and is removed after patient is positioned.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax



Patient Transfer Board with 2" Pad - 3080, 3080-R, 3085 SP

#BF204

Supports patient's legs during transfer to table. Recommended for leg support only. Inserts into IA Extender and is removed after patient is positioned. Includes a 2-inch removable foam pad. Pad attaches to transfer board with Velcro strips.

Table Compatibility: 3080, 3080-R, 3085 SP

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Replacement Pad(s): BF820

Special Note(s): Can also be used with the IA Extender (BF188) or 9" Urology Extension (BF578) attached to 5085SRT, 5085, 4085 and Cmax.



Bariatric Table Extender Pads Only; 4" Foot section - 3080, 3085SP

#BF559

This set of 4" pads for the Bariatric Table Extenders includes 2 pads for the Foot section.

Table Compatibility: 3080, 3085 SP

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: Bariatric; Obstetric Surgical Procedures

Fluid Collection System

#BF457

The Fluid Collection System includes a flexible drape support and fluid collection bag and it is ideal for use in cystoscopic or hysteroscopic procedures. The flexible drape support device quickly attaches to the end of the table seat section during set up, and is easily manipulated intra-procedurally to obtain maximum unobstructed surgical site access. Both the Sterile and Non-sterile Fluid Collector Bags fit either the Standard or Bariatric model drape support. The built-in flexible drain hose works with a kick bucket or a suction canister. Bariatric Fluid Collection System also available (BF556).

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, See Special Note(s)

Special Note(s): Also fits IA Extender and 9-inch Urology Extender. Set of 10 Fluid Collection Bags: Sterile (BF458) or Non-Sterile (BF459) are required (sold separately).



Bariatric Fluid Collection System

#BF56

The Bariatric Fluid Collection System includes adapters that extend to accommodate the additional width created by the Bariatric Table Extenders applied to the table. It is ideal for use in cystoscopic or hysteroscopic procedures, and comes with the flexible drape support and a fluid collection bag. The flexible drape support device quickly attaches to the Seat Section of the Bariatric Table Extender during set up, and is easily manipulated intra-procedurally to obtain maximum unobstructed surgical site access. Both the Sterile and Non-sterile Fluid Collector Bags fit either the standard or Bariatric model drape support. The built-in flexible drain hose works with a kick bucket or a suction canister.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Urologic; Cystoscopic; Hysteroscopic

Special Note(s): Also compatible with 3080 with Bariatric Table Extenders applied. Set of 10 Fluid Collection Bags: Sterile (BF458) or Non-sterile (BF459) are required (sold separately).



Fluid Collector Bags; Sterile Box of 10


#BF458

These sterile disposable bags slide over the flexible drape support. Each bag has a tissue retention filter as well as a flexible, expandable drain hose that can be connected to suction or used with a kick bucket. These bags fit both the Standard and Bariatric drape supports.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Urologic; Cystoscopic; Hysteroscopic

Special Note(s): Fluid Collection System: Standard (BF457) or Bariatric (BF56) is required (sold separately).



Fluid Collector Bags; Non-Sterile Box of 10


#BF459

These non-sterile disposable bags slide over the flexible drape support. Each bag has a tissue retention filter as well as a flexible, expandable drain hose that can be connected to suction or used with a kick bucket. These bags fit both the Standard and Bariatric drape supports.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax

Usage/Procedures: Urologic; Cystoscopic; Hysteroscopic

Special Note(s): Fluid Collection System: Standard (BF457) or Bariatric (BF56) is required (sold separately).



X-ray Top for 9 Inch Urology Extension

#BF579

Removable, radiolucent X-ray top mounts onto the 9 Inch Urology Extension. The 2-inch standoffs add elevation to the extension providing space to insert flat X-ray cassettes. Provides even surface when X-ray tops are mounted to the surgical table.

Table Compatibility: 4085, 5085, 5085 SRT, Cmax

Weight Capacity: Maximum weight capacity is the equivalent to that of the 9 Inch Urology Extension.

Usage/Procedures: Endourologic; KUB; Lithotomy

Required Accessories: BF578

Standard Siderail Adaptor (Each)

#TA00012

The Standard Siderail Adaptors provide attachment of accessories requiring standard siderails. Moveable and removable anywhere along side of table.

Table Compatibility: SurgiGraphic 1027, SurgiGraphic 6000, 36" Carbon Fiber Fluoro Extension



Standard Siderail Adaptor (Set of 3)

#TA00005

The Standard Siderail Adaptors provide attachment of accessories requiring standard siderails. Moveable and removable anywhere along side of table.

Table Compatibility: SurgiGraphic 1027, SurgiGraphic 6000, 36" Carbon Fiber Fluoro Extension



Patient Safety Straps (Pair)

#TA00001

Radiolucent strap wraps around patient and tabletop. Fits all size patients with Velcro securement.

Table Compatibility: SurgiGraphic 1027, SurgiGraphic 6000

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.



Catheter Extension

#TA00002

The Catheter Extension attaches to the foot end of table and moves with the table top during float, serving as a convenient sterile instrument surface.

Table Compatibility: SurgiGraphic 6000

Usage/Procedures: Endovascular Surgery

Special Note(s): Not a patient support.



Radial AV Fistula Board

#TA00003

This larger radiolucent platform is ideal for fistula repair with an uninterrupted fluoro scan for intraoperative flow studies from radial AV fistula to vena cava. Includes pad and pair of Side Rail Locks.

Table Compatibility: SurgiGraphic 1027, SurgiGraphic 6000, 36" Carbon Fiber Fluoro Extension

Weight Capacity: Maximum patient weight capacity is equivalent to that of the table.

Usage/Procedures: AV Fistula



Carbon Fiber Extension with Pad


#TA00006

Provides additional 12-inches to increase the normal 71-inch fluoro body scan to 83-inches of uninterrupted imaging coverage.

Table Compatibility: SurgiGraphic 6000

Weight Capacity: Maximum patient weight capacity is 400 lbs. (181 kg).

Replacement Pad(s): BF495



Side Rail Extender

#BF456

The Side Rail Extender extends the side rail of the SurgiGraphic 6000 table to allow for positioning of an armboard when the Carbon Fiber Extension (TA00006) is attached. This extender attaches to a standard side rail (or on the Standard Siderail Adaptor) and can be used on most surgical tables.

Table Compatibility: 3080, 3080-R, 3085 SP, 4085, 5085, 5085 SRT, Cmax, OrthoVision

Special Note(s): Also fits SurgiGraphic tables with use of Standard Siderail Adaptors (TA00012).



Description	Item	Page No.
3085 Shoulder Table Accessory Package	BF382	12
3085 Uro/Endo Accessory Package	BF402	11
36" Carbon Fiber Fluoro Extension	BF642	75
9 Inch Urology Extension	BF578	77
Adult Arm and Hand Table	BF435	17
A-Frame Accessory Cart	BF27	5
Andrews Frame Adaptor Kit	BF305	62
Anesthesia Armboard with 1" Pad	BF448	14
Anesthesia Armboard with 2" Pad	BF803	14
Anesthesia Armboard with 3" Pad	BF804	14
Anesthesia Armboard without Pad	BF431	14
AQUAGEL Arm Protector, Long	BF476	17, 38
AQUAGEL Arm Protector, Short	BF475	17, 37
AQUAGEL Head Protector	BF463	29, 35
AQUAGEL Heel Protector	BF477	34, 38
AQUAGEL Leg Protector	BF478	34, 38
AQUAGEL Patient Positioner	BF481	22, 36
AQUAGEL Sacral Protector, Small	BF474	22, 36
Arthroscopy Knee Holder	BF717	70
Bariatric Accessory Storage Cart	BF551	5, 48
Bariatric Fluid Collection System	BF56	58, 80
Bariatric Foot Extension	BF542	33, 49
Bariatric Power-lift Stirrups (Pair)	BF454	31, 48, 78
Bariatric Procedure Package with Split Leg	BF852	9
Bariatric Reverse Trendelenburg Strap (Pair)	BF563	19, 48
Bariatric Table Extender Pads Only - 4"; Leg Section - 5085 SRT, 5085, 4085, Cmax	BF659	55
Bariatric Table Extender Pads Only - 2"; Back and Seat Sections - 5085 SRT, 5085, 4085, Cmax	BF654	54
Bariatric Table Extender Pads Only - 2"; Head Section - 5085 SRT, 5085, 4085, Cmax	BF656	54
Bariatric Table Extender Pads Only - 2"; Leg Section - 5085 SRT, 5085, 4085, Cmax	BF658	55
Bariatric Table Extender Pads Only - 4"; Back and Seat Sections - 5085 SRT, 5085, 4085, Cmax	BF655	54
Bariatric Table Extender Pads Only - 4"; Complete Set - 5085 SRT, 5085, 4085, Cmax	BF653	53
Bariatric Table Extender Pads Only - 4"; Head Section - 5085 SRT, 5085, 4085, Cmax	BF657	54
Bariatric Table Extender Pads Only 2"; Complete Set - 5085 SRT, 5085, 4085, Cmax	BF652	53
Bariatric Table Extender Pads Only; 2" Back and Seat sections - 3080, 3085SP	BF544	57
Bariatric Table Extender Pads Only; 2" Complete Set - 3080, 3085SP	BF543	57
Bariatric Table Extender Pads Only; 2" Foot section - 3080, 3085SP	BF546	57
Bariatric Table Extender Pads Only; 2" Head section - 3080, 3085SP	BF545	57
Bariatric Table Extender Pads Only; 4" Back and Seat sections - 3080, 3085SP	BF557	56
Bariatric Table Extender Pads Only; 4" Complete Set - 3080, 3085SP	BF556	56
Bariatric Table Extender Pads Only; 4" Foot section - 3080, 3085SP	BF559	79
Bariatric Table Extender Pads Only; 4" Head section - 3080, 3085SP	BF558	57
Bariatric Table Extenders with 2" Pads; Back and Seat - 3080, 3085	BF553	55
Bariatric Table Extenders with 2" Pads; Back and Seat Sections - 5085SRT, 5085, 4085, Cmax	BF645	51
Bariatric Table Extenders with 2" Pads; Complete Set - 5085SRT, 5085, 4085, Cmax	BF643	50
Bariatric Table Extenders with 2" Pads; Complete Set - 3080, 3085	BF541	55
Bariatric Table Extenders with 2" Pads; Foot Section - 3080, 3085	BF555	56
Bariatric Table Extenders with 2" Pads; Head Section - 3080, 3085	BF554	56
Bariatric Table Extenders with 2" Pads; Head Section - 5085SRT, 5085, 4085, Cmax	BF647	52
Bariatric Table Extenders with 2" Pads; Leg Section - 5085SRT, 5085, 4085, Cmax	BF649	53
Bariatric Table Extenders with 4" Pads; Back and Seat Sections - 5085SRT, 5085, 4085, Cmax	BF646	51
Bariatric Table Extenders with 4" Pads; Complete Set - 5085SRT, 5085, 4085, Cmax	BF644	50
Bariatric Table Extenders with 4" Pads; Head Section - 5085SRT, 5085, 4085, Cmax	BF648	52
Bariatric Table Extenders with 4" Pads; Leg Section - 5085SRT, 5085, 4085, Cmax	BF651	53

Description	Item	Page No.
Bariatric Table Extenders without Pads; Back and Seat Sections - 5085SRT, 5085, 4085, Cmax	BF662	50
Bariatric Table Extenders without Pads; Complete Set - 5085SRT, 5085, 4085, Cmax	BF661	49
Bariatric Table Extenders without Pads; Head Section - 5085SRT, 5085, 4085, Cmax	BF663	51
Bariatric Table Extenders without Pads; Leg Section - 5085SRT, 5085, 4085, Cmax	BF664	52
Bariatric Urology Package	BF853	11
Cable Wrap	BF316	43
Carbon Fiber Extension with Pad	TA00006	82
Carbon Fiber Hip Rest with Pad	BF00072	71
Carbon Fiber Tibia/Lateral Countertraction Support	BF00062	70
Carbon Fiber T-Shaped Sacral Rest with Pad	BF00073	70
Carbon Fiber Wedge-Shaped Sacral Rest with Pad - Adult	BF00066	71
Cardiovascular Accessory Package	BF412	10
Cassette Tray with Handle - 14" x 17" Film	BF137	45, 76
Catheter Extension	TA00002	81
Clark Sockets (Pair)	BF083	24
Cone-Shaped Headrest	BF215	27, 59, 63
Drain Pan	BF186	78
Drain Pan Tubing	R003500835	26, 79
ELITE Flat Bottom Chest Roll, Large	BF483	23, 37
ELITE Flat Bottom Chest Roll, Medium	BF482	22, 37
ELITE Lateral Positioner	BF484	23, 37
ELITE Prone Headrest, Large	BF464	29, 35
ELITE Prone Headrest, Medium	BF465	29, 36
ELITE Prone Positioner	BF485	23, 37
ELITE Supine Headrest, Large	BF466	30, 36
ELITE Supine Headrest, Medium	BF467	30, 36
Eye-ENT-Neuro Headrest Adaptor - 3080, 3080-R, 3085	BF185	27, 59, 63
Eye-ENT-Neuro Headrest Adaptor - 5085 SRT	BF309	28, 60, 64
Eye-ENT-Neuro Headrest Adaptor - 5085, 4085, Cmax	BF173	27, 59, 63
Featherweight Leg Section - 5085 SRT	BF012	31
Featherweight Leg Section - 5085, 4085, Cmax	BF009	31
Fem / Pop Board	BF287	76
Fluid Collection System	BF457	79
Fluid Collector Bags; Non-Sterile Box of 10	BF459	26, 80
Fluid Collector Bags; Sterile Box of 10	BF458	26, 80
Foot Control - 3080 SP, 3080-R and 3080	BF187	43
Foot Control - 5085 SRT, 5085, 4085 and Cmax	BF586	42
Foot Control - SurgiGraphic 6000	TA00004	42
Foot Extension with Siderail Locks	BF2	33
Foot Rest; Left	MTPA09	32
Foot Rest; Right	MTPA10	32
Foot Traction Boot Pads with Straps - Adult	BF538	68
Foot Traction Boot Pads with Straps - Child	BF539	68
Foot Traction Boot Pads with Straps - Obese Adult	BF548	68
Foot Traction Boots - Obese Adult (Pair)	BF0008	73
Foot Traction Boots (Old Style) - Adult (Pair)	BF00023	72
Foot Traction Boots (Old Style) - Child (Pair)	BF00032	72
Frontal Sitting Attachment	BF01	28, 60
Gastric Bypass Package	BF851	9
GEL Lateral Package, Large	BF488	22, 35
GEL Lateral Package, Medium	BF489	22, 35
GEL Supine Package, Large	BF486	21, 35
GEL Supine Package, Medium	BF487	22, 35

Description	Item	Page No.
Hand Control - 3080 RC	BF201	41
Hand Control - 3085 SP, 3080 RL/SP	BF336	41
Hand Control - Amsco 3085 SP ACT-enabled	BF374	41
Hand Control - Cmax	BF61	40
Hand Control - SurgiGraphic 1027	BF346	42
Hand Control - SurgiGraphic 6000 (old style)	BF33	41
Hand Control - SurgiGraphic 6000 CE	BF35	41
Hand Control (Advanced) - 5085, 5085 SRT	BF329	40
Hand Control (Universal) - 5085, 4085, Cmax	BF319	40
Hand Control Bracket	BF318	43
Head Restraint Straps	BF039	66
Hip Spica Assembly - Adult	BF00031	67
Hip Spica Assembly - Child	BF00038	67
Hip Spica Assembly for Ortho Extension - Child	BF177	67
Horizontal Attachment #1800	BF008	28, 60
Horseshoe Headrest	BF036	28, 60
IA Extender with 2" Pad and X-ray Top	BF188	76, 77
IV Armboard with 2" Pad	BF363	15
IV Armboard with 3" Pad	BF364	15
IV Armboard without Pad	BF141	15
Knee Restraint Strap (Pair)	BF552	20, 48
Krause Arm Support	BF1	16
Lateral Brace Assemblies (Pair)	BF602	21
Lateral Safety Strap	BF562	20
Lithotomy Legholders - Adult (Pair)	BF054	32, 78
L-Shaped Traction Extension	BF0003	67
Micro-Surgical Arm and Hand Table	BF436	17
Morgan Pelvic Rest	BF135	61
Multi-Accessory Clamps (Pair)	BF43	24
Multi-axis Armboard, RC	BF71239	15
Multi-Poise Headrest	BF038	29, 61
Multi-Position Legholders (Pair)	BF433	33, 78
Multi-posture Armboard with 2-inch Pad	BF432	16
Narrow Headrest	BF592	27, 59, 63
OASIS Armboard, Long	BF427	18, 38
OASIS Armboard, Short	BF426	18, 38
OASIS Crutch Stirrup	BF48	34, 39
Orthopedic Extension Accessory Package	BF806	12, 65
Orthopedic Extension Cart	BF24	5
OrthoVision & Orthopedic Extension Pediatric Accessory Package	BF809	13
OrthoVision & Orthopedic Extension Supplemental Accessory Package	BF808	12
Palm Control - SurgiGraphic 6000 (old style)	BF347	42
Palm Control - Surgigraphic 6000 CE	BF348	42
Patient Safety Rail with Cover	BF092	19
Patient Safety Straps (Pair)	TA00001	20, 81
Patient Transfer Board with 2" Pad - 3080, 3080-R, 3085 SP	BF204	79
Patient Transfer Board with 2-inch Pad - OrthoVision	BF213	66
Patient Transfer Board with 4-inch Pad for Tables with Orthopedic Extension	BF214	66
Pediatric Arm and Hand Table	BF434	16
Perineal Cutout Filler Piece	BF205	21
Perineal Post Pad - Disposable	NO001	26, 72
Perineal Post with Pad - Adult	BF00065	71
Perineal Post without Pad - Child	BF41	71
Posterior Sitting Attachment #1801	BF011	29, 61
Power-lift / Bariatric Power-lift Stirrup Storage Cart	BF46	5
Power-lift Stirrups (Pair)	BF453	31, 77
Radial AV Fistula Board	TA00003	81

Description	Item	Page No.
Rail Clamp (Each)	BF462	24
Restraint Strap; Extra-long 126"	BF451	19
Restraint Strap; Standard Length 66-1/2"	BF45	19
Shoulder Chair	BF594	65
Shoulder Chair Package	BF571	11
Shoulder Table	BF437	66
Shoulder Table Accessory Cart	BF461	6
Shoulder Table Chin Cups	BF452	66
Side Rail Extender	BF456	82
Siderail Locks (Pair)	BF133	24
Solid Lateral Braces (Pair)	BF331	21
Split Leg Package	BF537	10
Split Leg Positioner	BF351	32, 49
Standard Accessory Package	BF807	13, 73
Standard Accessory Package with Foot Extension, 2-inch Armboard (to be used without X-ray tops), and Extra-long Restraint Strap 126"	BF567	7
Standard Accessory Package with Foot Extension, 2-inch Armboard (to be used without X-ray tops), and Standard Length Restraint Strap 66-1/2"	BF566	7
Standard Accessory Package with Foot Extension, 3-inch Armboard (to be used with X-ray tops), and Extra-long Restraint Strap 126"	BF565	7
Standard Accessory Package with Foot Extension, 3-inch Armboard (to be used with X-ray tops), and Standard Length Restraint Strap 66-1/2"	BF564	7
Standard Accessory Package with Foot Extension, and 2-inch Armboard Pad (to be used without X-ray tops)	BF099	9
Standard Accessory Package with Foot Extension, and 3-inch Armboard (to be used with X-ray tops)	BF093	9
Standard Accessory Package, 2-inch Armboard (to be used without X-ray tops), and Extra-long Restraint Strap 126"	BF425	8
Standard Accessory Package, 2-inch Armboard (to be used without X-ray tops), and Standard Length Restraint Strap 66-1/2"	BF424	8
Standard Accessory Package, 3-inch Armboard (to be used with X-ray tops), and Extra-long Restraint Strap 126"	BF414	8
Standard Accessory Package, 3-inch Armboard (to be used with X-ray tops), and Standard Length Restraint Strap 66-1/2"	BF413	8
Standard Reverse Trendelenburg Strap (Pair)	BF455	19
Standard Siderail Adaptor (Each)	TA00012	25, 81
Standard Siderail Adaptor (Set of 3)	TA00005	25, 81
Straight Traction Extension	BF268	67
Surgical Drape Support	BF00018	72
SurgiGraphic 1027 & 6000 Standard Accessory Package	TA005	13
SurgiGraphic 1027 Vascular Accessory Package	TA00502	13
SurgiGraphic 6000 Vascular Accessory Package	TA00501	13
Table Accessory Cart	BF7	5
TenZor™ Carbon Fiber Traction Unit Assembly	BF72	68
Traction Accessory Clamp	BF0002	70
Traction Block, Skeletal	BF00028	69
Traction Bow, Skeletal	BF724	69
Traction Bow/Block	BF371	69
Traction Bow/Block, Skeletal	BF725	68
Traction Unit Assembly	BF0006	69
Transducer Holder	BF263	47
Transfer Board	BF588	79
Trendelenburg Accessory Package	BF607	10
TrenStop	BF728	21
Universal Ether Screen with Clark Socket	BF074	47
Universal Legholder (Each)	BF266	33, 65
Uro/Cysto Accessory Package	BF574	10
Uro/Endo Accessory Package with Drain Pan	BF609	11

Description	Item	Page No.
Uro/Endo Accessory Package with Fluid Collection System	BF576	10
Wedge-Shaped Sacral Rest with Pad - Child	BF00048	71
Winged Ether Screen #1401	BF075	47
Wrist Rest for Narrow Headrest	BF593	30, 62, 64
Wristlets with Straps and Brackets (Pair)	BF111	20
X-ray Top for 9 Inch Urology Extension	BF579	46, 80
X-ray Tops Back Section - 5085 SRT, 5085, 4085, Cmax	BF175	44, 74
X-ray Tops Complete Set - 3080, 3080-R, 3085 SP	BF202	45, 75
X-ray Tops Complete Set - 5085SRT, 5085, 4085, Cmax	BF219	44, 74
X-ray Tops Head Section - 5085 SRT, 5085, 4085, Cmax	BF174	44, 74
X-ray Tops Leg Section - 5085 SRT, 5085, 4085, Cmax	BF744	45, 75
X-ray Tops Seat Section - 5085 SRT, 5085, 4085, Cmax	BF176	45, 75

Item	Description	Page No.
BF00018	Surgical Drape Support	72
BF0002	Traction Accessory Clamp	70
BF00023	Foot Traction Boots (Old Style) - Adult (Pair)	72
BF00028	Traction Block, Skeletal	69
BF0003	L-Shaped Traction Extension	67
BF00031	Hip Spica Assembly - Adult	67
BF00032	Foot Traction Boots (Old Style) - Child (Pair)	72
BF00038	Hip Spica Assembly - Child	67
BF00048	Wedge-Shaped Sacral Rest with Pad - Child	71
BF0006	Traction Unit Assembly	69
BF00062	Carbon Fiber Tibia/Lateral Countertraction Support	70
BF00065	Perineal Post with Pad - Adult	71
BF00066	Carbon Fiber Wedge-Shaped Sacral Rest with Pad - Adult	71
BF00072	Carbon Fiber Hip Rest with Pad	71
BF00073	Carbon Fiber T-Shaped Sacral Rest with Pad	70
BF0008	Foot Traction Boots - Obese Adult (Pair)	73
BF008	Horizontal Attachment #1800	28, 60
BF009	Featherweight Leg Section - 5085, 4085, Cmax	31
BF01	Frontal Sitting Attachment	28, 60
BF011	Posterior Sitting Attachment #1801	29, 61
BF012	Featherweight Leg Section - 5085 SRT	31
BF036	Horseshoe Headrest	28, 60
BF038	Multi-Poise Headrest	29, 61
BF039	Head Restraint Straps	66
BF054	Lithotomy Legholders - Adult (Pair)	32, 78
BF074	Universal Ether Screen with Clark Socket	47
BF075	Winged Ether Screen #1401	47
BF083	Clark Sockets (Pair)	24
BF092	Patient Safety Rail with Cover	19
BF093	Standard Accessory Package with Foot Extension, and 3-inch Armboard (to be used with X-ray tops)	9
BF099	Standard Accessory Package with Foot Extension, and 2-inch Armboard Pad (to be used without X-ray tops)	9
BF1	Krause Arm Support	16
BF111	Wristlets with Straps and Brackets (Pair)	20
BF133	Siderail Locks (Pair)	24
BF135	Morgan Pelvic Rest	61
BF137	Cassette Tray with Handle - 14" x 17" Film	45, 76
BF141	IV Armboard without Pad	15
BF173	Eye-ENT-Neuro Headrest Adaptor - 5085, 4085, Cmax	27, 59, 63
BF174	X-ray Tops Head Section - 5085 SRT, 5085, 4085, Cmax	44, 74
BF175	X-ray Tops Back Section - 5085 SRT, 5085, 4085, Cmax	44, 74
BF176	X-ray Tops Seat Section - 5085 SRT, 5085, 4085, Cmax	45, 75
BF177	Hip Spica Assembly for Ortho Extension - Child	67
BF185	Eye-ENT-Neuro Headrest Adaptor - 3080, 3080-R, 3085	27, 59, 63
BF186	Drain Pan	78
BF187	Foot Control - 3080 SP, 3080-R and 3080	43
BF188	IA Extender with 2" Pad and X-ray Top	76, 77
BF2	Foot Extension with Siderail Locks	33
BF201	Hand Control - 3080 RC	41
BF202	X-ray Tops Complete Set - 3080, 3080-R, 3085 SP	45, 75
BF204	Patient Transfer Board with 2" Pad - 3080, 3080-R, 3085 SP	79
BF205	Perineal Cutout Filler Piece	21
BF213	Patient Transfer Board with 2-inch Pad - OrthoVision	66
BF214	Patient Transfer Board with 4-inch Pad for Tables with Orthopedic Extension	66

Item	Description	Page No.
BF215	Cone-Shaped Headrest	27, 59, 63
BF219	X-ray Tops Complete Set - 5085SRT, 5085, 4085, Cmax	44, 74
BF24	Orthopedic Extension Cart	5
BF263	Transducer Holder	47
BF266	Universal Legholder (Each)	33, 65
BF268	Straight Traction Extension	67
BF27	A-Frame Accessory Cart	5
BF287	Fem / Pop Board	76
BF305	Andrews Frame Adaptor Kit	62
BF309	Eye-ENT-Neuro Headrest Adaptor - 5085 SRT	28, 60, 64
BF316	Cable Wrap	43
BF318	Hand Control Bracket	43
BF319	Hand Control (Universal) - 5085, 4085, Cmax	40
BF329	Hand Control (Advanced) - 5085, 5085 SRT	40
BF33	Hand Control - SurgiGraphic 6000 (old style)	41
BF331	Solid Lateral Braces (Pair)	21
BF336	Hand Control - 3085 SP, 3080 RL/SP	41
BF346	Hand Control - SurgiGraphic 1027	42
BF347	Palm Control - SurgiGraphic 6000 (old style)	42
BF348	Palm Control - Surgigraphic 6000 CE	42
BF35	Hand Control - SurgiGraphic 6000 CE	41
BF351	Split Leg Positioner	32, 49
BF363	IV Armboard with 2" Pad	15
BF364	IV Armboard with 3" Pad	15
BF371	Traction Bow/Block	69
BF374	Hand Control - Amsco 3085 SP ACT-enabled	41
BF382	3085 Shoulder Table Accessory Package	12
BF402	3085 Uro/Endo Accessory Package	11
BF41	Perineal Post without Pad - Child	71
BF412	Cardiovascular Accessory Package	10
BF413	Standard Accessory Package, 3-inch Armboard (to be used with X-ray tops), and Standard Length Restraint Strap 66-1/2"	8
BF414	Standard Accessory Package, 3-inch Armboard (to be used with X-ray tops), and Extra-long Restraint Strap 126"	8
BF424	Standard Accessory Package, 2-inch Armboard (to be used without X-ray tops), and Standard Length Restraint Strap 66-1/2"	8
BF425	Standard Accessory Package, 2-inch Armboard (to be used without X-ray tops), and Extra-long Restraint Strap 126"	8
BF426	OASIS Armboard, Short	18, 38
BF427	OASIS Armboard, Long	18, 38
BF43	Multi-Accessory Clamps (Pair)	24
BF431	Anesthesia Armboard without Pad	14
BF432	Multi-posture Armboard with 2-inch Pad	16
BF433	Multi-Position Legholders (Pair)	33, 78
BF434	Pediatric Arm and Hand Table	16
BF435	Adult Arm and Hand Table	17
BF436	Micro-Surgical Arm and Hand Table	17
BF437	Shoulder Table	66
BF448	Anesthesia Armboard with 1" Pad	14
BF45	Restraint Strap; Standard Length 66-1/2"	19
BF451	Restraint Strap; Extra-long 126"	19
BF452	Shoulder Table Chin Cups	66
BF453	Power-lift Stirrups (Pair)	31, 77
BF454	Bariatric Power-lift Stirrups (Pair)	31, 48, 78
BF455	Standard Reverse Trendelenburg Strap (Pair)	19

Item	Description	Page No.
BF456	Side Rail Extender	82
BF457	Fluid Collection System	79
BF458	Fluid Collector Bags; Sterile Box of 10	26, 80
BF459	Fluid Collector Bags; Non-Sterile Box of 10	26, 80
BF46	Power-lift / Bariatric Power-lift Stirrup Storage Cart	5
BF461	Shoulder Table Accessory Cart	6
BF462	Rail Clamp (Each)	24
BF463	AQUAGEL Head Protector	29, 35
BF464	ELITE Prone Headrest, Large	29, 35
BF465	ELITE Prone Headrest, Medium	29, 36
BF466	ELITE Supine Headrest, Large	30, 36
BF467	ELITE Supine Headrest, Medium	30, 36
BF474	AQUAGEL Sacral Protector, Small	22, 36
BF475	AQUAGEL Arm Protector, Short	17, 37
BF476	AQUAGEL Arm Protector, Long	17, 38
BF477	AQUAGEL Heel Protector	34, 38
BF478	AQUAGEL Leg Protector	34, 38
BF48	OASIS Crutch Stirrup	34, 39
BF481	AQUAGEL Patient Positioner	22, 36
BF482	ELITE Flat Bottom Chest Roll, Medium	22, 37
BF483	ELITE Flat Bottom Chest Roll, Large	23, 37
BF484	ELITE Lateral Positioner	23, 37
BF485	ELITE Prone Positioner	23, 37
BF486	GEL Supine Package, Large	21, 35
BF487	GEL Supine Package, Medium	22, 35
BF488	GEL Lateral Package, Large	22, 35
BF489	GEL Lateral Package, Medium	22, 35
BF537	Split Leg Package	10
BF538	Foot Traction Boot Pads with Straps - Adult	68
BF539	Foot Traction Boot Pads with Straps - Child	68
BF541	Bariatric Table Extenders with 2" Pads; Complete Set - 3080, 3085	55
BF542	Bariatric Foot Extension	33, 49
BF543	Bariatric Table Extender Pads Only; 2" Complete Set - 3080, 3085SP	57
BF544	Bariatric Table Extender Pads Only; 2" Back and Seat sections - 3080, 3085SP	57
BF545	Bariatric Table Extender Pads Only; 2" Head section - 3080, 3085SP	57
BF546	Bariatric Table Extender Pads Only; 2" Foot section - 3080, 3085SP	57
BF548	Foot Traction Boot Pads with Straps - Obese Adult	68
BF551	Bariatric Accessory Storage Cart	5, 48
BF552	Knee Restraint Strap (Pair)	20, 48
BF553	Bariatric Table Extenders with 2" Pads; Back and Seat - 3080, 3085	55
BF554	Bariatric Table Extenders with 2" Pads; Head Section - 3080, 3085	56
BF555	Bariatric Table Extenders with 2" Pads; Foot Section - 3080, 3085	56
BF556	Bariatric Table Extender Pads Only; 4" Complete Set - 3080, 3085SP	56
BF557	Bariatric Table Extender Pads Only; 4" Back and Seat sections - 3080, 3085SP	56
BF558	Bariatric Table Extender Pads Only; 4" Head section - 3080, 3085SP	57
BF559	Bariatric Table Extender Pads Only; 4" Foot section - 3080, 3085SP	79
BF56	Bariatric Fluid Collection System	58, 80
BF562	Lateral Safety Strap	20
BF563	Bariatric Reverse Trendelenburg Strap (Pair)	19, 48
BF564	Standard Accessory Package with Foot Extension, 3-inch Armboard (to be used with X-ray tops), and Standard Length Restraint Strap 66-1/2"	7
BF565	Standard Accessory Package with Foot Extension, 3-inch Armboard (to be used with X-ray tops), and Extra-long Restraint Strap 126"	7
BF566	Standard Accessory Package with Foot Extension, 2-inch Armboard (to be used without X-ray tops), and Standard Length Restraint Strap 66-1/2"	7
BF567	Standard Accessory Package with Foot Extension, 2-inch Armboard (to be used without X-ray tops), and Extra-long Restraint Strap 126"	7
BF571	Shoulder Chair Package	11

Item	Description	Page No.
BF574	Uro/Cysto Accessory Package	10
BF576	Uro/Endo Accessory Package with Fluid Collection System	10
BF578	9 Inch Urology Extension	77
BF579	X-ray Top for 9 Inch Urology Extension	46, 80
BF586	Foot Control - 5085 SRT, 5085, 4085 and Cmax	42
BF588	Transfer Board	79
BF592	Narrow Headrest	27, 59, 63
BF593	Wrist Rest for Narrow Headrest	30, 62, 64
BF594	Shoulder Chair	65
BF602	Lateral Brace Assemblies (Pair)	21
BF607	Trendelenburg Accessory Package	10
BF609	Uro/Endo Accessory Package with Drain Pan	11
BF61	Hand Control - Cmax	40
BF642	36" Carbon Fiber Fluoro Extension	75
BF643	Bariatric Table Extenders with 2" Pads; Complete Set - 5085SRT, 5085, 4085, Cmax	50
BF644	Bariatric Table Extenders with 4" Pads; Complete Set - 5085SRT, 5085, 4085, Cmax	50
BF645	Bariatric Table Extenders with 2" Pads; Back and Seat Sections - 5085SRT, 5085, 4085, Cmax	51
BF646	Bariatric Table Extenders with 4" Pads; Back and Seat Sections - 5085SRT, 5085, 4085, Cmax	51
BF647	Bariatric Table Extenders with 2" Pads; Head Section - 5085SRT, 5085, 4085, Cmax	52
BF648	Bariatric Table Extenders with 4" Pads; Head Section - 5085SRT, 5085, 4085, Cmax	52
BF649	Bariatric Table Extenders with 2" Pads; Leg Section - 5085SRT, 5085, 4085, Cmax	53
BF651	Bariatric Table Extenders with 4" Pads; Leg Section - 5085SRT, 5085, 4085, Cmax	53
BF652	Bariatric Table Extender Pads Only 2"; Complete Set - 5085 SRT, 5085, 4085, Cmax	53
BF653	Bariatric Table Extender Pads Only - 4"; Complete Set - 5085 SRT, 5085, 4085, Cmax	53
BF654	Bariatric Table Extender Pads Only - 2"; Back and Seat Sections - 5085 SRT, 5085, 4085, Cmax	54
BF655	Bariatric Table Extender Pads Only - 4"; Back and Seat Sections - 5085 SRT, 5085, 4085, Cmax	54
BF656	Bariatric Table Extender Pads Only - 2"; Head Section - 5085 SRT, 5085, 4085, Cmax	54
BF657	Bariatric Table Extender Pads Only - 4"; Head Section - 5085 SRT, 5085, 4085, Cmax	54
BF658	Bariatric Table Extender Pads Only - 2"; Leg Section - 5085 SRT, 5085, 4085, Cmax	55
BF659	Bariatric Table Extender Pads Only - 4"; Leg Section - 5085 SRT, 5085, 4085, Cmax	55
BF661	Bariatric Table Extenders without Pads; Complete Set - 5085SRT, 5085, 4085, Cmax	49
BF662	Bariatric Table Extenders without Pads; Back and Seat Sections - 5085SRT, 5085, 4085, Cmax	50
BF663	Bariatric Table Extenders without Pads; Head Section - 5085SRT, 5085, 4085, Cmax	51
BF664	Bariatric Table Extenders without Pads; Leg Section - 5085SRT, 5085, 4085, Cmax	52
BF7	Table Accessory Cart	5
BF71239	Multi-axis Armboard, RC	15
BF717	Arthroscopy Knee Holder	70
BF72	TenZor™ Carbon Fiber Traction Unit Assembly	68
BF724	Traction Bow, Skeletal	69
BF725	Traction Bow/Block, Skeletal	68
BF728	TrenStop	21
BF744	X-ray Tops Leg Section - 5085 SRT, 5085, 4085, Cmax	45, 75
BF803	Anesthesia Armboard with 2" Pad	14
BF804	Anesthesia Armboard with 3" Pad	14
BF806	Orthopedic Extension Accessory Package	12, 65
BF807	Standard Accessory Package	13, 73
BF808	OrthoVision & Orthopedic Extension Supplemental Accessory Package	12
BF809	OrthoVision & Orthopedic Extension Pediatric Accessory Package	13
BF851	Gastric Bypass Package	9
BF852	Bariatric Procedure Package with Split Leg	9
BF853	Bariatric Urology Package	11

Item	Description	Page No.
MTPA09	Foot Rest; Left	32
MTPA10	Foot Rest; Right	32
NO001	Perineal Post Pad - Disposable	26, 72
R003500835	Drain Pan Tubing	26, 79
TA00001	Patient Safety Straps (Pair)	20, 81
TA00002	Catheter Extension	81
TA00003	Radial AV Fistula Board	81
TA00004	Foot Control - SurgiGraphic 6000	42
TA00005	Standard Siderail Adaptor (Set of 3)	25, 81
TA00006	Carbon Fiber Extension with Pad	82
TA00012	Standard Siderail Adaptor (Each)	25, 81
TA005	SurgiGraphic 1027 & 6000 Standard Accessory Package	13
TA00501	SurgiGraphic 6000 Vascular Accessory Package	13
TA00502	SurgiGraphic 1027 Vascular Accessory Package	13

HERMES is a registered trademark of Intuitive
OR₁ is a registered trademark of Karl Storz Endoscopy, Inc.
*Velcro is a registered trademark of Velcro Corporation.



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 · USA
440-354-2600 · 800-548-4873
www.steris.com