

Valisafe **ScopeCheck**

& Valisafe **InstruCheck**

Medisafe[®]

Valisafe



The Valisafe **ScopeCheck** and Valisafe **InstruCheck** protein tests are a rapid test that have the ability to detect residual proteins left behind on the surfaces of washer-disinfectors, ultrasonic cleaners, endoscopes and other hard to clean surgical instruments. This protein test is based on a dye-binding solution used in clinical chemistry and can detect protein residues within 1 µg sensitivity.



Unexposed

0ug

1ug

2ug

5ug

10ug

Positive
control

- Conforms to BS EN ISO 15883
- Result within 10 seconds
- Easy to use - no incubation required
- Clear colour change
- Cost effective
- For use with:
 - Endoscopes
 - WD / Ultrasonic surfaces
 - Surgical instruments

Order#

V3504095 – Scope-Check / Instru-Check Standard Kit – 25 vials & 25 x standard swabs

V3505005 – Scope-Check 2.8mm – 25 vials & 25 x 2.8mm swabs

V3505010 – Scope-Check 1.9mm – 25 vials & 25 x 1.9mm swabs

V3505013 – Scope-Check 3.7mm – 25 vials & 25 x 3.7mm swabs

Valisafe is a division of Medisafe International

SafMed

Johannesburg Office:

Tel: 011 201 4300

Fax: 011 507 5550

Durban Office:

Tel: 082 464 3832

Fax: 031 561 5138

Cape Town Office:

Tel: 021 763 3280

Fax: 021 762 2760

E-mail: info@safmed.co.za

Website: www.safmed.co.za

www.facebook.com/Safmed.SA/

Twitter: @SafMedSA