

SafMed Tyvek Sterilization Reels



PMS Tyvek sterilization reels and pouches for use in H₂O₂ sterilization are constructed from transparent multilayer PET/PE copolymer film and uncoated Tyvek web in compliance with EN 868-9. A water based, non-toxic process indicator for hydrogen peroxide sterilization complying with ISO 11140-1 is applied on the Tyvek surface and helps to differentiate between processed and unprocessed packages.

Features:

- Excellent microbial barrier properties
- Easy peel and aseptic presentation
- Wide sealing temperature window
- Easy seal with old sealing machines
- Clear and accurate color change

A wide range of standard sizes allow optimum choice of correct sized packaging of each item. Go to www.safmed.co.za/tyvek-reels.html to see the full range of sizes.