

SafMed Sterilization Reels



PMSSteripack sterilization reels and pouches are constructed from transparent PET/PP multi-layer copolymer film and medical grade paper. Water based, non-toxic process indicators for steam and ethylene oxide sterilization complying with ISO 11140-1 are applied on the paper surface and help to differentiate between processed and unprocessed packages.

Features:

- Superior barrier with 60gsm or 70gsm medical grade paper conforming EN 868-3
- Transparent, reinforced multilayer co-polymer PET/PP film
- Three independent seal lines
- Excellent gas permeability and microbial barrier properties
- Can be easily sealed with all sealing machines
- Manufactured under validated production process

Flat and gusseted type and wide range of standard sizes allow optimum choice of correct sized packaging of each item. [Click here to view the available sizes.](#)