

SafMed Reliance 440 SMMMS

Tear and burst resistance combined with antistatic treatment, low memory and low linting makes SafMed Reliance 440 one of the safest protective wraps.

The key to the barrier protection of SafMed Reliance 440 SMMMS is the three Meltblown layers of microscopic fibers, delivering a tortuous path for a consistent bacteriological barrier to maintain sterility.

SafMed Reliance 440 conforms to global and European sterile barrier system standards for terminally sterilized medical devices ISO 11607-1, EN 868-2, SANS 11607-1:2020 and SANS 11607-2:2020.

It is available in various sizes to satisfy different user requirements.

[Click here to see the various sizes available.](#)

Features & Benefits

- Excellent bacterial barrier
- Liquid repellent and antistatic
- High strength to avoid rips, tears or punctures
- Low linting



STEAM



EO



FORM



H2O2

