

SafMed

Surgical Gel Pads & Positioners



SafMed

Your solutions partner for **30 YEARS**

We specialise in providing solutions and products that improve the productivity and effectiveness in the CSSD, Operating Theatre & Endoscopy Unit. From expert advice on new theatre or CSSD designs to providing the latest digital operating theatre and instrument tracking technology, we have the products and experience to offer meaningful value to our customers.



Process Improvement



Track, Trace & Asset Management



Digital Theatre Systems



CSSD & Operating Theatre Design Advice



Online & Classroom Based Education



Sterility Assurance



Decontamination



Customized Service Programs

**Gauteng, Mpumalanga,
Limpopo & North West**
Tel: 011 201 4300

Kwa-Zulu Natal
Tel: 082 464 3832

Western & Eastern Cape
Tel: 021 763 3280

Free State & Northern Cape
Tel: 082 686 8075

E-mail: info@safmed.co.za
Website: www.safmed.co.za
www.facebook.com/Safmed.SA/
Twitter: @SafMedSA

Blue Diamond® Gel Pads and Positioners

The Blue Diamond® range of Gel Pads and Positioners were designed to help reduce patient pressure injuries.

FEATURES:

- 100% Latex-Free
- Radiolucent
- Reusable
- Contains no water, silicones or plasticisers that tend to degrade over time
- The urethane film is formulated to provide optimal elongation and wear and tear.
- A ¼" continuous seal to the outer edge of the product to eliminate the chance for contaminants to collect between the skins
- Washable with everyday theatre cleaning products

GEL HEAD DONUT and HORSESHOE HEAD DONUT			
	<p>NEONATAL GEL HORSESHOE DONUT</p> <p>Use: Neonates and Premies</p> <p>Designed to support the head in prone positions</p> <p>Procedures: ENT, Ophthalmic, plastics and general surgery</p> <p>Dimensions Height: 1.9 cm Diameter: 8.25 cm Hollow Keyhole Centre</p> <p><i>Reorder number: BD2165</i></p>		<p>NEONATAL HEAD DONUT</p> <p>Use: Neonates and Premies</p> <p>Designed to support the head in prone positions</p> <p>Procedures: ENT, Ophthalmic, plastics and general surgery</p> <p>Dimensions Height: 1.9 cm Diameter: 8.25 cm Hollow Centre</p> <p><i>Reorder number: BD2155</i></p>
	<p>PEADIATRIC GEL HORSESHOE DONUT</p> <p>Designed to support the head in prone positions</p> <p>Procedures: ENT, Ophthalmic, plastics and general surgery</p> <p>Dimensions Height: 3 cm Diameter: 14 cm Hollow Keyhole Centre</p> <p><i>Reorder number: BD2175</i></p>		<p>PAEDIACTRIC HEAD DONUT</p> <p>Designed to support the head in prone positions</p> <p>Procedures: ENT, Ophthalmic, plastics and general surgery</p> <p>Dimensions Height: 3 cm Diameter: 15 cm Hollow Centre</p> <p><i>Reorder number: BD2160</i></p>
	<p>ADULT GEL HORSESHOE DONUT</p> <p>Designed to support the head in prone positions</p> <p>Procedures: ENT, Ophthalmic, plastics and general surgery</p> <p>Dimensions Height: 5 cm Diameter: 20 cm Hollow Keyhole Centre</p> <p><i>Reorder number: BD2170</i></p>		<p>ADULT HEAD DONUT</p> <p>Designed to support the head in prone positions</p> <p>Procedures: ENT, Ophthalmic, plastics and general surgery</p> <p>Dimensions Height: 5 cm Diameter: 20 cm Hollow Centre</p> <p><i>Reorder number: BD2150</i></p>










GEL PRONE POSITIONER			
	<p>PRONE HEADREST / RICHARDS HEADREST</p> <p>Gel Prone Headrest with rigid bottom Designed for the head while in prone position Designed with an open centre</p> <p>Dimensions: Standard: 23 cm L x 20 cm W x 7.5 cm H Large: 23 cm L x 20 cm W x 11.5 cm H</p> <p><i>Reorder number:</i> Standard: <i>BD2191</i> Large: <i>BD2192</i></p>		<p>PRONE HEADREST for Anaesthesia with Rigid Bottom</p> <p>Gel Prone Headrest with rigid bottom Designed for the head while in prone position Designed to provide easy access for ET Tubes</p> <p>Dimensions: 28 cm L x 23 cm W x 23 cm H</p> <p><i>Reorder number:</i> <i>BD2225</i></p>
EXTREMITY PROTECTORS			
	<p>CLAM SHELL HEEL PROTECTOR</p> <p>Gel fully encapsulates heel</p> <p>Adjustable Velcro® hook and loop straps</p> <p>Dimensions: Heel to toe: 32 cm Heel to ankle: 17 cm Thickness: 1 cm</p> <p><i>Reorder number:</i> <i>BD2290 each</i></p>		<p>OVAL ULNAR NERVE PROTECTOR</p> <p>Ideal for arm and elbow Oval Ulnar Protector</p> <p>Adjustable Velcro® hook and loop straps</p> <p>Dimensions: 33 cm x 18 cm x 12cm</p> <p><i>Reorder number:</i> <i>BD2270</i></p>
	<p>HEEL CUPS (PAIR)</p> <p>Heel cup that supports heels Notch for the heel to be cradled Does not elevate the leg</p> <p>Dimensions: 15 cm W x 7.5 cm L x 4 cm H</p> <p><i>Reorder number:</i> <i>BD4502P (pair)</i></p>		

TABLE GEL PADS

	<p>TABLE PAD</p> <p>Designed to evenly distribute patients' weight Gel pad can be rolled as a positioning aid Sheer reduction Gel pad for large surfaces</p> <p>Dimensions:</p> <p>¾ length: 51 cm x 117 cm x 1 cm thick Full length: 51 cm x 183 cm x 1 cm thick</p> <p><i>Reorder number:</i> ¾ Lengths: <i>BD2100</i> Full length: <i>BD2110</i></p>		<p>3 PIECE OPERATING TABLE SET WITH CUTOUT</p> <p>3 Piece set Fits most tables Designed to evenly distribute patients' weight Gel pad can be rolled as a positioning aid Sheer reduction Foot section has cut-out</p> <p>Dimensions:</p> <p>Head section: 51 cm x 23 cm x 1cm Body section: 51 cm x 102 cm x 1cm Foot section: 51 cm x 51 cm x 1cm</p> <p><i>Reorder number:</i> <i>BD2130</i></p>
<p>MISCELLANEOUS</p>			
	<p>CHEST ROLLS / AXILLIARY ROLLS / FLAT BOTTOMS</p> <p>Various sizes and shapes</p> <p>Discuss with your representative your specific needs</p>		<p>SANDBAG POSITIONERS</p> <p>Various sizes and shapes</p> <p>Discuss with your representative your specific needs</p>

ABOUT SafMed

- B-BBEE Level 1 Contributor
- ISO 13485 Accredited
- Member of SAMED
- Medical Device Licence Numbers:
 - Cape Town - 00000243MD (Manufacture/Distribute)
 - Johannesburg - 00000227MD (Manufacture/Distribute)
 - Cape Town - 0000680MD (Wholesale)
 - Johannesburg - 00000643MD (Wholesale)



Safmed.SA



SafMed



@SafMedSA



SafMed South Africa



www.safmed.co.za