

SafMed Prolystica 2X Enzymatic Cleaner



PROLYSTICA 2X Concentrate Enzymatic Cleaner is a versatile enzymatic cleaning chemistry, effective on all types of soils and compatible with surgical instruments & flexible endoscopes. It has a dual enzyme formula which provides exceptional proteolytic performance within a broad range of water qualities and types. It is engineered to provide enhanced cleaning of fatty orthopaedic and protein soils in multiple applications.

Features:

- Low foaming
- Neutral pH
- Pleasant Fragrance
- Concentrated Low Dose (1-4ml/L)
- Enhanced corrosion inhibitor system
- Exceptional cleaning performance with speed in the removal of blood, fat, mucus and other surgical soils from medical devices in all water types