
The logo for SafMed features the word "SafMed" in a bold, sans-serif font. The letters "Saf" are blue, and the letters "Med" are red. The letter "M" is split vertically, with its left half being blue and its right half being red.

Your solutions partner for **30 YEARS**

SafMed Aseptor Bags

Instructions for use

Product description

Paper bags are the most economical type of sterilization packaging used worldwide. Bags are made in a range of sizes suitable for the packaging of most items. They conform to internationally recognized standards for paper quality, strength of side seals and process indicator ink performance. We certify that the Bags are Latex Free.

Printing:

Printing must show:

Bag Size; Manufacturing Date/Lot Number; Manufacturers name; sterilization process indicator (see below) and instructions for interpretation of colour change. The process indicator must be lead free ink which, under normal steam autoclave conditions change from a pink colour, to a brown colour.

Paper:

The base paper used is a bleached wet strength paper designed for conversion into sterilization bags and complies with ISO 11607-1, and SANS 50868-3:2010. Test reports from accredited Test Laboratories are provided.

Substance 60 +/- 3gsm

Microbial Barrier Approved DIN 58953-6

Water Repellency 20 seconds minimum

Finished bags:

Manufactured and tested to SANS 50868-3:2010 under ISO9001:2008. Test reports are provided.

Packing:

Bags must be packed in rigid cardboard cartons sealed at all edges and joins.

Indications for use

Chose the correct size of bag from the list shown below. Insert the contents so that sufficient space is left at the open end of the bag to allow the bag to be closed properly. To close the bag make a double fold of the open end and secure with Autoclave Tape.

Please ensure that items to be sterilized are double packed with an approved sterilization wrapping material.

SafMed Product Code	Tender Lot / Item	SafMed Product Description	Pack Size
BAGS001	NA	Aseptor Bag No 1 89x51x102mm	1000
BAGS002	NA	Aseptor Bag No 2 89x229mm	1000
BAGS003	NA	Aseptor Bag No 3 114x36x190mm	1000
BAGS004	Lot 10 / Item 53	Aseptor Bag No 4 127x48x260mm	1000
BAGS005	NA	Aseptor Bag No 5 127x85x229mm	1000
BAGS006	Lot 10 / Item 54	Aseptor Bag No 6 192x62x318mm	1000
BAGS007	Lot 10 / Item 55	Aseptor Bag No 7 178x74x395mm	1000
BAGS008	Lot 10 / Item 56	Aseptor Bag No 8 64x34x483mm	1000
BAGS009	NA	Aseptor Bag No 9 120x135mm	1000
BAGS011	Lot 10 / Item 51	Aseptor Bag No11 292x127x444mm	500
BAGS012	Lot 10 / Item 52	Aseptor Bag No12 330x152x559mm	500
BAGS013	NA	Aseptor Bag No13 257x108x381mm	1000
BAGS025	NA	Aseptor Bag No 25 64x36x152mm	1000
BAGS1020	NA	Aseptor Bag No1020 64x36x203mm	1000
BAGS3050	NA	Aseptor Bag No3050 64x36x254mm	1000