


# KEY SURGICAL INSTRUMENT CLEANING BRUSHES

## BOTTLE BRUSHES

- BR-06-062 (BR-3050)
- BR-18-062 (BR-3100)
- BR-24-070 (BR-3125)
- BR-24-080 (BR-3115)
- BR-18-087 (BR-3120)
- BR-18-105 (BR-3130)
- BR-18-125 (BR-3140)
- BR-18-250 (BR-3150)
- BR-24-312 (BR-3155)
- BR-06-375 (BR-3170)
- BR-24-440 (BR-3175)
- BR-18-880 (BR-3160)
- BR-18-500 (BR-4100)
- BR-06-438 (BR-2300)

## BURR BRUSHES



BR-1010

BR-1031

BR-1033

## TOOTHBRUSH STYLE BRUSHES




N-2121

N-2100

S-2100

ALL BRUSH DETAILS AND TIPS ARE SHOWN AS ACTUAL SIZE

## EXTRA-LONG SPECIALITY BRUSHES



BR-3308

BR-3310

BR-3320

BR-3330

BR-3030

BR-3040

BR-3410

BR-3420

BR-3300

## BRUSH SIZES AND DESCRIPTION

Instrument Cleaning Brushes (Gottle Brushes)	Brush Width/mm	Length Brush/mm	End-to-End Length/mm
NEW CODE (OLD CODE)			
Nylon Brushes	2	13	150
BR-06-062 (BR-3050)	2	13	450
BR-18-062 (BR-3100)	2	13	600
BR-24-070 (BR-3125)	2	25	800
BR-24-080 (BR-3115)	2	13	450
BR-18-087 (BR-3120)	2	13	450
BR-18-105 (BR-3130)	3	25	450
BR-18-125 (BR-3140)	3	50	450
BR-18-250 (BR-3150)	6	50	450
BR-24-312 (BR-3155)	6	50	450
BR-06-375 (BR-3170)	10	38	150
BR-24-440 (BR-3175)	10	75	500
BR-18-880 (BR-3160)	11	89	450
BR-18-500 (BR-4100)	13	75	450
BR-06-438 (BR-2300)	11	50	150

### BURR CLEANING BRUSHES

Autoclavable brushes with stainless steel wire bristles, specially designed for cleaning burrs, rasps, files or saws

BR-1010	Round Burr Brush	16	13	64
BR-1031	Single Row Burr Brush	3	35	165
BR-1033	Three-Row Burr Brush	16	35	165

### TOOTHBRUSH STYLE CLEANING BRUSHES

N-2121	Double Ended Nylon Bristles	6	29	130
N-2100	Three-Row Nylon Bristles	6	130	190
S-2100	Three-Row Stainless Steel Bristles	6	29	190

### EXTRA-LONG SPECIALITY CLEANING BRUSHES

BR-3308	Extralong brush with nylon bristle tip attached to rigid tubing	3.5	10	200
BR-3310	Bristles on one end only	3.5	20	200
BR-3320	Bristles on both ends	3.5	10	200
BR-3330	Bristles on both ends	3.5	20	200

BR-3030	Delicate pull-through cleaning brush, Nylon bristle tip attached to stainless steel wire	1	10	650
BR-3040	Bristles on one end only	2	10	1200
BR-3410	Bristles on one end only	3	10	1200
BR-3420	Bristles on one end only	3	10	1800

BR-3300	Extralong, Nylon bristles attached to stainless steel coil wire, loop handle at one end and angled tip at bristle end	4.6	10	2200
---------	---	-----	----	------



**SafMed**  
SAFMED (Pty) Ltd

JOHANNESBURG OFFICE: Cape Town Office: Durban Office:  
Tel: 011 201 4300 Tel: 021 763 3280 Cell: 082 464 3832  
Fax: 011 507 5550 Fax: 021 762 2760 Fax: 031 561 5138