

# SafMed CSSD Documentation System



The Printex Label Gun can print 2 print lines and uses labels that are 37x19mm. Each print line can be numerical or alphanumeric. Printex labels for CSSD documentation are made in double adhesive support and with Type 1 process indicators according to EN ISO 11140-1. The labels are produced with indicators for **DECONTAMINATION, STEAM, EO and PLASMA sterilization.**

The following information can be printed on the tracking labels:

- **Sterilization date**
- **Expiry date**
- **Operator**
- **Sterilization cycle number**
- **Batch/Lot number**
- **Product number**
- **Sterilizer code**