

**SafMed Trayliners**

SafMed Trayliners absorb and rapidly disperse condensate during sterilization. This greatly facilitates the production of dry packs and assists the sterilization process. When placed directly on transport trolley shelves they can help in the management of condensate. Trayliners can be also be placed under trays to protect wrap from sharp edges, this can reduce or prevent wrap from tearing that sometimes occurs during loading and unloading. Produced from pure cellulose wood pulp, manufactured with low chloride and sulphate content and free from optical brighteners.

**Features:**

- Facilitates the production of dry packs
- Assists the sterilization process
- Safe to use with delicate instruments
- Absorbs and rapidly disperses condensate
- Available in rolls and pre-cut sheets

To view all available sizes go to [www.safmed.co.za/absorbent-tray-liners.html](http://www.safmed.co.za/absorbent-tray-liners.html)