



**SafMed Self-Seal Sterilization Pouches** are pre-folded and assure accurate and fast closing without the requirement of a heat sealing device. They are especially suited for use in small general practices, dental surgeries and by users who want to avoid costs of additional sealing equipment.

Product Code	Description	Features
S1	90mm x 190mm (515)	<ul style="list-style-type: none"> <li>○ Pack Size: 200</li> <li>○ Constructed from transparent PET/PP multi-layers copolymer film, medical grade paper and adhesive strip which provides excellent resistance to over-exposure of steam and clean opening.</li> <li>○ Tear, fibre-free opening and aseptic presentation</li> <li>○ Triple band seal for higher package integrity</li> </ul>
S2	90mm x 230mm (520)	
S4	130mm x 250mm (525)	
S6	130mm x 360mm (530)	
S7	200mm x 330mm (545)	
S8	300mm x 400mm (555)	