

SafMed PROLYSTICA Multi-Enzymatic Cleaner



Prolystica Multi-Enzymatic Cleaner is a low foaming enzymatic detergent that contains 8 enzymes for use in manual, ultrasonic, and automated cleaning applications. This product may be used on surgical instruments, endoscopes, and other medical devices for the removal of blood, protein, carbohydrates, mucus, and lipids. Prolystica Multi-Enzymatic Cleaner is compatible with stainless steel, plastics, elastomers, and soft metals, including aluminum.

Prolystica Multi-Enzymatic Cleaner is uniquely formulated with biodegradable ingredients delivered to you within packaging made from 98 % recycled materials.

Features:

- Eight Enzyme Formulation
- Low Foaming
- Concentrated Low Dose (2mL to 7.5mL per of 1L of warm water)
- Enhanced corrosion inhibitor system
- Fast and powerful cleaning of soils
- Visual confirmation - Blue dye shows chemistry presence in the sink
- Better visibility of instruments under water with this low-foaming formulation