

SafMed Hinge Free Instrument Lubricant



Hinge-Free Instrument Lubricant is a non-silicone-based surgical instrument lubricant, sometimes known as instrument milk, designed for use with metal instruments. Regular use of Hinge-Free will keep surgical instruments operating freely and easily. Hinge-Free Instrument Lubricant contains an antimicrobial as a preservative, ensuring that the product will not support the growth of microorganisms.

The neutral pH of this surgical instrument lubricant provides compatibility with stainless steel, tungsten carbide, aluminium, and titanium.

Features:

- Excellent lubricating ability
- Non-silicone-based formula
- Contains corrosion inhibitors
- Neutral pH
- Contains an antimicrobial as a preservative
- Pack Size: 3.78 Litre Bottle & 5 Litre "Short" Bottle (Automated Washing)

[Click here to view the Instructions for Use.](#)