



SafMed Disposable Aprons are manufactured from “Virgin” Low Density Polyethylene (LDPE) and all aprons have a minimum thickness of 20 Microns. SafMed disposable aprons are intended for use in all relevant departments of a hospital. The aprons are provided on a cardboard hanger and are all sold in quantities of 200. There are 25 aprons attached to each cardboard hanger & 8 cardboard hangers per quantity of 200.



Various Colours and Sizes Available (Pack Size : 200)

- Econopron - White, 990mm Length
- Medium (Blue, Red, Green, Yellow) - 1100mm Length
- Medium White - 1100mm Length
- Long White - 1500mm Length