

# **HOSPI STERILIZERS**

## A DIVISION OF **SafMed**

# Product Catalogue

---



## Hospi JSD160 (160L) Autoclave with onboard 18kW boiler



The Hospi JSD160 is a 160 litre, trans atmospheric steam sterilizer available in a single door and double door configuration. Manufactured in South Africa, the Hospi JSD 160 offers a dedicated Bowie-Dick and vacuum leak test cycle as well as water saving functionality with the onboard water & energy saving system which utilizes recycled water. The Hospi JSD160 uses a water ejector to create a vacuum.

### Features:

- Trans atmospheric sterilization system – Class B.
- Chamber Size - 440 x 440 x 760.
- Sterilizes at both 121°C and 134°C.
- Dedicated Bowie & Dick test cycle.
- Vacuum leak test cycle - To test pressure vessel integrity.
- Electronic control system - PLC with touch screen.
- Water saving tank on-board - Recycled water.
- Displays cycle progress information.
- Printing capability.
- Four sterilization cycles to choose from:
  - Packs - 134 degrees
  - Flash - 134 degrees
  - Liquids - 121 degrees
  - Rubbers - 121 degrees

## Hospi JSD200 (200L) Autoclave with onboard 18kW boiler



The Hospi JSD200 is a 200 litre, trans atmospheric steam sterilizer available in a single door and double door configuration. Manufactured in South Africa, the Hospi JSD200 offers a dedicated Bowie-Dick and vacuum leak test cycle as well as water saving functionality with the onboard water & energy saving system which utilizes recycled water. The Hospi JSD200 uses a water ejector to create a vacuum.

### Features:

- Trans atmospheric sterilization system – Class B.
- Chamber Size - 460 x 460 x 960.
- Sterilizes at both 121°C and 134°C.
- Dedicated Bowie & Dick test cycle.
- Vacuum leak test cycle - To test pressure vessel integrity.
- Electronic control system - PLC with touch screen.
- Water saving tank on-board - Recycled water.
- Displays cycle progress information.
- Printing capability.
- Four sterilization cycles to choose from:
  - Packs - 134 degrees
  - Flash - 134 degrees
  - Liquids - 121 degrees
  - Rubbers - 121 degrees

## Hospi JSD400 (400L) Autoclave with onboard 36kW boiler



The Hospi JSD400 is a 400 litre, trans atmospheric steam sterilizer available in a single door and double door configuration. Manufactured in South Africa, the Hospi JSD400 offers a dedicated Bowie-Dick and vacuum leak test cycle as well as water saving functionality with the onboard water & energy saving system which utilizes recycled water. The Hospi JSD400 uses a liquid ring vacuum pump to create a vacuum.

### Features:

- Trans atmospheric sterilization system – Class B.
- Chamber Size - 660 x 660 x 915.
- Sterilizes at both 121°C and 134°C.
- Dedicated Bowie & Dick test cycle.
- Vacuum leak test cycle - To test pressure vessel integrity.
- Electronic control system - PLC with touch screen.
- Water saving tank on-board - Recycled water.
- Displays cycle progress information.
- Printing capability.
- Four sterilization cycles to choose from:
  - Packs - 134 degrees
  - Flash - 134 degrees
  - Liquids - 121 degrees
  - Rubbers - 121 degrees

## Hospi JSD600 (600L) Autoclave with onboard 36kW boiler



The Hospi JSD600 is a 600 litre, trans atmospheric steam sterilizer available in a single door and double door configuration. Manufactured in South Africa, the Hospi JSD600 offers a dedicated Bowie-Dick and vacuum leak test cycle as well as water saving functionality with the onboard water & energy saving system which utilizes recycled water. The Hospi JSD600 uses a liquid ring vacuum pump to create a vacuum.

### Features:

- Trans atmospheric sterilization system – Class B.
- Chamber Size - 660 x 660 x 1380.
- Sterilizes at both 121°C and 134°C.
- Dedicated Bowie & Dick test cycle.
- Vacuum leak test cycle - To test pressure vessel integrity.
- Electronic control system - PLC with touch screen.
- Water saving tank on-board - Recycled water.
- Displays cycle progress information.
- Printing capability.
- Four sterilization cycles to choose from:
  - Packs - 134 degrees
  - Flash - 134 degrees
  - Liquids - 121 degrees
  - Rubbers - 121 degrees

## Stainless Steel Products



### Stainless Steel Packing Table

- Packing table with bottom shelf
- 1800Lx 800W x 900H



### Stainless Steel Storage

- 1850H x 500W x 1000L
- 5 or 6 tier available
- Fixed and mobile options available



### Stainless Steel Working Tables

- Working Tables with bottom shelf
- Various sizes
- Splash back available on request



### Stainless Steel Shelves

- Various sizes available

## Water Filters & Softeners

**SPT 20 UNIT – Water softening system**



**Big Blue water filter system**

