



Operating Table
XHT-D


Electro-hydraulic Operating Table XHT-D

XHT-D is a high performance modular and motorized hydraulic operating table with outstanding features for any surgical disciplines. The complete structure, column and base cover are made of high quality 304 stainless steel, ensuring the safety and stability needed for complex demands of surgery.




Lowest Table Height

With lowest height up to 520mm for neurosurgery, surgeons can sit and perform the surgery with great comfort, maximizing the work efficiency and reducing whole patients recovery time.



Large Height Lift

A 510mm height range supports a variety of procedures, whether the surgeon is seated or standing. This will reduce strain and the risk of workplace injury.



One-button to leveling, Flex/Reflex

XHT-D is designed with one-button to leveling, flex/reflex to simplify medical staff's work load. Every seconds in OR is counted, which requires all instruments inside are user-friendly and time-saved.



Quick connection Interface

With the help of imported gas spring, the easy click interfaces of XHT-D surgical table allow simple and quick module changes. Various modules can be exchanged in no time at all.



Options

 Kidney board, motorized operation: 100mm



 Longitudinal slide: up to 300mm



 Carbon fiber tabletop



 Foot switch control



 Orthopedic Traction Frame



 Support accessories



Basic motion parameters

