

SCOPE BUDDY™

Endoscope Flushing Aid



CONSISTENT, HANDS-FREE FLUSHING OF ENDOSCOPES

REDUCES THE RISK OF REPETITIVE MOTION INJURIES

- Replaces repetitive and labour-intensive syringing of channels
- Simple, hands-free operation reduces muscle fatigue
- Provides constant fluid delivery through all endoscope channels

PROVIDES CONSISTENT FLUSHING

- Continuous, timed fluid delivery ensures consistent, repeatable endoscope flushing
- Audible and visual alerts upon completion
- Compatible with all major endoscope brands
- Works with your existing sink; compact size for countertop or wall mounting
- Meets or exceeds endoscope manufacturers requirements for channel flushing volumes

SHORTENS FLUSHING TIME

- Thoroughly flushes most endoscopes in 35 seconds
- Flushes multiple endoscope channels simultaneously
- Quickly and quietly delivers fluids through endoscope channels

INTERCEPT™ DETERGENT PROVIDES SUPERIOR ENDOSCOPE CLEANING

- Fast, one-minute contact time
- Low foaming, neutral pH formulation safe for endoscope use
- Cost-effective concentrated formula reduces per-use costs



THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, **only Cantel delivers the Complete Circle of Protection**, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimise your success.

MANUAL CLEAN Manual cleaning is essential in removing bioburden and preparing the endoscope for high-level disinfection. For optimal results, Cantel's manual cleaning products and detergents are designed to be a critical first step to efficiently and effectively eliminating bacteria on endoscopes prior to reprocessing.



Visual display
of endoscope
flushing time

User-selectable
time setting

Audible and
visual alerts
upon completion

One-touch
start/stop operation

Simple and quick
tubing connections

*SCOPE BUDDY™ is only
intended to be used as
a Flushing Aid during
Manual Cleaning step of
Endoscope Reprocessing.
It is not intended to
manually clean endoscopes.

Complete endoscope connection kit



Specifications

Unit dimensions	24.1 H x 15.2 W x 20.3 D (cm) (9.5 H x 6 W x 8 D (inches))
Unit weight	4.5 kg (10 lbs)
Electrical requirements	110 – 230 V, 50/60 Hz
Power requirements	230v 50/60Hz
Electrical safety certifications	UL Standard 61010-1 CANADA/USA C22.2 No. 61010-1 ETL and cETL listed
Fluid delivery rate	1 litre/minute
Flush time	0 – 99 minutes
Timer control	Digital display with touchpad controls
System pressure	Built-in safety bypass
Cycle completion indicators	Audible and visual alerts
Unit placement	Counter-top or wall-mounting (bracket included)
Warranty	One year

SCOPE BUDDY™ Flushing Aid information

PART NUMBER	DESCRIPTION
EFA-IN-G	SCOPE BUDDY™ Flushing Aid

Accessory ordering information

PART NUMBER	DESCRIPTION
78399-669	Complete endoscope connection kit

Detergent ordering information

PART NUMBER	DESCRIPTION
ML02-0134	INTERCEPT™ Detergent (4 x 3.8 litre)



SCOPE BUDDY™ and INTERCEPT™ are trademarks of Medivators Inc.
© 2017-2021 Cantel. All Rights Reserved. The Complete Circle of Protection diagram is a trademark owned by Cantel Medical Corp.
and is subject to copyright protection owned by Cantel Medical Corp.

www.cantelmedical.eu

Medivators BV | t: +31 (0)45 547 1471 | e: csdistributors@medivators.com

© 2021 Cantel Medical Corp. LIT-00372/A-ENEU

