



» Sterisystem Perfo-Safe® wire baskets

Our Sterisystem highly polished stainless steel Perfo-Safe® wire basket, which is available in both ISO and DIN sizes, is a safe and hygienic means to transport and store sterile packs.

As a sterilisation, transport and storage component, these stainless steel Sterisystem Perfo-Safe® wire baskets will prove to be an invaluable tool in your instrument set sterilisation process.

The Sterisystem Perfo-Safe® wire baskets can be used in all our transport carts and storage systems.

Sterisystem Perfo-Safe® wire baskets are invaluable if you want to ensure that your sterile sets of instruments are stored safely. Our Sterisystem Perfo-Safe® wire baskets form the basis of the CSSD basket-in-basket system.



» Stackable and nestable

Without the addition of extra handles, the baskets are stackable and nestable, making them more user-friendly and safer (no fear of fingers getting jammed).

For hygiene purposes, the upper edge features one single wire only.



» Outstanding quality

The baskets have been manufactured in accordance with the most stringent quality standards. Baskets with welding burrs, gaps, wire distortions or discolorations in the welding work never see the light of day.

Our baskets are devoid of any sharp edges or welding burrs which could cause tears in packaging.



» Safety

These Sterisystem Perfo-Safe® wire baskets are ideal if you wish to follow the basket-in-basket system to store your packaged Sterisystem Dry-Base® instrument trays. As the packaged instrument net is stored in a wire basket, direct contact with sterile instrument set is avoided.

» Sterisystem Perfo-Safe® wire baskets




» Sterisystem Perfo-Safe® wire baskets ISO

Status	Article Number	Name	Ext. Dimensions	Description
A	30005 00002	Anna	600 x 400 x 100 mm	-
A	30005 00004	Alice	600 x 400 x 130 mm	-
A	30005 00018	Amanda	400 x 300 x 100 mm	-
C	30005 00007	Annetta	600 x 400 x 200 mm	Upon request



» Sterisystem Perfo-Safe® wire baskets DIN

Status	Article Number	Name	Ext. Dimensions
A	30005 00001	Cindy	600 x 300 x 50 mm



» Silicone pad for Sterisystem Perfo-Safe® wire baskets and Sterisystem Perfo-Safe® wire shelves

Status	Article Number	Description	Thickness	Ext. Dimensions
A	30010 00004	Standard	3 mm	540 x 280 mm
A	30010 00006	Heavy duty	6 mm	600 x 324 mm