

## KEN IWD 2314

### Washer-disinfector for hospitals and clinics



- Machine for autoclave baskets
- Elegant and clean-cut Danish design
- Optimal hygiene level

For more information consult – [www.ken.dk](http://www.ken.dk)



#### Extra-spacious machine with limited outer measurements

KEN IWD 2314 is identical with IWD 2311, but with its extra 4 cm of depth it is also possible to wash and disinfect autoclave baskets and instrument containers at, for example, sterile centres. The machine is available both as a pass-through machine and with one door. The pass-through model ensures effective separation of an unclean and clean side, thus optimising the

level of hygiene and reducing the risk of hospital infections. The machine is elegant and clean-lined Danish design. KEN IWD 2314 can wash at up to 3 levels and, with its extra 4 cm of depth and large range of flexible equipment, it is a good choice for disinfecting washes at e.g. surgical wards, clinics for out-patients and sterile centres.



## Advantages:

- User-friendly control panel ensures correct operation
- Optional KEN NetCom or XML connect software solutions
- Individual programme selection possible. Four programme positions are standard
- Process Indication Bar shows how far the programme has come
- Double monitoring of temperature (EN/ISO 15883)
- Possible selection of intelligent, maintenance-free drying module
- Possible connection of up to 4 types of chemicals - 2 as standard
- Able to have printer connected for optimum process documentation
- Lockable lower compartment for keeping chemicals - up to 4 x 5 l containers



## Service

The machine has built-in KEN Quality control for optimum reliability. It demands a quality assurance check every 24 hours as well as servicing after 2,500 washes. The construction is service-friendly, with easy access to all components.

All KEN machines conform to present EU Directives, and all hospital machines have been constructed on the basis of EN ISO 15883, including, for example, the testing of its cleaning capacity according to SIS-TR 3:2001 (formerly SPRI).



## Accessories

Above are examples of accessories/equipment for this machine. Look for equipment with the aid of the equipment selector at [www.ken.dk](http://www.ken.dk), or contact KEN for information about further possible equipment, including customized special solutions.

## Contact a KEN specialist

KEN will gladly help you find the right solution. Consult [www.ken.dk](http://www.ken.dk) to find your KEN dealer.

Tlf.: +45 62 63 10 91

e-mail: [ken@ken.dk](mailto:ken@ken.dk)

Internet: [www.ken.dk](http://www.ken.dk)

## Specifications

Power connection	3 x 400 V,N,J, 50Hz or multiple choices available
Connected load/pre-fuse	10kW/16A
Water connection	Cold untreated water: 1 pcs., 3/4" RG, 5-30°C, 0-30°dH
Water connection	Hot soft water: 1 pcs., 3/4" RG, 30-70°C, 0-3°dH
Water connection	Demin. water: (Extra) 1 pcs., 3/4" RG 5-70°C
Water pressure	200-800 kPa
Outlet valve	Ø 50
Water amount (approx.) per. intake	13-14 l
Ventilation (if supplied with drying module)	Ø 70, 100 m³/hour
Pump capacity	370 l/min.
Free depth 2 doors	620 mm
Free depth 1 door	630 mm
Free width without central level	510 mm
Noise level	approx. 59-64 dBa (sound pressure)
Sealing class	IP 21 (not scavenger-proof)



**We make a difference – both locally and globally**