

Saff Control of Contr

Product Catalogue





Johannesburg Office: Tel: 011 201 4300 Fax: 011 507 5550

Durban Office: Tel: 082 464 3832 Fax: 031 561 5138 Cape Town Office:Tel:021 763 3280Fax:021 762 2760

E-mail: info@safmed.co.za Website: www.safmed.co.za www.facebook.com/Safmed.SA/ Twitter: @SafMedSA





We specialise in providing solutions and products that improve the productivity and effectiveness in the CSSD, Operating Theatre & Endoscopy Unit.

From expert advice on new theatre or CSSD designs to providing the latest digital operating theatre and instrument tracking technology, we have the products and experience to offer meaningful value to our customers.





Johannesburg Office:Tel:011 201 4300Fax:011 507 5550

Durban Office: Tel: 082 464 3832 Fax: 031 561 5138 Cape Town Office: Tel: 021 763 3280 Fax: 021 762 2760 E-mail: info@safmed.co.za Website: www.safmed.co.za www.facebook.com/Safmed.SA/ Twitter: @SafMedSA



Reliance 310 Sensitive - Crepe Paper

Item Code	Size (mm)	Pack Size
42281808-00011	1200 x 1700	100

SafMed Reliance crepe paper is made with more than 97% of cellulosic fibers and wet strength and sizing agents.

Cellulosic fibers create a tortuous path for excellent sterility maintenance and bacterial barrier. Excellent bacterial barrier. with Compatible EO and formaldehyde. Liquid Repellent. Cost sterilization steam, by effective. Softest grade of crepe paper available. Reliance 310 Sensitive is intended for single use applications and be considered as a Medical Device Class 1 Accessory wrapping fabric. Can be used as an inner and outer wrap for trays and packs. Manufactured to comply with ISO 11607 and EN 868.







Johannesburg Office: Durban Office: Tel: 011 201 4300 Fax: 011 507 5550

Tel: 082 464 3832

Cape Town Office: Tel: 021 763 3280 Fax: 031 561 5138 Fax: 021 762 2760

E-mail: info@safmed.co.za Website: www.safmed.co.za www.facebook.com/Safmed.SA/ Twitter: @SafMedSA



Reliance 440 - SMMMS Sterile Wrap Product

Item Code	Size (mm)	Pack Size
42281902-000334	880 x 1500	100
42281902-000335	1160 x 1200	100
42281902-000336	1160 x 1400	100

SafMed Reliance 440 SMS is 100% polypropylene and treated to resist static charge build-up. It does not contain natural fibers (cellulose) or synthetic binder.

Reliance 440 is intended for single use applications and be considered as a Medical Device Class 1 Accessory wrapping fabric. It could be used an inner or as an external layer. Compatible with steam (121-134C), as Ethylene Oxide, low temperature formaldehyde and Vapour Hydrogen Peroxide sterilization methods.







Johannesburg Office: Durban Office: Tel: 011 201 4300 Fax: 011 507 5550

Tel: 082 464 3832

Cape Town Office: Tel: 021 763 3280 Fax: 031 561 5138 Fax: 021 762 2760

E-mail: info@safmed.co.za Website: www.safmed.co.za www.facebook.com/Safmed.SA/ Twitter: @SafMedSA




Reliance Tandem 2 - Interleaved Sequential Wrap

Product	Item Code	Size (mm)	Pack Size
SMS Regular	42281902-00035	1160 x 1200	100 (50 x 2 sheets)
SMMS Regular	42281902-000355	1160 x 1200	100 (50 x 2 sheets)
SMMS Medium	42281902-000358	1160 x 1200	100 (50 x 2 sheets)
SMMS Heavy	42281902-000361	1160 x 1200	100 (50 x 2 sheets)
SMS Medium	42281902-00040	1160 x 1200	100 (50 x 2 sheets)
SMS Heavy	42281902-00045	1160 x 1200	100 (50 x 2 sheets)
SMMS Regular	42281902-000353	580 x 600	100 (50 x 2 sheets)
SMMS Medium	42281902-000356	580 x 600	100 (50 x 2 sheets)
SMMS Heavy	42281902-000359	580 x 600	100 (50 x 2 sheets)
SMMS Regular	42281902-000354	880 x 900	100 (50 x 2 sheets)
SMMS Medium	42281902-000357	880 x 900	100 (50 x 2 sheets)
SMMS Heavy	42281902-000360	880 x 900	100 (50 x 2 sheets)
SMS Heavy	42281902-00044	880 x 900	100 (50 x 2 sheets)

SafMed's Reliance Tandem 2 is a combination of two interleaved SMMMS sheets of nonwoven Reliance 435 Green and Reliance 435 Blue, providing a superior sterile barrier system. Both nonwovens deliver the barrier protection and sterilant permeability required for optimal sterilization of medical devices

SafMeds Reliance Tandem 2 is effectively lint and memory-free and is soft and strong. Reliance Tandem 2 is intended for single use applications and be considered as a Medical Device Class 1 Accessory wrapping fabric. Tandem 2 has been validated for use as a 2 layer sterile barrier system. The layers MUST NOT BE used separately.





Tel: 011 201 4300 Fax: 011 507 5550

Tel: 082 464 3832

Tel: 021 763 3280 Fax: 031 561 5138 Fax: 021 762 2760 Website: www.safmed.co.za www.facebook.com/Safmed.SA/ Twitter: @SafMedSA



Samed Aseptor Bags

Product	Item Code	Size (mm)	Pack Size
BAGS001	42281916-00014	89 x 51 x 102	1000
BAGS002	42281916-00015	89 x 229	1000
BAGS003	42281916-00016	114 x 36 x 190	1000
BAGS004	42281916-00017	127 x 48 x 260	1000
BAGS005	42281916-00018	127 x 85 x 229	1000
BAGS006	42281916-00019	192 x 62 x 318	1000
BAGS007	42281916-00020	178 x 74 x 380	1000
BAGS008	42281916-00021	64 x 34 x 483	1000
BAGS009	42281916-00022	120 x 135	1000
BAGS011	42281916-00023	292 x 127 x 444	500
BAGS012	42281916-00024	330 x 152 x 559	500
BAGS1020	42281916-00025	864 x 36 x 203	1000

SafMed Aseptor Bags are locally manufactured from a special medical grade paper which forms a bacteria barrier and are used to enclose items to be sterilized. It has a distinctive process indicator which changes colour when exposed to steam cycles. All SafMed Aseptor Bags conform to internationally recognized standards ISO 11607-1 and tested to SANS 50868-3:2010. SafMed Aseptor Bags are meant for use in Steam Sterilization Cycles only.

Features:

- High paper quality & strength of side seals
- Process indicator ink performance.
- Manufactured in South Africa
- Available in a wide range of sizes



Johannesburg Office: Durban Office: Tel: 011 201 4300 Fax: 011 507 5550

Tel: 082 464 3832

Cape Town Office: Tel: 021 763 3280 Fax: 031 561 5138 Fax: 021 762 2760 E-mail: info@safmed.co.za Website: www.safmed.co.za www.facebook.com/Safmed.SA/ Twitter: @SafMedSA

RT28-1-2019 Catalogue

Saffed

Samed Aseptor Bags









8

















Johannesburg Office: Durban Office: Tel: 011 201 4300 Fax: 011 507 5550

Tel: 082 464 3832 Fax: 031 561 5138

Cape Town Office: Tel: 021 763 3280 Fax: 021 762 2760

E-mail: info@safmed.co.za Website: www.safmed.co.za www.facebook.com/Safmed.SA/ Twitter: @SafMedSA



Samed CSSD Documentation System

Item Code	Item Description	Pack Size
42281902-000338	Labelling Device for use with steam & LTS Sterilization	1
42281902-000342	Documentation Envelopes	100
42281902-000339	Tracking labels for Low Temp Sterilization	500
42281902-000340	Tracking labels for Steam Sterilization	650
42281902-000341	Tracking label for Decontamination for Washer	650



The Printex Label Gun can print 2 print lines and uses labels that are 37x19mm. Each print line can be numercial or alphanumerical. Printex labels for CSSD documentation are made in double adhesive support and with Type 1 process indicators according to EN ISO 11140-1. The labels are produced with indicators for DECONTAMINATION, STEAM, EO and PLASMA sterilization.

The following information can be printed on the tracking labels:

- Sterilization date
- Expiry date
- Operator
- Sterilization cycle number
- Batch/Lot number
- **Product number**
- Sterilizer code



Johannesburg Office: Durban Office: Tel: 011 201 4300 Fax: 011 507 5550

Tel: 082 464 3832 Tel: 021 763 3280 Fax: 031 561 5138 Fax: 021 762 2760

Cape Town Office:

E-mail: info@safmed.co.za Website: www.safmed.co.za www.facebook.com/Safmed.SA/ Twitter: @SafMedSA



ABOUT Samed

- B-BBEE Level 1 Contributor ISO
- > 13485:2016 Accredited
- Member of SAMED
- Medical Device Licence Numbers:
 - Cape Town 00000243MD (Manufacture/Distribute)
 - Johannesburg 00000227MD (Manufacture/Distribute)
 - Cape Town 0000680MD (Wholesale)
 - Johannesburg 00000643MD (Wholesale)

